INSTRUCTIONS FOR USE

XENIOS COMPACT HOLDER





CE

XENIOS


A FRESenius
MEDICAL CARE
COMPANY

SYMBOL GLOSSARY


| | |
|--|---|
|  | <p>Consult Instructions for Use Gebrauchsanweisung beachten استخدام بالتعليمات الرجوع إلى Čtěte návod k použití Se brugsanvisningen Consultar las instrucciones de uso Katso käyttöohjeet Consulter la notice d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Olvassa el a használati utasítást Consultare le istruzioni per l'uso 사용 설명서 참고 Lees de gebruiksaanwijzing Przestrzegaj instrukcji użycia Consultar as instruções de utilização Consultați instrucțiunile de utilizare Обратитесь к инструкции по эксплуатации Se bruksanvisningen Kullanma Talimatlarına Bakın 參考使用說明</p> |
|  | <p>Serial number Seriennummer الرقم التسلسلي Sériové číslo Seriennummer Sarjanumero Numéro de série Σειριακός αριθμός Sorozatszám Numero di serie 일련 번호 Seriennummer Numer seryjny Número de série Număr de serie Серийный номер Seriennummer Seri numarası 序號</p> |
|  | <p>Manufacturer Hersteller الشركة المصنعة Výrobce Producent Número de serie Fabricante Valmistaja Fabricant Κατασκευαστής Gyártó Costruttore 제조자 Fabrikant Producent Fabricante Producător Производитель Tillverkare Üretici 製造商</p> |
|  | <p>Date of manufacture Herstellungsdatum تاريخ التصنيع Datum výroby Fremstillingsdato Fecha de fabricación Valmistuspäivä Date de fabrication Ημερομηνία κατασκευής Gyártás dátuma Data di produzione 제조연월 Productiedatum Data produkcji Data de fabrico Data de fabricație Дата изготовления Tillverkningsdatum Üretim tarihi 生產日期</p> |
|  | <p>Catalog number Artikelnummer رقم الكتالوج Datum výroby Katalogové číslo Katalognummer Referencia Tuotenumero Référence produit Katalógusszám Numero di catalogo 주문번호 Catalogusnummer Numer katalogowy Referencia de catálogo Număr de catalog Номер по каталогу Artikelnummer Katalog numarası 目錄編號</p> |
|  | <p>Keep dry Trocken aufbewahren يُحفظ جافاً Uchovávejte v suchu Beskyttes mod fugt Mantener seco Säilytettävä kuivassa Conserver au sec Αριθμός καταλόγου Να διατηρείται στεγνό Száras helyen tárolandó Conservare in luogo asciutto 건조하게 유지 Droog bewaren Chronić przed wilgocią Manter em local seco A se depozita într-un spațiu uscat Беречь от влаги Förvaras torr Kuru yerde saklayın 保持乾燥</p> |


6



7

8

LANGUAGES INCLUDED IN THIS INSTRUCTIONS FOR USE

| | |
|----------------------|-----|
| English | 8 |
| Deutsch | 18 |
| العربية | 26 |
| Čeština | 34 |
| Dansk | 42 |
| Español | 50 |
| Suomi | 58 |
| Français | 66 |
| ΕΛΛΗΝΙΚΑ | 74 |
| Magyar | 82 |
| Italiano | 90 |
| 한국어 | 98 |
| Nederlands | 106 |
| Polski | 114 |
| Português | 122 |
| Română | 132 |
| Русский | 140 |
| Svenska | 148 |
| Türkçe | 156 |
| 繁體中文 | 164 |

Copyright

The copyright to these Instructions for Use is held by Xenios AG. All rights reserved, including for translation. No part of this information may be reproduced or electronically processed, copied, translated or disseminated in any form (printed, photocopy, microfilm or by a different process) either in whole or in part without the written permission of Xenios AG. All violations involving copyright can result in legal action.

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Item number bracket: 38350070

CONTENTS

| | | |
|-------|--|----|
| 1. | Introduction | 11 |
| 1.1. | Explanation of Symbols and Notices. | 11 |
| 2. | Intended Use | 11 |
| 3. | Contraindications | 12 |
| 4. | Safety Instructions | 12 |
| 5. | Set up | 13 |
| 5.1. | Display Safety Plate. | 13 |
| 5.1.1 | Attaching the Display Safety Plate | 13 |
| 6. | Application | 14 |
| 6.1. | Handling the Clamping Screw | 14 |
| 6.2. | Attach the Compact Holder | 14 |
| 6.2.1 | Attach the Compact Holder to the Rear of the Console | 14 |
| 6.3. | Attach the Compact Holder to a Standard Medical Device Rail. | 15 |
| 6.4. | Fitting the Pin or Bracket | 15 |
| 6.5. | Mounting the Gas Exchanger | 16 |
| 6.5.1 | Mounting the Oxygenator of the iLA active kits / Novalung kits / Support sets, i-cor puls TC kit | 16 |
| 6.5.2 | Mounting the iLA / i-cor Membrane Lung | 16 |
| 6.6. | Fastening the Pump Drive | 16 |
| 6.7. | Adjusting the Pump Drive Inclination Angle | 16 |
| 6.8. | Attaching the Sensor Box | 17 |
| 7. | Cleaning and Maintenance | 17 |
| 8. | Packaging and Storage. | 17 |

1. INTRODUCTION

These Instructions for Use are provided for all medical staff required to work with the Xenios compact holder and should be read prior to use.

The compact holder may only be used by specially trained and qualified medical staff. The compact holder is a Class I device according to the Medical Device Directive 93/42/EEC.

If not listed separately, all information in the Instruction for Use referring to the console also apply to the deltaxstream MDC console, iLA active console, i-cor console and Xenios console.

1.1. EXPLANATION OF SYMBOLS AND NOTICES

DANGER DANGER indicates a hazardous situation which, if not avoided, will result in death or serious injury.

WARNING WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION CAUTION indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

2. INTENDED USE

The compact holder is used for modular attachment of different types of compatible gas exchanger (see section 4), pump drives and sensor boxes to the trolley's rail or to an appropriate standard medical device rail.

When the console is used without the trolley (e.g. for inner-hospital transport) the compact holder can be attached to the back of the console. In this setup, the system can also be used as a table unit.

CONTRAINDICATIONS

See figure **1**:

| ID | Item | Weight | Item number |
|----|-----------------------|--------|-------------|
| 1 | Xenios compact holder | 1520 g | 30000145 |
| 2 | Pin (standard) | 260 g | 38350041 |
| 3 | Bracket (optional) | 500 g | 38350070 |

3. CONTRAINDICATIONS

No contraindications are known.

4. SAFETY INSTRUCTIONS

- The compact holder may only be used to mount the following equipment: pump drive, gas exchanger, pin, bracket, sensor box.
- The compact holder may only be used for the following products with the respective compact holder:

| Pin | Bracket |
|---|--|
| Novalung kits <ul style="list-style-type: none"> – XLung – XLung kit 230 – MiniLung kit – MiniLung petite kit | iLA product family <ul style="list-style-type: none"> – iLA active iLA kit – iLA active iLA kit IPS – iLA membrane ventilator – i-cor puls kit |
| iLA active kits <ul style="list-style-type: none"> – iLA active XLung kit – iLA active MiniLung kit – iLA active MiniLung petite kit | |
| Support sets | |
| i-cor puls TC kit | |

- The user bears sole liability for the use of the compact holder with any other than these approved sets.

- Ensure that the compact holder and the individual components are mounted and fixed securely and screws are tightened. All screws shall be tightened and loosened by hand only.
- Attach components in a way that prevents kinking of the tubing.
- Do not expose the compact holder to high mechanical stress.
- Do not use the compact holder if it is damaged (e.g. damages on screws, and holder frame).
- If the compact holder is dropped or falls down, the compact holder and all fastening components must be checked for damage. The damaged compact holder or components should not be used further and must be exchanged.

5. SET UP

5.1. DISPLAY SAFETY PLATE

Before the compact holder is used, the display safety plate must be attached to the console.

WARNING The compact holder can only be attached to a console that is equipped with a display safety plate.

CAUTION The safety plate must be securely and firmly fastened and must no longer be detachable by hand.

5.1.1 ATTACHING THE DISPLAY SAFETY PLATE

CAUTION The safety plate must be securely and firmly fastened and must no longer be detachable by hand.

To prevent electric shock, disconnect power supplies prior to cleaning the device.

The display safety plate must be attached to the top of the console.

Proceed as follows [figure 2]:

Using either a mild detergent solution or a common disinfectant (e.g. Isopropanol) to remove any dirt and grease residue from the area to which the display safety plate will be attached.

Peel the paper off of the adhesive surface on the reverse side of the display safety plate. Affix the display safety plate to the cleaned area. Ensure that the display safety plate is firmly fastened to the console.

6. APPLICATION

6.1. HANDLING THE CLAMPING SCREW

After the clamping screw has been tightened, it can be returned to its starting position and tightened further if necessary. To do this, pull the handle of the clamping screw out and return it to its starting position [figure 3].

6.2. ATTACH THE COMPACT HOLDER

CAUTION The compact holder must only be fastened to an appropriate standard medical device rail or the rear of the console.

6.2.1 ATTACH THE COMPACT HOLDER TO THE REAR OF THE CONSOLE

Hook the compact holder onto the handle extension [figure 4] while connecting it to the rail on the rear of the power supply. Make sure that it is securely connected to the rail. Tighten the clamping screw without tools. Make sure that the compact holder is securely attached to the power supply [figure 5].

WARNING

Always place the console on a stable and level surface. Be aware that the weight of the fully equipped compact holder on the back of the console is sufficient to cause the console to tilt or fall if positioned incorrectly. In case of tilting or falling the user and/or the patient and/or (disposable) components may be harmed.

6.3. ATTACH THE COMPACT HOLDER TO A STANDARD MEDICAL DEVICE RAIL

Attach the compact holder to a standard medical device mounting rail (10 mm x 25 mm) such as is located on the trolley [figure 6] using the orange coloured clamp located on the rear side of the compact holder. Once the compact holder is positioned correctly, tighten the clamping screw so that the compact holder can no longer be moved. Prior to attempting subsequent assembly steps, ensure that the compact holder is securely attached to the trolley or standard medical device rail.

CAUTION

Before attaching the compact holder, make sure the standard medical device rail is stable.

6.4. FITTING THE PIN OR BRACKET

Insert the pin or bracket into the base carrier and fasten it with the clamping screw [figure 7]. Make sure that the screw of the bracket is visible from the front.

6.5. MOUNTING THE GAS EXCHANGER

WARNING

Mount the gas exchanger without excessive mechanical stress, as described in these Instructions for Use of the specific disposable.

Mount the gas exchanger onto the holder in a way that the gas inlet is at the highest and the gas outlet at the lowest point.

Position the gas exchanger in a way that the blood inlet of the gas exchanger is higher than the blood outlet of the pump head.

6.5.1 MOUNTING THE OXYGENATOR OF THE ILA ACTIVE KITS / NOVALUNG KITS / SUPPORT SETS, I-COR PULS TC KIT

Place the oxygenator onto the pin, lock the oxygenator into place as far down as possible on the pin and ensure the tubing is free of kinks. After tightening the screws of the compact holder, the oxygenator should no longer rotate on the pin [figure 8].

6.5.2 MOUNTING THE ILA / I-COR MEMBRANE LUNG

Turn the screw to open the bracket until the membrane lung can be inserted. To fix the membrane lung onto the bracket close the screw [figure 9]. Ensure that the bracket engages the ribs of the membrane lung [figure 10]. For in-hospital transport, the membrane lung can be lowered via the integrated height adjuster of the compact holder for better protection of the oxygenator.

6.6. FASTENING THE PUMP DRIVE

Attach the pump drive to the top of the compact holder. Make sure that the pump drive is fully seated in the compact holder and then tighten the clamping screw by hand [figure 11].

6.7. ADJUSTING THE PUMP DRIVE INCLINATION ANGLE

The angle of inclination of the pump drive can be adjusted. To do this, open the clamping screw and bring the pump drive into the desired position. Fix the desired position by tightening the clamping screw [figure 12].

6.8. ATTACHING THE SENSOR BOX

Attach the sensor box to the outer right side of the compact holder. To do this, use the fastener on the rear of the sensor box. Ensure that the sensor box is attached in a way that it cannot slip off the side of the compact holder.

7. CLEANING AND MAINTENANCE

The compact holder is used for modular attachment of different types of compatible oxygenators (see section 4), pump drives and sensor boxes. To maintain reliable functionality and prevent infection, contamination, or malfunction, the compact holder must be cleaned and checked on intactness before each use. For disinfection, use a medical grade disinfectant (e.g. Bode Chemie GmbH Kohrsolin, Sterilium, Bacilol Plus). For further approved cleaning solutions refer to console Instructions for Use.

Maintenance of the compact holder is recommended as part of the console's maintenance interval. If any problem occurs, contact the manufacturer.

8. PACKAGING AND STORAGE

Store and transport the compact holder in clean and dry condition.

The manufacturer shall not be held liable for damage caused by incorrect or improper storage. Only use the original packaging for transportation.

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