
INSTRUCTIONS FOR USE

MULTISUPPORT GROUND

CE

XENIOS | A FRESENIUS
MEDICAL CARE
COMPANY

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IFU Version:
CE_003_2021_02

Item number IFU:
38350331



Should questions arise or problems/malfunctions occur, note the affected serial number and contact Xenios AG.





Technical Service group +49 7131 2706-234
technicalsupport@xenios-ag.com





Customer Service group +49 7131 2706-100
customerservice@xenios-ag.com

If the MultiSupport GROUND is used incorrectly or for purposes other than those defined in these Instructions for Use, the sole responsibility is with the user.

These Instructions for Use are provided for all qualified medical staff required to work with the MultiSupport GROUND. For safe and proper use, follow these instructions. Keep them for future reference.

SYMBOL GLOSSARY

	<p>Consult Instructions for Use Gebrauchsanweisung beachten راجع تعليمات الاستخدام Čtěte návod k použití Se brugervejledningene Συμβουλευτείτε τις οδηγίες χρήσης Consultar las instrucciones de uso Katso käyttöohjeet Consulter la notice d'utilisation Olvassa el a használati utasítást Consultare le istruzioni per l'uso 사용 설명서 참고 Lees de gebruiksaanwijzing Patrz Instrukcja użycia Consultar as Instruções de utilização Consultați instrucțiunile de utilizare Обратитесь к инструкциям по применению Läs bruksanvisningen Kullanim Talimatlarına bakın 參考使用說明</p>
	<p>Follow Instructions for Use Gebrauchsanweisung befolgen اتبع تعليمات الاستخدام Postupujte dle návodu k použití Følg brugervejledningene Εφαρμόζετε τις οδηγίες χρήσης Seguir las instrucciones de uso Noudata käyttöohjeita Suivre la notice d'utilisation Kövesse a használati utasítást Seguire le istruzioni per l'uso 사용 설명서에 따를 것 Volg de gebruiksaanwijzing Należy postępować zgodnie z Instrukcją użycia Seguir as Instruções de utilização Urmați instrucțiunile de utilizare Соблюдайте инструкции по применению Följ bruksanvisningen Kullanim Talimatlarına uyun 遵守使用說明</p>
	<p>Serial number Seriennummer الرقم التسلسلي Sériové číslo Seriennummer Σειριακός αριθμός Número de serie Sarjanumero Numéro de série Sorozatszám Numero di serie 일련 번호 Seriennummer Numer seryjny Número de série Număr de serie Серийный номер Seriennummer Seri numarası 序號</p>
	<p>Manufacturer Hersteller الجهة المصنعة Výrobce Producent Κατασκευαστής Fabricante Valmistaja Fabricant Gyártó Produttore 제조자 Fabrikant Producent Fabricante Producător Производитель Tillverkare Üretici 製造商</p>

	<p>Date of manufacture Herstellungsdatum تاريخ التصنيع Datum výroby Fremstillingsdato Ημερομηνία κατασκευής Fecha de fabricación Valmistuspäivämäärä Date de fabrication Gyártási dátum Data di fabbricazione 제조일 Productiedatum Data produkcji Data de fabrico Data fabricației Дата изготовления Tillverkningsdatum Üretim tarihi 生產日期</p>
	<p>Catalog number Artikelnummer رقم الكتالوج Katalogové číslo Varenummer Αριθμός καταλόγου Número de referencia Luettelonumero Référence produit Cikkszám Numero di catalogo 주문번호 Catalogusnummer Numer katalogowy Referência de catálogo Număr de catalog Номер по каталогу Artikelnummer Katalog numarası 目錄編號</p>
	<p>Keep dry Trocken aufbewahren يُحفظ جافاً Uchovávejte v suchu Beskyttes mod fugt Να διατηρείται στεγνό Mantener seco Säilytettävä kuivassa Conserver au sec Száraz helyen tárolandó Proteggere dall'umidità 건조한 곳에 보관 Droog bewaren Chronić przed wilgocią Guardar em local seco Păstrați produsul uscat Хранить в сухом месте Förvaras torrt Kuru tutun 保持乾燥</p>
	<p>Atmospheric pressure limitation Umgebungsdruckbegrenzung حد الضغط الجوي Omezení atmosférického tlaku Grænse for atmosfærisk tryk Περιορισμός ατμοσφαιρικής πίεσης Límite de presión atmosférica Ilmanpainerajoitus Limite de pression atmosphérique Léggöri nyomás tartománya Limitazione della pressione atmosferica 기압 제한 Limiet atmosferische druk Ograniczenie ciśnienia atmosferycznego Limite de pressão atmosférica Limitarea presiunii atmosferice Ограничение атмосферного давления Atmosfärisk tryckbegränsning Atmosferik basınç sınırlaması 大氣壓力限制</p>

Humidity limitation | Luftfeuchte, Begrenzung | حد الرطوبة | Omezení vlhkosti | Fugtighedsgrænse | Περιορισμός υγρασίας | Límite de humedad | Kosteusrajoitus | Limite d'humidité | Páratartalom-tartomány | Limitazione dell'umidità | 습도 제한 | Toelaatbare luchtvochtigheid | Ograniczenie wilgotności | Limite de humidade | Limitarea umidității | Ограничения по влажности окружающей среды | Luftfuktighetsgräns | Nem sınırlaması | 溼度限制

Temperature limit | Temperaturbegrenzung | حد درجة الحرارة | Teplotní limit | Temperaturgrænse | Όριο θερμοκρασίας | Límite de temperatura | Lämpötilaraja | Limite de température | Hőmérsékleti határérték | Limite di temperatura | 온도 한계 | Temperatuurlimiet | Ograniczenie temperatury | Limite de temperatura | Limită de temperatură | Температурные ограничения | Temperaturgräns | Sıcaklık limiti | 溫度限制

Inter-hospital transport (ground-based) | Interklinischer Transport (bodengebunden) | النقل بين المستشفيات (الأرضي) | Přeprava mezi nemocnicemi (pözemní) | Transport mellem hospitaler (landbaseret) | Μεταφορά μεταξύ νοσοκομείων (επίγεια) | Transporte interhospitalario (por tierra) | Sairaaloiden välinen kuljetus (maakuljetus) | Transport inter hospitalier (au sol) | Kórházak közötti szállítás (szárazföldön) | Trasporto interospedaliero (a terra) | 병원 간 이동(지상 기반) | Transport tussen ziekenhuizen (over de grond) | Transport pomiędzy szpitalami (naziemny) | Transporte inter-hospitalar (terrestre) | Transport interspitalicesc (pe sol) | Межбольничное средство транспортировки (наземное) | Transport mellan sjukhus (markbaserad) | Hastaneler arası taşıma (kara yolu ile) | 醫院間運輸(地面)

LANGUAGES INCLUDED IN THIS INSTRUCTIONS FOR USE

English	10
Deutsch	56
العربية	102
Čeština	150
Dansk	196
Ελληνικά	242
Español	288
Suomi	334
Français	380
Magyar	426
Italiano	472
한국어	518
Nederlands	562
Polski	608
Português	656
Română	702
Русский	748
Svenska	796
Türkçe	842
繁體中文	888

EN

DE

AR

CS

DA

EL

ES

FI

FR

HU

IT

KO

NL

PL

PT

RO

RU

SV

TR

ZH

CONTENTS

1.	Explanation of symbols	12
1.1.	Glossary	12
2.	General safety	14
2.1.	Handling the MSG Holder	16
2.2.	Handling the MSG Crossbar	17
2.3.	Handling the MSG Base Plate	17
3.	Contents liability	18
4.	Target group	18
4.1.	Requirements	18
5.	Intended use	19
5.1.	Indications	20
5.2.	Contraindications	20
6.	Scope of delivery	21
7.	Packaging and storage	21
8.	Overview of the MSG Components	22
8.1.	MSG Holder	22
8.2.	MSG Base Plate	24
8.3.	MSG Crossbar	24
8.4.	MSG Backup Module	25
8.5.	External Battery Charger	25
9.	Application	26
9.1.	Basic information on patient transport	26
9.2.	Safety during patient transport	27
9.3.	Mounting the MSG Holder	28
9.4.	Mounting the MSG Pin or MSG Bracket	29
9.4.1.	Fastening the pin or bracket	29
9.4.2.	Removing the pin or bracket	29
9.5.	Mounting the sensor box	30
9.6.	Mounting the pump drive	30

9.7.	Installing patient kits on the MSG Holder	31
9.8.	Mounting the MSG Cossbar	31
9.8.1.	Mounting the console to the MSG Crossbar	32
9.8.2.	Removing the console from the MSG Crossbar	32
9.9.	Mounting the MSG Base Plate.	33
9.9.1.	Mounting the MSG Base Plate on two airline rails (recommended).	33
9.9.2.	Mounting the MSG Base Plate on one airline rail (optional)	33
9.9.3.	Mounting the console on the MSG Base Plate	34
9.9.4.	Removing the console from the MSG Base Plate	34
9.10.	Connecting the console to the AC power connection on board of the ambulance	34
9.11.	Charging the battery packs	36
9.12.	Deploying the MSG Backup Module	36
9.12.1.	Before use	37
9.12.2.	After use	37
10.	Cleaning	37
11.	Maintenance	38
12.	Troubleshooting	39
12.1.	Ensure backup operation	39
12.2.	Alarms - measures	40
13.	Disposal	51
14.	Technical data	51
14.1.	MSG Holder	51
14.2.	MSG Crossbar	52
14.3.	MSG Base Plate	52
14.4.	MSG Backup Module	52
14.5.	MSG ambient conditions.	53
15.	Checklist for transport	54

1. EXPLANATION OF SYMBOLS

DANGER

DANGER indicates a hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

CAUTION indicates a hazardous situation which, if not avoided, can result in minor or moderate injury.

NOTICE

ATTENTION indicates a situation which, if not avoided, can result in damage to materials and impaired function of the product.

Note: Supplementary information on the operation of the product.

1.1. GLOSSARY

Term	Explanation
Console	The term console includes the following components: <ul style="list-style-type: none"> — Pump drive — Sensor box — Control panel — Power supply with two integrated battery packs
MSG	MultiSupport GROUND transport device for the transport of the console with its components and a patient kit.

Term	Explanation
Patient kit	<p>A patient kit consists of the following components:</p> <ul style="list-style-type: none"> – Gas exchanger – Pump head – Tubing lines <p>Patient kits are not reusable. A patient kit shall only be used for a single patient.</p>
Patient stretcher	e.g. Intensive-Transport-Unit ITU1, or a comparable product for the transport of intensive care patients in an intensive care unit
Pump drive	Pump drive of the console
Sensor box	Sensor box (without ECG function) of the console
System	The system described in this manual consists of a console, patient kit as well as the components of the MSG.
AC	Alternating electrical current (Symbol: ~)
DC	Direct electrical current (Symbol: ==)
Ground transport	Ground-based inter-hospital transport
Stationary application	Intra-hospital transport
ECMO/ECLS	Extracorporeal membrane oxygenation / extracorporeal life support

2. GENERAL SAFETY

Before each use of the components (crossbar, holder and/or base plate) the user must check the components for damage.

Danger to life due to console failure.

The console can fail in cases of moisture ingress.

Cover the console with the MSG Rain Cover as soon as the transport takes place in an uncovered outdoor area. The console must not be operated for more than 30 minutes under the Rain Cover.

DANGER

Danger to life due to failure of console components.

If the display, pump drive or battery pack fails, the patient's supply can be compromised.

To maintain gas exchange and extracorporeal circuit support in an emergency, always carry the MSG Backup Module with you during transport and keep it within easy reach of the console.

Reassemble the complete MSG Backup Module after each use.

Charge the battery packs after each use with the MSG Battery Charger.

WARNING

Danger to life due to console failure.

Unfavorable influencing factors, such as those that can occur during airborne transport, can lead to system failure (complete failure of the console).

The MSG is not approved for airborne transport.

Only use the MSG in ground transport.

Danger to life due to incompatible patient kits.

Incompatible patient kits can interfere with or cause the failure of the gas exchange or circulatory support for the patient.

Only use the following MSG patient kits for ground-based inter-hospital transport:

- XLung kit 230
- XLung kit
- MiniLung kit
- MiniLung petite kit

Other patient kits may not be used for ground-based inter-hospital transport.

WARNING**Risk of fatal injury due to damaged tubing.**

Damaged or kinked tubing can interrupt therapy and endanger the patient.

The transport of patients under ECMO/ECLS generally leads to an increased risk for the treatment and for the patient himself. Before transport begins:

- Carefully secure the vascular accesses (fixation of the cannulas).
- Secure the tubing lines of the extracorporeal circuit in such a way that they cannot get caught in other objects during transport.
- Secure the tubing lines of the extracorporeal circuit in such a way that no kinks can occur.
- Position the MSG Holder, including gas exchanger and pump drive, below the patient's heart level.
- Correctly position and secure the MSG Holder as well as the console at all times during transport.

CAUTION

Risk of injury from falling components.

The components can loosen and fall off if they are exposed to strong vibrations.

Secure the components of the MSG during transport.

This applies in particular to:

- the MSG Holder on a standard rail of a patient stretcher
- the MSG Crossbar on the infusion poles of a patient stretcher
- the MSG Base Plate on the airline rails

2.1. HANDLING THE MSG HOLDER

WARNING

Danger to life due to system component failure.

Inadequately secured components can interfere with or cause the failure of the gas exchange or circulatory support for the patient.

During transport, the MSG Holder must be attached to the standard rail of the patient stretcher.

Danger to life from damaged cables and tubing.

Damaged or kinked tubing can interrupt therapy and endanger the patient.

Secure the tubing lines and cables for transport in such a way that they cannot be damaged.

CAUTION

Risk of crushing due to falling parts.

The pump drive can fall out of its dedicated holder when the adjustable clamping handle is opened.

Hold the pump drive in place when the adjustable clamping handle of the pump drive holder is opened.

2.2. HANDLING THE MSG CROSSBAR

DANGER

Risk of crushing due to falling parts.

The MSG Crossbar can slide down or fall off if subjected to strong vibrations.

Use the MSG Crossbar exclusively for intra-hospital transport (from the intensive care unit to the ambulance and back).

The console must not be mounted on the MSG Crossbar during transport in the ambulance.

2.3. HANDLING THE MSG BASE PLATE

WARNING

Danger to life due to console failure.

The console can fail if subjected to strong vibrations.

During transport in an ambulance, the console must be fastened to the base plate.

Ensure that the console attached to the floor of the ambulance is not damaged by accidental kicks. A strong jolt or vibration can cause the following problems:

- the console display is damaged
- the treatment parameters are not displayed
- the treatment cannot be continued

CAUTION

Risk of crushing due to uncontrolled movements of the components.

The console can move uncontrolled in the ambulance if it is not mounted correctly.

Fasten the console to the MSG Base Plate during transport in the ambulance. To do this, fasten the console properly with the clamping bracket. The locking bolt automatically secures the clamping bracket in the closed position.

3. CONTENTS LIABILITY

The manufacturer accepts no liability for damage in the following cases:

- Assembly or connection error
- Damage to the product due to mechanical influences and incorrect supply voltage
- Modifications to the product without the express permission of the manufacturer
- Use for purposes other than those described in the Instructions for Use

Changes or repairs to the device may only be carried out by technicians certified by the manufacturer.

4. TARGET GROUP

The Instructions for Use are intended for all personnel working with the MultiSupport GROUND (MSG) and must be read carefully before use. The MSG shall only be used by specially trained, qualified medical personnel. Only medical professionals (qualified medical staff) may be in charge of the device during transport preparation and during transport itself.

The MSG is a class I medical device in accordance with the Regulation (EU) 2017/745 on medical devices (MDR).

4.1. REQUIREMENTS

Users of the MSG must be familiar with and trained in the use of extracorporeal gas exchange procedures and circulatory therapies under the guidance of a physician.

In addition, users must:

- Have knowledge of the methodology of extracorporeal gas exchange procedures and circulatory therapies as well as the relevant physiology.
- Be familiar with the contents of these Instructions for Use.
- Be fully trained and qualified in the operation of the system, and be able to distinguish between normal and abnormal operation.

5. INTENDED USE

The MSG is intended for intra- and inter-hospital transport of intensive care patients who require extracorporeal gas exchange or extracorporeal circulatory support.

The MSG is intended to fix the following components:

- Xenios console
- DP3 pump drive
- Xenios sensor box
- Patient kits approved for ground based inter-hospital transport:
 - XLung kit 230
 - XLung kit
 - MiniLung kit
 - MiniLung petite kit

The MSG is only approved for ground transport and accepts the patient kit components (oxygenator, DP3 pump head) as well as the technical components (pump drive, sensor box and console) during transport. (According to the standards IEC 60601-1-12, IEC 60068-2-27 / IEC 60068-2-64).

The MSG is intended for use in a hospital environment and for transport in emergency services.

Examples of possible locations:

- Intensive care units
- Catheter laboratories
- Operating rooms
- Ambulance
- Intensive care transport

Note: Detailed information can be found in the Instructions for Use of the console and/or the respective patient kits.

Note: Device category according to IEC 60601-1-12: Fixed in ambulance and transportable within the hospital environment.

Intended use also includes:

- the use according to the specifications
- compliance with the maintenance intervals
- correct storage

WARNING

Danger to life due to console failure.

Do not use the MSG for airborne transport.
Only use the MSG for ground transport.

Therapy failure due to incompatible components.

Incompatible components endanger patient care.
Only consoles with the following symbol on the console and sensor box may be used for inter-hospital transport:

Therapy failure due to incompatible patient kits.

Incompatible patient kits endanger patient care.
Only use the approved patient kits for inter-hospital transport (see WARNING “Danger to life due to incompatible patient kits.” on page 15)

5.1. INDICATIONS

The MSG is generally indicated for use for neonatal, pediatric and adult patients with cardiac and/or respiratory failure.

The MSG is also indicated for use in adult patients in cardiogenic shock or requiring high-risk percutaneous coronary intervention (PCI)

The duration of use depends on the patient kit selected for the therapy required. Please read and comply with the relevant Instructions for Use.

5.2. CONTRAINDICATIONS

No contraindications are known with regard to the use of MSG components.

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

6. SCOPE OF DELIVERY

Number	Part name
1	MSG Holder
1	MSG Base Plate (for inter-hospital transport)
1	MSG Crossbar (for intra-hospital transport)
1	MSG Backup Module

7. PACKAGING AND STORAGE

Use the original packaging for shipping if the unit must be returned to the manufacturer for service or repair.

Only store MSG components in a completely dry condition.

Protect the components from dust and high humidity during storage.

The manufacturer is not liable for damage resulting from incorrect or improper storage.

DANGER

In order to guarantee a constant power supply to the unit, the charge state of the battery packs must be checked before each use. The battery packs must be connected to an external power supply for 6 hours before being put into operation, so that sufficient capacity is guaranteed!

When operating the console, e.g. while transporting the patient, there must be a sufficient number of fully charged battery packs (MSG Backup Module) close at hand. Battery packs not in the Power Supply may only be stored when fully charged.

Battery packs may only be stored when fully charged.

Do not short-circuit a battery pack.

Store the battery packs in a dry and clean environment.

Prevent contact of metal objects with the terminals of the battery packs.

8. OVERVIEW OF THE MSG COMPONENTS

MSG Holder and MSG Backup Module must be available for each transport. Depending on the situation, one of the two following components is added:

- Transport by ambulance: MSG Base Plate
- Transport in the hospital: MSG Crossbar

8.1. MSG HOLDER

The MSG Holder [fig. 1] fixes the pump drive, sensor box and gas exchanger. The MSG Holder is standard equipped with an MSG Pin 7000 as standard.

When using the MiniLung kit or MiniLung petite kit with the MSG Holder, the MSG Pin 2400 is required (optionally available).

When using the iLA active iLA kit IPS or the i-cor puls kit with the MSG Holder, the MSG Bracket is required (optionally available).

CAUTION

Danger of components breaking off.

The MSG Bracket is not approved for ground transport.

The MSG Bracket shall only be used for inpatient use in hospitals including intra-hospital transport.

EN

Pos. in fig. 1:

Components of the MSG Holder

1	MSG Pin 7000 (standard): Only for the XLung kit oxygenators
2	MSG Pin 2400: Only for MiniLung kit or MiniLung petite kit (optionally available)
3	MSG Bracket: Only for iLA active iLA kit IPS or i-cor puls kit (optionally available)
4	U-shaped handle
5	Adjustable clamping handle for the pump drive holder
6	Pump drive holder
7	Sensor box carrier
8	Locking bolt of the sensor box carrier
9	Cable guides
10	Pedestal of the gas exchanger holder
11	Locking button of the pedestal of the gas exchanger holder
12	MSG Handle Fixations

8.2. MSG BASE PLATE

The MSG Base Plate [fig. 2] fixes the console during transport in an ambulance.

Pos. in fig. 2: Components of the MSG Base Plate	
1	Folding handle
2	Threaded hole for MSG Distance Screws
3	MSG Fittings (4x)
4	Lobed knob of the MSG Fitting
5	Counter plate of the MSG Fitting
6	Screw of the MSG Fitting
7	Slot for MSG Fittings
8	Airline rails
9	MSG Distance Screw (4x)
10	Locking bolt of the clamping brackets
11	Clamping bracket
12	Console mount for the base plate

8.3. MSG CROSSBAR

The MSG Crossbar [fig. 3] fixes the console during transport in the hospital on a patient stretcher (e.g. Intensive-Transport-Unit 1 [ITU1]). For this purpose, the MSG Crossbar is attached to two infusion poles positioned at opposite corners at the foot of the patient stretcher.

Pos. in fig. 3: Components of the MSG Crossbar	
1	Screw clamps
2	MSG Distance Screw

Pos. in fig. 3: Components of the MSG Crossbar

3	Release button
4	Console mount

8.4. MSG BACKUP MODULE

The MSG Backup Module (case) [fig. 5] is designed for emergencies. In order to maintain extracorporeal circulation support in an emergency, a backup pump drive and two fully charged battery packs must always be within immediate reach. This is why the MSG Backup Module must always be carried directly next to the console during transport.

The backup pump drive must be inserted into the MSG Backup Module before transport begins. The backup pump drive must be taken from the same console used for transport.

Pos. in fig. 5: Components of the MSG Backup Module

1	MSG Backup Module (case)
2	Recess for a pump drive
3	Four (4) tubing clamps
4	MSG Rain Cover
5	Battery Pack
6	MSG Battery Case with inserted battery pack
7	MSG Instructions for Use
8	MSG Power Cable

8.5. EXTERNAL BATTERY CHARGER

The external battery charger [fig. 4] charges the battery packs in the MSG Backup Module.

Only one battery pack can be charged at a time.

9. APPLICATION

9.1. BASIC INFORMATION ON PATIENT TRANSPORT

The transport of a patient with ongoing ECMO/ECLS treatment should generally be weighed very carefully. The patient must be in a stable condition for transport.

The intra- and inter-hospital transport of a patient under extracorporeal circulation (gas exchange or circulatory support) shall only be carried out by medical personnel who are trained and experienced in the field of ECLS and in handling the console and the components of the MSG.

During the transport of a patient with an extracorporeal system, the vascular accesses (cannulas) in particular, must be fastened and carefully secured in such a way that they can be checked at any time by the accompanying specialist personnel. In addition, the patient kit tubing system must be positioned so that it cannot be accidentally pinched, compressed or kinked. The tubing system must be secured during transport in such a way that there is no pull on the cannulas and that it cannot be damaged.

For transport, the activated alarms, the alarm limits and the activated functions of the console must be checked again.

If the patient is transferred (from the patient's bed to a stretcher, from patient's bed to an examination table, e.g. CT), care must be taken to ensure that the safety of the patient and the treatment is not affected at any time.

Components (e.g. the MSG Holder, the console) must always be positioned or secured during the transfer process in such a way that the patient kit and/or the cannulas are not exposed to any mechanical stress or even danger of pulling and/or dislocation.

After the transfer process has been carried out, the components (e. g. the MSG Holder, the console) must be secured immediately (e. g. attaching and securing the MSG Holder to a standard rail) near the patient.

9.2. SAFETY DURING PATIENT TRANSPORT

DANGER

Danger to life due to failing components.

If a complete MSG Backup Module is not available, the system cannot be safely operated in the event of an emergency or failure of the console, the pump drive or battery packs.

Before using the MSG, insert the appropriate backup pump drive into the MSG Backup Module.

Before use, check the charge status of the battery packs in the console and in the MSG Backup Module and replace them with fully charged battery packs if necessary.

Always keep the MSG Backup Module close to the console during transport.

WARNING

Therapy failure due to damaged cables and tubing.

Damaged or kinked tubing can interrupt therapy and endanger the patient.

The console and the MSG Holder must be positioned around the patient so that the patient cannot touch the components.

Place the components around the patient so that the extracorporeal circulation can be safely maintained.

WARNING

Patient hazard due to inadequate care.

The patient can be endangered due to undersupply (e.g. insufficient circulatory support).

Monitor the patient's vital parameters during transport.

During transport, monitor flow, P1 and P3 closely and set the respective alarm limits taking into consideration the patient's condition and the transport situation.

CAUTION**Danger of heat loss for the patient.**

If there is no active temperature regulation of the patient's blood during transport, the patient must be protected against heat loss.

9.3. MOUNTING THE MSG HOLDER

The MSG Holder accommodates the pump drive, the sensor box and the associated gas exchanger. If transport is considered prior to treatment, the MSG Holder should be used from the start of treatment. This eliminates the need to convert components of the patient kit from another holder to the MSG Holder.

Note: When transferring the patient kit from another holder (e.g. a compact holder) to the MSG Holder, the pump drive and pump head must be repositioned, which inevitably interrupts therapy. Therefore, moving the patient kits from another holder to the MSG Holder should be avoided.

The MSG Holder can be attached to the standard rail of the following devices:

- Xenios Trolley N
- Standard rail of a patient stretcher

WARNING**Inadequate therapy due to incorrectly positioned components.**

If the MSG Holder is positioned above the patient's heart level, the pump must apply a greater suction force, which can impair therapy.

When mounting the MSG Holder to the patient stretcher, make sure that the gas exchanger and the pump drive are positioned at the patient's heart level, preferably slightly below it.

1. Clean the MSG Holder (see section 10).
2. Hang the MSG Holder on the standard rail.

3. Move the adjustable clamping handles of the MSG Handle Fixations upwards [fig. 6].
4. Make sure that the adjustable clamping handles are correctly seated on the standard rail.
 - ✓ The MSG Holder is locked onto the standard rail.

9.4. MOUNTING THE MSG PIN OR MSG BRACKET

Note: MSG Pins and the MSG Bracket are not compatible with other holders (from the manufacturer). MSG Pins and the MSG Bracket are only suitable for the MSG Holder!

The following components fit together:

XLung kit:	MSG Pin 7000 (standard)
MiniLung kit	MSG Pin 2400 (optionally available)
MiniLung petite kit	
iLA active iLA kit	MSG Bracket (optionally available)
iLA active iLA kit IPS	
i-cor puls kit	

Note: Information on the installation/equipment of the patient kits can be found in the respective Instructions for Use of the corresponding patient kit.

9.4.1. FASTENING THE PIN OR BRACKET

1. Press the two locking buttons [fig. 7, pos. 1].
2. Insert the pin or bracket. Position the stud in the recess of the pin or bracket [fig. 7, pos. 2].
3. Release the two locking buttons.
4. Check that the pin or bracket is firmly seated and secured.

9.4.2. REMOVING THE PIN OR BRACKET

1. Firmly hold the pin or bracket.
2. Press the two locking buttons [fig. 7, pos. 1].

3. Remove the pin or bracket.
4. Release the two locking buttons.

9.5. MOUNTING THE SENSOR BOX

1. Pull the locking bolt and swivel the sensor box carrier forward [fig. 8, 1. and 2.].
2. Mount the sensor box on the sensor box carrier [fig. 8, 3.].
3. Pull the locking bolt and push the sensor box carrier backwards until it is completely in contact with its rear side in the holder [fig. 8, 4. and 5.].
4. Release the locking bolt.
 - ✓ The sensor box is locked.

9.6. MOUNTING THE PUMP DRIVE

1. Release the adjustable clamping handle [fig. 10, 1.].
2. Insert the pump drive into the holding fixture [fig. 10, 2.]. The pump drive is in the correct position if the following conditions are met:
 - The motor is located to the left of the sensor box.
 - The nose on the handle of the pump drive is within the holding fixture.
 - The distance between the handle of the pump drive and the U-shaped handle of the MSG Holder is approx. 20 mm [fig. 10].
 - The pump drive **does not** touch the sensor box.
3. Close the adjustable clamping handle [fig. 10, 3.]. The pump drive is correctly seated in the holding fixture when the "nose" of the handle touches the holding fixture of the pump drive completely and is aligned horizontally.

Note: When the pump drive is inserted into the holding fixture for the pump drive, the adjustable clamping handle can move slightly upwards. This is normal and does not represent a malfunction of the adjustable clamping handle.

4. Pass the pump drive cable through the cable guides on the right side of the MSG Holder [fig. 9].
5. Wind up excess cable behind the sensor box. Make sure that the cable is not exposed to tensile stress.
6. Plug the pump drive cable into the sensor box.

9.7. INSTALLING PATIENT KITS ON THE MSG HOLDER

1. Install the patient kit according to the Instructions for Use of the patient kit (placing the gas exchanger on the pin, etc.).
 - ✓ The patient kit is fully prepared for connection to the patient.
2. Position the cables of the integrated pressure sensors and the flow sensor and fasten them using the clips of the MSG Holder [fig. 1, pos. 9]. Make sure that the tubing lines are still movable. The tubing lines must not be obstructed by the cables.

CAUTION

Patient hazards due to damaged cables and tubing.

Damaged or kinked tubing can interrupt therapy and endanger the patient.

Make sure that the tubing between the gas exchanger and the pump head is not kinked. When positioning the gas exchanger, also make sure that the tubing section between the pump head and the gas exchanger is free of kinks.

3. During transport, move the U-shaped handle forward 90° to protect the essential components within the MSG Holder.

9.8. MOUNTING THE MSG COSSBAR

CAUTION

Risk of crushing due to falling parts.

The MSG Crossbar can fall down when attached to infusion poles that cannot carry the load.

Only attach the MSG Crossbar together with the console to patient carriers with infusion poles that can support a weight of at least 13 kg.

The MSG Crossbar is intended exclusively for mounting the console on a patient stretcher for intra-hospital transport (from the intensive care unit to the ambulance and back). For this purpose, the MSG Crossbar is attached to two infusion poles positioned at opposite corners at the foot of the patient stretcher.

1. Loosen the MSG Distance Screw.
2. Position and tighten the screw clamps of the MSG Crossbar by hand on the infusion poles of the patient stretcher [fig. 3, 1]. Make sure that the MSG Crossbar is aligned horizontally.
3. Adjust the distance between the infusion poles with the MSG Distance Screw.
4. Tighten the MSG Distance Screw by hand [fig. 3, 2].
5. Ensure that the MSG Crossbar is securely attached to the infusion poles of the patient stretcher.
 - ✓ The console can be attached to the console mount of the MSG Crossbar.

9.8.1. MOUNTING THE CONSOLE TO THE MSG CROSSBAR

Note: Before starting the transport, make sure that the console is firmly and securely seated on the console mount of the MSG Crossbar.

Note: The display of the console should be oriented in such a way that it can be clearly seen by the physician in charge during transport.

1. Place the console onto the console mount of the MSG Crossbar.
2. Push the console sideways until it snaps into the console mount.
 - ✓ The console is placed and secured onto the console mount of the MSG Crossbar.

9.8.2. REMOVING THE CONSOLE FROM THE MSG CROSSBAR

1. Ensure that the console is held securely by hand.
2. Pull the release button [fig. 3, 3].
3. Push the console sideways and remove it from the console mount of the MSG Crossbar.

9.9. MOUNTING THE MSG BASE PLATE

With the help of the MSG Base Plate, the console is securely fastened in the ambulance during inter-hospital transport. The MSG Base Plate can be mounted as follows:

- on two (2) airline rails [fig. **11 B**] (recommended)
- on one (1) [fig. **11 A**] (optional)

Note: The MSG Base Plate can be mounted on airline rails with both metric and inch dimensions.

9.9.1. MOUNTING THE MSG BASE PLATE ON TWO AIRLINE RAILS (RECOMMENDED)

Note: When mounting the MSG Base Plate on two airline rails, the MSG Distance Screws are not required and can be removed.

1. Insert the screw parts of the MSG Fittings into the airline rail. Note the distance between the holes [fig. **11**, 1. and 2.].
2. Fix the screw parts with the counter plates [fig. **11**, 3.].
3. Place the MSG Base Plate on the MSG Fittings in the airline rails [fig. **11 B**, 4.].
4. Fix the MSG Base Plate manually with the four (4) lobed knobs and tighten hand-tight [fig. **11 B**, 5.].
5. Make sure that the MSG Base Plate is properly fixed.

9.9.2. MOUNTING THE MSG BASE PLATE ON ONE AIRLINE RAIL (OPTIONAL)

1. Insert the screw parts of the MSG Fittings into the airline rails. Note the distance between the holes [fig. **11**, 1. and 2.]. Use two (2) MSG Fittings per side.
2. Fix the screw parts with the counter plates [fig. **11**, 3.].
3. **Only if the airline rail is above floor level of the ambulance:** Insert the MSG Distance Screws into the threaded holes of the MSG Base Plate [fig. **11 A**, 4. and 5.].
4. Place the MSG Base Plate on the MSG Fittings in the airline rails [fig. **11 A**, 6.].
5. Fix the MSG Base Plate manually with the four (4) lobed knobs and tighten hand-tight [fig. **11 A**, 7.].

6. Using the MSG Distance Screws, position the MSG Base Plate so that it rests solidly on the floor of the ambulance.
7. Secure the MSG Distance Screws using the lobed knob.
8. Make sure that the MSG Base Plate is properly fixed.

9.9.3. MOUNTING THE CONSOLE ON THE MSG BASE PLATE

1. Attach the console to the console mount of the MSG Base Plate [fig. 2, 12]. The center knob of the console should point to the middle of the clamping bracket. The console display should point in the direction of the accompanying emergency physician.
2. Press down the clamping bracket of the MSG Base Plate [fig. 2, 11].
 - ✓ The locking bolt automatically secures the clamping bracket when closing [fig. 2, 10].
 - ✓ The console is fixed.
3. Ensure that the console is properly mounted on the MSG Base Plate.

9.9.4. REMOVING THE CONSOLE FROM THE MSG BASE PLATE

1. Ensure that the console is held securely by hand.
2. Pull the locking bolt [fig. 2, 10] and lift the clamping bracket.
3. Pull the console towards the clamping bracket and remove it from the MSG Base Plate.

9.10. CONNECTING THE CONSOLE TO THE AC POWER CONNECTION ON BOARD OF THE AMBULANCE

WARNING

Therapy failure due to system or pump failure.

The system cannot be safely operated without the backup pump drive in the event of a malfunction or failure of the primary pump drive.

Before using the MSG, insert the backup pump drive into its designated recess in the MSG Backup Module.

Always carry the MSG Backup Module with the MSG Holder during transport.

CAUTION**Patient hazards due to console failure.**

Incorrect mains voltage can damage the console resulting in an interruption to therapy and endangering the patient.

The console alternating electrical current (AC) terminal must **not** be connected to the ambulance alternating electrical current (AC) supply.

To connect the console to the ambulance's AC power supply, only the supplied MSG Power Cable may be used (located in the MSG Backup Module).

The console can be powered via the MSG Power Cable by the AC power supply on board the ambulance to ensure a continuous voltage to the console during transport.

1. Plug the MSG Power Cable into the AC terminal of the ambulance and into the direct electrical current (DC) terminal of the console.
2. To unplug the MSG Power Cable, press the button on the DC plug to pull the DC plug out of the console's DC terminal.

9.11. CHARGING THE BATTERY PACKS

DANGER

In order to guarantee a constant power supply to the unit, the charge state of the battery packs must be checked before each use. The battery packs must be connected to an external power supply for 6 hours before being put into operation, so that sufficient capacity is guaranteed!

When operating the console, e.g. while transporting the patient, there must be a sufficient number of fully charged battery packs (MSG Backup Module) close at hand. battery packs not in the Power Supply may only be stored when fully charged.

Battery packs may only be stored when fully charged.

Do not short-circuit a battery pack.

Prevent contact of metal objects with the terminals of the battery packs.

1. Connect the external battery charger to the mains.
2. Connect the battery pack to the external battery charger.
 - ✓ The status light of the external battery charger flashes green rapidly.
 - ✓ The status light of the external battery charger flashes green slowly when the battery pack is fully charged.
3. Press the status button of the battery pack to check its charging status.
4. Disconnect the battery pack from the external battery charger.

9.12. DEPLOYING THE MSG BACKUP MODULE

The MSG Backup Module comprises the replacement components needed to maintain emergency operation during intra- and inter-hospital transport.

9.12.1. BEFORE USE

1.  Press the status button on each battery pack to check the charging status of the battery packs.
2. Only place fully charged battery packs in the MSG Backup Module.
3. Place a suitable backup pump drive in the recess for the pump drive.
4. Check the MSG Backup Module for completeness and refill if necessary (see section 8.4).

9.12.2. AFTER USE

1. Charge the battery packs regularly with the external MSG Battery Charger.
2. Check the MSG Backup Module for completeness and refill if necessary (see section 8.4).

10. CLEANING

WARNING

Danger of electric shock.

To prevent electric shock, disconnect power supplies prior to cleaning the device.

Before each use:

1. Clean the MSG Holder (see below). The cleaning maintains the functionality and prevents infections, contaminations or malfunctions.

If necessary:

1. Clean the components with a mild detergent (e.g. dishwashing liquid).
2. For disinfection, use a medical grade disinfectant (e.g. Bode Chemie GmbH Kohrsolin, Sterilium, Bacilol Plus). For further approved cleaning solutions refer to the console Instructions for Use. Observe the safety instructions of the disinfectant manufacturer.
 - ✓ After cleaning and disinfection, allow the MSG components to dry completely.

11. MAINTENANCE

DANGER

Maintenance work may only be done by qualified maintenance staff who have been certified by the manufacturer. Incorrectly performed maintenance activities may endanger patient lives.

Xenios AG only guarantees the function and safety of the Transport Holder System MSG if the specified maintenance is done on time and in full.

Generally, the Transport Holder System MSG must undergo annual maintenance that is precisely defined by Xenios AG. Regular and correct maintenance and also technical safety checks are requirements for any warranty services and liability on the part of Xenios AG.

Maintenance work must be carried out regularly and properly. The operator must keep records of the maintenance work.

Damaged components must not be used. Damaged components may only be repaired or replaced by a Xenios authorized service technician.

Inspection Interval

Type of inspection	Interval
Visual inspection	12 months
Electrical Safety Test	12 months
Functional test	12 months

Service Interval

Type of service action	Interval
Replacement of battery packs	36 months

12. TROUBLESHOOTING

12.1. ENSURE BACKUP OPERATION

In order to be able to act quickly in alarm situations that require the use of the backup operation, the correct procedure for the backup operation is described here.

1. Clamp tubing lines (clamps in the MSG Backup Module if required).
2. Remove the battery pack with the MSG Battery Case from the MSG Backup Module.
3. Remove backup pump drive from MSG Backup Module.
4. Insert the backup pump drive directly into the connection socket of the battery pack (align the red markings with each other).
 - ✓ After a few seconds you will hear a sequence of beeps.
 - ✓ The pump drive is ready for operation.
5. Press the locking clamps on the defective pump drive and remove the pump head.
6. Pull pump drive's cable out of the sensor box.
7. Remove the defective pump drive from the MSG Holder.
8. Mount the backup pump drive on the MSG Holder.
9. Connect the pump head to the backup pump drive.
10. Press the  key.
 - ✓ The pump drive starts.
11. Set the speed (rpm) to the target value:

Symbol	Press key	Press and hold
	Increase speed by 100 rpm	Increase speed once by 3000 rpm
	Reduce speed by 100 rpm	Reduce speed once by 3000 rpm

12. Remove the clamps from the tubing.
13. Attach the battery pack to the U-shaped handle of the MSG Holder or to the infusion pole of the patient stretcher.
14. Check the efficacy of the extracorporeal circulation in backup operation and verify it on the basis of the patient's vital parameters (e.g. SpO₂).

12.2. ALARMS - MEASURES

This section describes important alarms that can occur during transport and appear on the display of the console. It also describes the measures to be taken to eliminate the cause of an alarm during transport.

Alarm priority: The alarms are prioritized according to the following levels:

- [', ] : "high"
- [', ] : "medium"
- [', ] : "low"

Alarm No.:	Alarm Priority	Alarm message	Measures to be taken during transport
No. 007	•	Technical failure, power supply	<ul style="list-style-type: none"> – Start backup operation (see section 12.1 on page 39). – Monitor the patient's vital parameters closely!
No. 008	•	Technical failure, sensor box	
No. 009	•	Technical failure control panel	
No. 00A	•	Technical failure, pump drive	

Alarm message		Measures to be taken during transport	
Alarm No.:	Alarm priority		
No. 103	•	<p>Failure pump head (Thrombus?):</p> <ul style="list-style-type: none"> – Set SETPOINT SPEED cannot be reached. – Increased resistance in pump head (Thrombus?) 	<p>⚠ Check patient vital parameter (e.g. SpO₂)</p> <ul style="list-style-type: none"> – Check speed (rpm) and blood flow! – If the speed (rpm) and thus the required blood flow (L/min) cannot be reached, start backup operation (see section 12.1 on page 39). – Monitor the patient's vital parameters closely!
No. 109	•	<p>Technical failure, pump drive</p>	<ul style="list-style-type: none"> – Start backup operation (see section 12.1 on page 39). – Monitor the patient's vital parameters closely!

Alarm message		Measures to be taken during transport	
Alarm No.:	Priority	Air bubbles detected – zero flow activated	<ul style="list-style-type: none"> – Check the tubing lines after the gas exchanger for air. 1. If only small air bubbles are visible (diameter of a pin head): <ul style="list-style-type: none"> – Turn the center knob clockwise to adjust the speed according to the previous setting. 2. If a lot of air is visible: <ul style="list-style-type: none"> – Check the pump head and gas exchanger for air. If there is air in the pump head and gas exchanger, the patient kit must be replaced. ⚠ The patient kit cannot be exchanged during transport! <ul style="list-style-type: none"> – Clamp the cannulas using the clamps contained in the MSG Backup Module, using cross-clamping technique. – Adjust ventilation settings to the patient's needs and/or take life support measures. – If possible, flush the cannulas with sterile NaCl 0.9% to avoid blood clots. – Get to a hospital as fast as possible! – Monitor the patient's vital parameters closely.
No. 112	•	Backflow detected – zero flow activated:	<ul style="list-style-type: none"> – Negative flow measured – P1 limiter possibly disabled
No. 114	•	Zero flow is activated manually!	<ul style="list-style-type: none"> – Check the patient's vital parameters. – Check the tubing system and the condition to the patient. – Turn the center knob clockwise to increase the speed until the flow is adjusted to the patient's needs.
No. 117	•		<ul style="list-style-type: none"> – Turn the center knob immediately clockwise to adjust the speed according to the previous setting.

Alarm message		Measures to be taken during transport	
Alarm No.:			
Alarm priority:	•		
No.:	210	Air bubbles detected	<ul style="list-style-type: none"> – Check the tubing lines after the gas exchanger for air. – 1. If only small air bubbles are visible (diameter of a pin head): <ul style="list-style-type: none"> – Turn the center knob clockwise to adjust the speed according to the previous setting. – 2. If a lot of air is visible: <ul style="list-style-type: none"> – Check the pump head and gas exchanger for air. If there is air in the pump head and gas exchanger, the patient kit must be replaced. <p>⚠ The patient kit cannot be exchanged during transport!</p> <ul style="list-style-type: none"> – Clamp the cannulas using the clamps contained in the MSG Backup Module, using cross-clamping technique. – Adjust ventilation settings to the patient's needs and/or take life support measures. – If possible, flush the cannulas with sterile NaCl 0.9% to avoid blood clots. – Get to a hospital as fast as possible! – Monitor the patient's vital parameters closely. <p>Note: Consider cause of air entry.</p>

Alarm message Measures to be taken during transport


Alarm No.:
Priority

No.
214

Pressure P1 too low:


- Below limiting setpoint value


- Check access line for kinks, clots and obstructions.
- Check the access cannula for correct positioning.
- If the cannula is correctly in place and the blood flow is within the normal range, adjust the pressure limit:
- Touch the P1 pressure window:

- Touch the alarm limit (in the pressure limits window):
- Set a suitable limit value with the center knob and confirm by pressing the green check mark.
- Close the window with the button **X**.
- If the patient's arterial blood pressure is low, infusion of volume should be considered.
- If the cannula is displaced:
- Clamp the cannulas using the clamps contained in the MSG Backup Module, using cross-clamping technique.
- Adjust ventilation settings to the patient's needs and/or take life support measures.
- If possible, flush the cannulas with sterile NaCl 0.9% to avoid blood clots.
- Get to a hospital as fast as possible!
- Monitor the patient's vital parameters closely.

Alarm message Measures to be taken during transport	
Alarm No.: 21A	<p>Alarm priority: •</p> <p>Pressure P3 too low:</p> <ul style="list-style-type: none"> – Below limiting setpoint value
	<ul style="list-style-type: none"> – Check the return cannula. – If the cannula is properly in place, check the patient's vital parameters and the alarm limit and adjust if required: – Tap the P3 pressure window:  <ul style="list-style-type: none"> – Touch the lower P3 alarm limit (in the pressure limit window): – Set a suitable limit value with the center knob and confirm by pressing the green check mark. – Close the window with the button X.

Alarm message		Measures to be taken during transport	
Alarm No.:			
Alarm priority:			
No. 221	<ul style="list-style-type: none"> • Blood flow too low: <ul style="list-style-type: none"> – Below limiting setpoint value 	<ul style="list-style-type: none"> – Check access line for kinks, clots and obstructions. – Check the patient's vital parameters. – If necessary, adjust the blood flow. – Turn the center knob immediately clockwise until the flow is adjusted to the patient's needs. – Adjust the lower alarm limit: <ul style="list-style-type: none"> – Tap the flow window:  – Touch the lower alarm limit (in the flow limit window):  – Set a suitable limit value with the center knob and confirm by pressing the green check mark. – Close the window with the button . 	

Alarm message		Measures to be taken during transport	
Alarm No.:	Priority	No. 222	<ul style="list-style-type: none"> • <p>Blood flow too high:</p> <ul style="list-style-type: none"> – Exceeds limiting setpoint value  <ul style="list-style-type: none"> – Check the patient's vital parameters – If necessary, adjust the blood flow. – Turn the center knob counter clockwise until the flow is adapted to the patient's needs. – Adjust the upper alarm limit: <ul style="list-style-type: none"> – Tap the flow window: – Touch the upper alarm limit button (in the flow limit window):  (0.0 limit) – Set a suitable limit value with the center knob and confirm by pressing the green check mark. – Close the window with the button .
Alarm No.:	Priority	No. 225	<ul style="list-style-type: none"> • <p>Failure pump head:</p> <ul style="list-style-type: none"> – Speed deviation detected. There may be a defect in the speed in the speed monitoring. <p>⚠ Check patient vital parameter (e.g. SpO₂)</p> <ul style="list-style-type: none"> – Check speed (rpm) and Blood flow! – If the speed (rpm) is thus the required flow (L/min) cannot be reached, start backup operation (see section 12.1 on page 39). – Monitor the patient's vital parameters closely!

Alarm message Measures to be taken during transport

Alarm No.:	Alarm Priority	Alarm message
No. 225 and No. 103		<p>Failure pump head:</p> <ul style="list-style-type: none"> Speed deviation detected. There may be a defect in the speed monitoring. There may be increased resistance in the pump head. <p>⚠ Extracorporeal circulation can no longer be maintained due to pump head embolism. The patient kit cannot be exchanged during transport!</p> <ul style="list-style-type: none"> Clamp the cannulas using the clamps contained in the MSG Backup Module, using cross-clamping technique. Adjust ventilation settings to the patient's needs and/or take life support measures. If possible, flush the cannulas with sterile NaCl 0.9% to avoid blood clots. Get to a hospital as fast as possible! Monitor the patient's vital parameters closely.
No. 231	●	<p>Pressure P3 too high:</p> <ul style="list-style-type: none"> Set limit value exceeded <p>– Check the return cannula. – If the cannula is properly in place, adjust the alarm limit: – Tap the P3 pressure window:</p>  <p>– Touch the upper P3 alarm limit (in the pressure limit window).</p> <p>– Set a suitable limit value with the center knob and confirm by pressing the green check mark. – Close the window with the button X.</p>

Alarm message Measures to be taken during transport

Alarm No.:	Alarm priority	Alarm message	Measures to be taken during transport
No. 306	•	Left battery, reduced capacity Capacity on start $\leq 85\%$!	– Provide power to the console from the on-board power supply of the ambulance (see section 9.10).
No. 307	•	Left battery capacity $\leq 50\%$:	
No. 308	•	Left battery charge status $\leq 20\%$:	
No. 30C	•	Right battery, reduced capacity Capacity on start $\leq 85\%$!	
No. 30D	•	Right battery capacity $\leq 50\%$:	
No. 30E	•	Right battery charge status $\leq 20\%$:	
No. 309	•	Left battery empty	
No. 30F	•	Right battery empty	

Alarm message		Measures to be taken during transport	
Alarm No.:	Alarm priority		
No. 312	•	<p>Change to battery operation</p>	<ul style="list-style-type: none"> - The system is disconnected from the main power supply. - Reconnect the system to a main power supply as soon as possible.
No. 31E	•	<p>Battery capacity depleted: The console will be shut down in a few minutes</p>	<ul style="list-style-type: none"> - Immediately ensure power is supplied to the console via the on-board power supply of the ambulance (see section 9.10). - If this is not possible, proceed as follows: <ul style="list-style-type: none"> - Clamp the tubing lines (clamps in the MSG Backup Module). - Remove the Battery Pack with the MSG Battery Case from the MSG Backup Module. - Unplug the cable of the pump drive from the sensor box. - Start the backup operation (see section 12.1 on page 39, step 10 and 11). - Remove the clamp from the tubing lines. - Attach the battery pack to the U-shaped handle of the MSG Holder or to the infusion pole of the patient stretcher. - Monitor the patient's vital parameters.
No. 321	•	<p>Failure DC power supply: Protection shutdown</p>	<p>Power can only be supplied via the battery packs. Note: Monitor batteries charging status closely</p>

13. DISPOSAL

The MSG must be disposed of in accordance with the local waste and disposal guidelines. The manufacturer shall not be held liable for violations of local guidelines.

14. TECHNICAL DATA

IP67	Protection class IP67: Complete protection against contact, protection against ingress of dust; protection against ingress of water during temporary immersion
IP33	Protection class IP33: Protection against contact with instruments; protection against obliquely falling water

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please visit this page: www.xenios-ag.com/SVHC.

14.1. MSG HOLDER

Length	542 mm
Width	415 mm
Depth	230 mm
Weight of the component	4.7 kg
Total weight (complete equipment)	9.5 kg
Dimensions for standard rail	10 x 25 mm
Minimum required length of the standard rail	350 mm

14.2. MSG CROSSBAR

Length	740 mm
Width	80 mm
Depth	220 mm
Weight of the component	2.9 kg
Maximum distance between infusion poles	570 mm
Minimum distance between infusion poles	500 mm
Total weight (with mounted console)	13 kg

14.3. MSG BASE PLATE

Length	560 mm
Width	145 mm
Depth	390 mm
Weight	7.0 kg
Maximum distance between the airline rails	505 mm
Minimum distance between the airline rails	375 mm
Total weight (with mounted console)	17.0 kg

14.4. MSG BACKUP MODULE

Length	478 mm
Width	390 mm
Depth	194 mm
Weight	4.0 kg
Total weight (complete equipment)	7.5 kg

14.5. MSG AMBIENT CONDITIONS

Operating temperature (standard temp. conditions)	+18 °C to +40 °C
Operating temperature (20 min. extended temp. conditions)*	-20 °C to +50 °C
Relative humidity	20 % to 90 %
Atmospheric pressure	700 hPa to 1060 hPa
Patient population	The MSG is intended for intra- and inter-hospital transport of intensive care patients (newborns to adults) who require extracorporeal gas exchange and/or circulatory support therapy.

*The MultiSupport Ground including console can be used for 20 minutes in a temperature range of -20 °C to +50 °C. This allows patient transport also in cold or hot days since 20 minutes are sufficient to bring the patient from the hospital to the mobile intensive care unit

15. CHECKLIST FOR TRANSPORT

What is required for transport?		
Check the console:	Yes	No
– System status and settings (flow, speed, sweep gas flow, alarm limits) checked		
– Both battery packs fully charged		
MSG Holder:	Yes	No
– Pump drive properly mounted on MSG Holder		
– Sensor box properly mounted on MSG Holder		
– Gas exchanger properly mounted on MSG Pin		
– Cable securely fixed with clips		
– MSG Handle Fixations functional		
MSG Crossbar:	Yes	No
– MSG Crossbar available		
– Screws of the screw clamps function perfectly		
– MSG Distance Screw present and functioning correctly		
– Console mount intact		
MSG Base Plate:	Yes	No
– Clamping bracket and locking bolt function		
– Four (4) MSG Fittings available		
– Mounting on airline rails clarified (see Instructions for Use, section 9.9.1/9.9.2)		
– Four (4) MSG Distance Screws (for mounting the MSG Base Plate on only one [1] airline rail)		
MSG Backup Module:	Yes	No
– Two (2) battery packs fully charged and available		
– A battery pack already inserted in the MSG Battery Case		
– Backup pump drive available		
– Four (4) tubing clamps available		
– MSG Power Cable / AC/DC-Adapter available		
– MSG Rain Cover available		
– MSG Instructions for Use available		

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