# Handling Devices for PD Accessories

## INSTRUCTIONS FOR USE

Please read the following instructions carefully.

## GENERAL INFORMATION

### General description of the product

Handling devices for PD accessories hold components or products to facilitate peritoneal dialysis (PD) treatment. Products are not sterilized.

Handling devices for PD accessories comprise the following products:

- *stay•safe*® Organizer (called "Organizer" in the following) 2842561N
- *stay•safe*® Holder for Organizer (called "Holder" in the following) 2842571N
- Clip for *stay•safe*® Organizer (called "Clip" in the following) F00006845

## COMPOSITION

Product Name	Material
<i>stay•safe</i> ® Organizer	Body: polyphenylsulfone (PPSU) / brass Suction cups: polyvinyl chloride (PVC) / galvanized steel
<i>stay•safe</i> ® Holder for Organizer	Aluminium / steel / polyamid (PA)
Clip for <i>stay•safe</i> ® Organizer	Acrylonitrile butadiene styrene (ABS)

## INTENDED PURPOSE

Handling devices for PD accessories are intended for facilitating handling of other PD products during continuous ambulatory peritoneal dialysis (CAPD) or automated peritoneal dialysis (APD) treatment in combination with CAPD systems, APD systems and PD-Paed Plus from Fresenius Medical Care.

### Medical indication

The handling devices for PD accessories are indicated to support the treatment of kidney failure requiring CAPD or APD treatment.

### Intended patient population

Adult and pediatric patients suffering from kidney failure.

### Intended user / environment

Handling devices for PD accessories must only be used by individuals with the appropriate training, knowledge and experience and who are certified to have been trained. Handling devices for PD accessories have been specified by the manufacturer for application in rooms suitable for PD located in professional health care facilities, or for home environment.

## SIDE EFFECTS

Side effects are not applicable for Handling devices for PD accessories under their intended purpose.

### Reporting of serious incidents

If any serious incident occurs in relation to the device, including those not listed in this leaflet, the treating physician shall be informed immediately. Within the EU the user must report any serious incident that has occurred in relation to the device to the manufacturer according to labelling (📄) and the competent authority of the EU Member State in which the user is established.

A serious incident can be any incident that directly or indirectly leads to the death of a patient, user or other person; to the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health threat.

## CONTRAINDICATIONS

There are no contraindications known for Handling devices for PD accessories under their intended purpose.

## PERFORMANCE CHARACTERISTICS AND CLINICAL BENEFITS

The Organizer guides the patient during the connection and disconnection procedure in CAPD and APD treatment to improve convenience and to contribute to reduced touch contamination risk. The Holder allows the Organizer to be fixed onto the *sleep•safe*® porter or infusion pole (optional).

Only for APD and PD-Paed Plus treatment, the Clip is inserted into the Organizer for fixation of patient connector.

## METHOD OF ADMINISTRATION

### Handling Instructions

#### Before use

In order to prepare the CAPD systems, APD systems and PD-Paed Plus from Fresenius Medical Care, clean the Organizer and the Clip with a wipe using an appropriate disinfectant.

## Handling

### *stay•safe*® Organizer

Fix the Organizer either on a table by means of its suction cups or fix it on the Holder.


The Organizer presents three areas where components can be inserted during preparation or completion of treatment: The middle of the Organizer can hold a DISC from a CAPD system or the Clip for APD treatments. The two holders on left and right side of Organizer front can hold the *stay•safe*® Catheter Extension and *stay•safe*® Disinfection Cap or *stay•safe*® Closure Cap. For the procedure of preparation and completion of treatment please refer to the instructions from the training or to the instructions for use (IFU) of CAPD systems, APD systems or PD-Paed Plus from Fresenius Medical Care.

### *stay•safe*® Holder for Organizer

Push the tongue of the Holder into the groove in the back of Organizer until stop. Afterwards fix the assembly by the screw with black knob.

### Clip for *stay•safe*® Organizer

In APD and PD-Paed Plus treatment the Clip is part of the Organizer. Only both together enable the preparation and completion of treatment during APD treatment. Insert the Clip into the middle of the Organizer body as shown in the following figure:

Now the patient connector of the APD tubing system or PD-Paed Plus of Fresenius Medical Care can be inserted into the Clip during preparation or completion of treatment. The two holders on the left and right side of the Organizer fix *stay•safe*® Catheter Extension and *stay•safe*® Disinfection Cap or *stay•safe*® Closure Cap during preparation or completion of treatment. For cleaning and disinfection purposes, remove the Clip from the Organizer by pushing the Clip out from the Organizer's back side.

### Intended application time

During the respective PD treatment, the handling devices are used repeatedly during the day, depending on the prescription of the treating physician of the patient (continuous treatment):

– For CAPD treatments, the intended application time is 30 – 60 min

– For APD treatments, the intended application time is 7 – 10 h

– For treatments with PD-Paed Plus, the maximum intended application time is 24 h

### Termination

For the procedure of preparation and completion of treatment, please refer to the instructions from the training or to the IFU of CAPD systems, APD systems and PD-Paed Plus from Fresenius Medical Care.

## WARNINGS AND PRECAUTIONS

### Warnings

The application of wrong disinfectants can lead to damages on products. Fixation of *stay•safe*® Catheter extension and *stay•safe*® Disinfection Cap or *stay•safe*® Closure Cap can be impaired by the damaged two holders on left and right side of the Organizer.

Before start of treatment, make sure that there are no visible damages on Handling Devices for PD Accessories. Do not perform a treatment with defective device. If there are visible damages patient or user needs to inform responsible organization or service support.

### Particular notes on materials and substances

There is no hypersensitivity to the handling devices for PD accessories or any of the materials known.

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please use this page: [www.freseniusmedicalcare.com/en/svhc](http://www.freseniusmedicalcare.com/en/svhc).

## SPECIAL PRECAUTIONS FOR STORAGE

Keep dry.

Keep away from sharp objects such as knives and scissors.

## DISPOSAL

The accredited organization is liable for the safe disposal of the product. Any unused product or waste material shall be disposed of in accordance with local requirements and in accordance with local requirements for potentially contaminated materials.

## SYMBOLS USED ON LABELS

	Catalogue number
	Batch code
	Use-by date
	Medical Device
	Unique Device Identifier
	Caution
	Manufacturer
	Date of manufacture
	Consult instructions for use
	Units
	Patient information website

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**Fresenius Medical Care AG**