# multiFiltratePRO paed HD

<b>en</b>	INSTRUCTIONS FOR USE multiFiltratePRO paed HD.....	3	<b>de</b>	GEBRAUCHSANWEISUNG multiFiltratePRO paed HD.....	4
<b>ar</b>	تعليمات الاستخدام multiFiltratePRO paed HD جهاز .....	5	<b>bg</b>	ИНСТРУКЦИИ ЗА УПОТРЕБА multiFiltratePRO paed HD.....	6
<b>bs</b>	UPUTSTVA ZA KORIŠTENJE multiFiltratePRO paed HD.....	7	<b>cn</b>	使用说明书 multiFiltratePRO paed HD.....	8
<b>cs</b>	NÁVOD K POUŽITÍ multiFiltratePRO paed HD.....	9	<b>da</b>	BRUGSANVISNING multiFiltratePRO paed HD.....	10
<b>el</b>	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ multiFiltratePRO paed HD.....	11	<b>es</b>	MANUAL DE USUARIO multiFiltratePRO paed HD.....	12
<b>et</b>	KASUTUSJUHEND multiFiltratePRO paed HD.....	13	<b>fi</b>	KÄYTTÖOHJEET multiFiltratePRO paed HD.....	14
<b>fr</b>	MODE D'EMPLOI multiFiltratePRO paed HD.....	15	<b>hr</b>	UPUTE ZA UPOTREBU multiFiltratePRO paed HD.....	16
<b>hu</b>	HASZNÁLATI ÚTMUTATÓ multiFiltratePRO paed HD.....	17	<b>it</b>	ISTRUZIONI PER L'USO multiFiltratePRO paed HD.....	18
<b>kk</b>	қолдану жөніндегі нұсқаулық multiFiltratePRO paed HD магистральдар жиынтығы.....	19	<b>ko</b>	사용 설명서 multiFiltratePRO paed HD.....	20
<b>lt</b>	NAUDOJIMO INSTRUKCIJOS multiFiltratePRO paed HD.....	21	<b>lv</b>	LIETOŠANAS INSTRUKCIJA multiFiltratePRO paed HD.....	22
<b>nl</b>	GEBRUIKSAANWIJZING multiFiltratePRO paed HD.....	23	<b>no</b>	BRUKSANVISNING multiFiltratePRO paed HD.....	24
<b>pl</b>	INSTRUKCJA UŻYCIA multiFiltratePRO paed HD.....	25	<b>pt</b>	INSTRUÇÕES DE UTILIZAÇÃO multiFiltratePRO paed HD.....	26
<b>ro</b>	INSTRUCȚIUNI DE UTILIZARE multiFiltratePRO paed HD.....	27	<b>ru</b>	ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ Комплект магистралей multiFiltratePRO paed HD.....	28
<b>sk</b>	NÁVOD NA POUŽITIE multiFiltratePRO paed HD.....	29	<b>sl</b>	NAVODILA ZA UPORABO multiFiltratePRO paed HD.....	30
<b>sr</b>	UPUTSTVO ZA UPOTREBU multiFiltratePRO paed HD.....	31	<b>sv</b>	BRUKSANVISNING multiFiltratePRO paed HD.....	32
<b>tr</b>	KULLANIM TALİMATLARI multiFiltratePRO paed HD.....	33	<b>tw</b>	使用說明 multiFiltratePRO paed HD.....	34

**FRESENIUS  
MEDICAL CARE**

CE 0123

**Fresenius Medical Care AG**

Else-Kröner-Str. 1, 61352 Bad Homburg  
GERMANY ☎ +49 6172 609-0

[www.freseniusmedicalcare.com/en/product-information](http://www.freseniusmedicalcare.com/en/product-information)

8258447/C (07/2023)

# multiFiltratePRO paed HD

## INSTRUCTIONS FOR USE

Please read the following instructions carefully.

### GENERAL INFORMATION

#### General description of the product

The multiFiltratePRO paed HD is intended to be used and is compatible with an extracorporeal blood purification system multiFiltratePRO device.

The multiFiltratePRO paed HD is intended for single use.

The multiFiltratePRO paed HD is designed to provide extracorporeal blood treatment to patients suffering from renal insufficiency.

The multiFiltratePRO paed HD is composed of the arterial line / access line (RED line), venous line / return line (BLUE line), filter line (YELLOW line), dialysate line (GREEN line) and heparin (as anticoagulant) line.

The multiFiltratePRO paed HD is pyrogen-free.

#### Sterilization

The multiFiltratePRO paed HD is ETO sterilised.

### COMPOSITION

Tubes: medical grade soft PVC.

Connectors and other components: Acrylonitrile butadiene styrene (ABS), Polyamide (PA), Polybutylene Terephthalate (PBT), Polycarbonate (PC), Polyethylene (PE), Polyethylene Terephthalate (PET), Polypropylene (PP), Polytetrafluoroethylene (PTFE), Polyvinylchloride (PVC), Silicone, Isoprene Rubber (IR), Thermoplastic Elastomer (TPE), Polyethylene High-Density (PEHD), Methyl Methacrylate Acrylonitrile Butadiene Styrene (MABS), Acrylic Hydrophobic Polyethyl Methacrylate (PMMA), Polyethylene Low-Density (PELD).

### INTENDED PURPOSE AND RELATED DEFINITIONS

#### Intended purpose

Channelling of blood and fluid in an extracorporeal treatment.

#### Medical indication

Acute kidney injury/disease and volume overload requiring CV-VHD.

#### Intended patient population

The multiFiltratePRO paed HD has been specified by the manufacturer for the purpose of treating paediatric patients with a body weight between 8 and 40 kg, under consideration of the specified technical data of the device and the single-use items used (delivery rates, prime volumes).

#### Intended user group and intended environment

The disposable must only be used by individuals with the appropriate training, knowledge and experience of the proper operation and handling and for whom proof of instruction can be shown.

Operations in rooms suitable for dialysis located in professional healthcare facilities.

Normative and local regulations must be observed.

### SIDE EFFECTS

Occasional occurrence of the following side effects is reported in current literature:

- Hypotension
- Blood loss (e.g. hypovolemia and hypovolemic shock, hypotension, anaemia, cardiac arrhythmia, cardiac arrest)
- Coagulation disorders (e.g. thrombocytopenia)
- Embolism (air-, thrombo-)
- Bacterial infection
- Hypothermia
- Haemolysis
- Hypersensitivity or hypersensitivity-like reactions (e.g. dyspnoea, cardiopulmonary arrest, hypotension, urticaria, flushing, itching, abdominal pain, nausea, convulsions and unconsciousness).

Additional side effects might be specific to other devices and drugs used in the therapy.

Refer to Instructions for Use the multiFiltratePRO device for more information on treatment-related side effects.

#### Reporting of serious incidents

If any serious incident occurs in relation to the device, including those not listed in this leaflet, the treating physician shall be informed immediately. Within the EU the user must report any serious incident that has occurred in relation to the device to the manufacturer according to labelling (📄) and the competent authority of the EU Member State in which the treatment is performed.

A serious incident can be any incident that directly or indirectly leads to the death of a patient, user or other person; to the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health treat.

### CONTRAINDICATIONS

#### Product related contraindications

Hypersensitivity to any of the material of the tubings as listed under composition or to ETO residues.

In severe cases the treatment must be discontinued, the blood should not be returned to the patient, and appropriate emergency medical treatment should be initiated.

#### Therapy related contraindication

Refer to the Instructions for use multiFiltratePRO device Instructions for Use for more information of general contraindications for extracorporeal blood purification.

### PERFORMANCE CHARACTERISTICS

Refer to the Instructions for Use multiFiltratePRO device for essential performance parameters.

### METHOD OF ADMINISTRATION

#### Handling Instructions

The multiFiltratePRO paed HD tubing can be used only in combination with the specific paediatric treatment mode.

Refer to Instructions for Use the multiFiltratePRO device of the regarding how to handle the disposable during set-up (e.g.: selection of the treatment type), priming, treatment and reinfusion. Also, the Instructions for Use of other disposables such as the haemofilter must be taken into account.

#### Preparation

• Refer to Instructions for Use the multiFiltratePRO device for details on any preparatory handling of the disposable before it is ready for use or during its use and information to verify if the disposable is properly installed.

- The multiFiltratePRO paed HD is intended to be used and compatible with multiFiltratePRO device (refer to label) and shall be used only after appropriate instruction or training.
- The infusion administration lines should be clamped, except when needed.
- Unpack and connect the disposable aseptically without touching open connectors.
- Colour codes should be followed and used in line with the corresponding markings on the device
- Handle pumps and pump clips in accordance to the Instructions for Use of the device.
- Tighten all closure caps and ensure that all connectors and caps are secured.
- Fill and rinse the disposable in accordance with the Instructions for Use of the device or supplemental Instructions for Use for additional options, and training, if applicable.
- The Fresenius Medical Care disposable is designed to withstand the maximum and minimum manufacturer's recommended pressures and flow rates generated in use with the respective multiFiltratePRO device. All the other relevant specifications of the device apply.

#### Treatment

- In cases where the extracorporeal volume of the circuit consists of a significant amount of the circulating blood volume of the patient, as well as in anaemic or hemodynamically unstable patients, priming of the circuit with blood or blood substitution fluid may be considered.
- To ensure correct functionality of the bubble catchers, fill to about 1 cm below the top. This avoids fluid contacting the hydrophobic filter as well as air entering the extracorporeal circuit. In the event that fluid reaches and comes into contact with the hydrophobic filter within the Transducer Protector (TP) the disposable must be replaced.
- The disposable has an additional connection intended to be connected to the pressure measurement system. This can be used to connect a new pressure measurement line (available as an accessory in the Fresenius Medical Care product range).
- Do not push the fluid the pressure measurement line by using a syringe. This might damage the membrane in the TP and thus lead to contamination. In the event that fluid might have passed through the TP, check the device for contamination after treatment is completed. If contamination has occurred, the device must be taken out of service and disinfected according to the manufacturer's recommendations before further use.
- After manipulation of lines or use of components during treatment, check and, if necessary, restore the correct position of the lines and components.
- Disinfect the corresponding access sites without protective cap prior to connection with other products.

#### Reinfusion

- Refer to Instructions for Use the multiFiltratePRO device for the termination of the treatment and also to the "Disposal" Section of the present Instructions for Use.
- Close all clamps on the multiFiltratePRO Paed HD before removing the disposable to reduce the risk of fluid leakage.

### WARNINGS AND PRECAUTIONS

#### Warnings

- The disposable must only be used by individuals with the appropriate training, knowledge, and experience of the proper operation and handling and for whom proof of instruction can be shown.
- The multiFiltratePRO Paed HD is intended for single use. The correct function of all interfaces is ensured only for single use. Re-use may be hazardous to both patient and operator (e.g. impaired performances, contamination).
- Do not use after use-by-date (refer to label).
- If the carton is damaged, check the products contained carefully. Do not use if the sterile package is damaged, if the protective or closure caps are not in place, or if there is any visible damage on the finished products (e.g. kinked lines).
- Make sure all lines and chambers are correctly inserted into the respective holders. Avoid kinking or occluding the disposable in order to avoid mechanical and chemical damage to cellular blood constituents.
- For hygienic and functional reasons, it is recommended that the selected multiFiltratePRO Paed HD is only inserted immediately prior to the treatment time, thereby keeping the preparation and circulation times nearer to the beginning of the treatment and in compliance with applicable guidelines.
- Connect the multiFiltratePRO Paed HD aseptically without touching open connectors.
- Disinfect with 70% alcohol access sites without protective cap prior to the connection with other products to reduce the risk of infection and let it dry before the connection.
- The compatibility of disinfectants (other than those recommended) with access sites shall be determined prior to clinical use.
- Excessive negative pressure may cause partial collapse of the blood pump segment resulting in an actual blood flow substantially less than indicated on the device.
- Minimum temperature of use of the disposable is 18° C (64° F).
- Do not use needles with a beveled tip and a diameter larger than 20 gauge (outer needle diameter to be 0.9 mm or less) to puncture injection sites (if present).
- To ensure a secure connection between patient access and bloodline, hold and screw the colored (blue, red) coupling nut on the bloodline only. Do not apply the screwing torque to the inner part of the connector. After the connection, check that the components are firmly screwed together.
- Make sure that the components/tubes in direct interface with the device are fitting properly and ensure that all caps and all screw-lock joints are secure and properly tightened (particularly those of the patient connections, the dialyser connections, the device and the solution bags connections).
- Inspect the disposable for kinking and leaks during the priming and treatment phase, taking corrective measures (e.g. tightening Luer-Lock connection) or exchanging it as necessary.
- To avoid air embolism, make sure that the return patient line is correctly inserted into the machine's air bubble detector.
- Careful attention must be given to the installation and blood level of the bubble catchers. Make sure that due to pressure fluctuation neither the blood contacts the hydrophobic filter nor air enters the extracorporeal circuit.
- Cleansing solutions and disinfectants may damage materials employed for the multiFiltratePRO paed HD. Safety and performance of use can no longer be guaranteed, and the manufacturer assumes no liability.

- The plastics used can be incompatible with drugs or disinfectants (e.g. connectors made of polycarbonate can develop cracks when in contact with aqueous solutions with pH > 10).
- If nutritional solutions are administered into the bloodline, wetting of the Luer lock connection with lipidic fluids can weaken the properties of the plastic material used. Make sure that during the connection, the infusion line close to the Luer lock connection site remains completely free of nutrition solution.

#### Operating time

The maximum application time is 72 hours. The disposable must be replaced after the maximum usage time or respectively after a specific blood volume pumped as indicated on the primary packaging.

#### Particular notes on materials and substances

CMR substances and endocrine-disrupting substances  
For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please use this page: [www.freseniusmedicalcare.com/en/svhc](http://www.freseniusmedicalcare.com/en/svhc)

#### Special precautions for storage

Follow the indication of the product label. Protect from moisture, freezing and excessive heat.

#### Disposal

Ensure safe disposal of any unused product or waste material in accordance with local regulations.

Materials that have been in contact with blood or other material of human origin may be infectious. Dispose of such materials by taking the necessary precautionary measures and in accordance with local regulations for (potentially) infectious materials.

### SYMBOLS USED ON LABELS

	Medical Device
	Unique Device Identifier
	Patient information website
	Latex-free
	Sterilized using ethylene oxide
	Exchange the bloodline after maximum usage time or pumped blood volume as indicated on the primary package
	Pump segment diameter / length
	Blood Priming Volume
	Units
	Temperature limit
	Do not use if package is damaged
	Sterile barrier system

### DATE OF REVISION OF TEXT

11/2021

For the electronic version of instruction for use (e-IFU), please use this page:

[www.freseniusmedicalcare.com/en/product-information](http://www.freseniusmedicalcare.com/en/product-information).

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