multiFiltratePRO

en	INSTRUCTIONS FOR USE multiFiltratePRO.....	3	de	GEBRAUCHSANWEISUNG multiFiltratePRO.....	4
ar	تعليمات الاستخدام multiFiltratePRO جهاز	5	bg	ИНСТРУКЦИИ ЗА УПОТРЕБА multiFiltratePRO.....	6
bs	UPUTSTVA ZA KORIŠTENJE multiFiltratePRO.....	7	cn	使用说明书 multiFiltratePRO.....	8
cs	NÁVOD K POUŽITÍ multiFiltratePRO.....	9	da	BRUGSANVISNING multiFiltratePRO.....	10
el	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ multiFiltratePRO.....	11	es	MANUAL DE USUARIO multiFiltratePRO.....	12
et	KASUTUSJUHEND multiFiltratePRO.....	13	fi	KÄYTTÖOHJEET multiFiltratePRO.....	14
fr	MODE D'EMPLOI multiFiltratePRO.....	15	hr	UPUTE ZA UPOTREBU multiFiltratePRO.....	16
hu	HASZNÁLATI ÚTMUTATÓ multiFiltratePRO.....	17	it	ISTRUZIONI PER L'USO multiFiltratePRO.....	18
kk	қолдану жөніндегі нұсқаулық multiFiltratePRO магистральдар жиынтығы.....	19	ko	사용 설명서 multiFiltratePRO.....	20
lt	NAUDOJIMO INSTRUKCIJOS multiFiltratePRO.....	21	lv	LIETOŠANAS INSTRUKCIJA multiFiltratePRO.....	22
nl	GEBRUIKSAANWIJZING multiFiltratePRO.....	23	no	BRUKSANVISNING multiFiltratePRO.....	24
pl	INSTRUKCJA UŻYCIA multiFiltratePRO.....	25	pt	INSTRUÇÕES DE UTILIZAÇÃO multiFiltratePRO.....	26
ro	INSTRUCȚIUNI DE UTILIZARE multiFiltratePRO.....	27	ru	ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ Комплект магистралей multiFiltratePRO.....	28
sk	NÁVOD NA POUŽITIE multiFiltratePRO.....	29	sl	NAVODILA ZA UPORABO multiFiltratePRO.....	30
sr	UPUTSTVO ZA UPOTREBU multiFiltratePRO.....	31	sv	BRUKSANVISNING multiFiltratePRO.....	32
tr	KULLANIM TALİMATLARI multiFiltratePRO.....	33	tw	使用說明 multiFiltratePRO.....	34

**FRESENIUS
MEDICAL CARE**

CE 0123

Fresenius Medical Care AG
Else-Kröner-Str. 1, 61352 Bad Homburg
GERMANY ☎ +49 6172 609-0

www.freseniusmedicalcare.com/en/product-information

8258481/B (07/2023)

multiFiltratePRO

INSTRUCTIONS FOR USE

Please read the following instructions carefully.

GENERAL INFORMATION

General description of the product

The multiFiltratePRO tubings are intended to be used and compatible with an extracorporeal blood purification system multiFiltratePRO device.

The multiFiltratePRO tubings are intended for single use.

The multiFiltratePRO tubings are designed to provide extracorporeal blood treatment of patients suffering from renal insufficiency. multiFiltratePRO tubings are composed of the arterial line / access line (RED line), venous line / return line (BLUE line), filtrate line (YELLOW line), dialysate line (GREEN line), substitute line / plasma line (WHITE line) and citrate line (as anticoagulant), calcium line and heparin line (as anticoagulant).

The variants of multiFiltratePRO tubings are:

- multiFiltratePRO Ci-Ca HDF
- multiFiltratePRO Ci-Ca HD
- multiFiltratePRO HDF
- multiFiltratePRO TPE
- multiFiltratePRO SecuCas Ci-Ca HD
- multiFiltratePRO SecuCas Ci-Ca HDF

multiFiltratePRO tubings can be used in combination with ILA active Medical Device.

The multiFiltratePRO tubings are pyrogen free.

Sterilization

The multiFiltratePRO tubing are ethylene oxide sterilized (EO).

COMPOSITION

Tubes: medical-grade soft PVC.

Connectors and other components: Acrylonitrile butadiene styrene (ABS), Polybutylene Terephthalate (PBT), Polycarbonate (PC), Polyamide (PA), Polyethylene (PE), Polyethylene Terephthalate Glycol-modified (PETG), Polypropylene (PP), Polytetrafluoroethylene (PTFE), Polyvinylchloride (PVC), Silicone, Isoprene Rubber (IR), Thermoplastic Elastomer (TPE), Polyethylene High-Density (PEHD), Methyl Methacrylate Acrylonitrile Butadiene Styrene (MABS), Hydrophobic Acrylic—Polymethyl Methacrylate—(PMMA), Polyethylene Low-Density (PELD), Polyethylene Terephthalate (PET).

INTENDED PURPOSE AND RELATED DEFINITIONS

Intended purpose

Channeling of blood and fluid in an extracorporeal treatment.

Medical indication

- Acute kidney injury requiring continuous renal replacement therapy
- Volume overload requiring continuous renal replacement therapy
- Intoxications requiring continuous renal replacement therapy
- Conditions requiring the exchange of blood plasma via TPE
- Conditions requiring CRRT combined with hemoperfusion to eliminate further pathogens from blood
- Conditions requiring CRRT combined with ECCO2R for additional CO2 removal

Intended patient population

The devices have been specified by the manufacturer for the purpose of treating patients with a body weight of 40 kg and more, irrespective of their age, under consideration of the specified technical data of the device and the single-use items used (e.g. delivery rates, fill volumes).

Intended user group and intended environment

The disposable must only be used by individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown.

Operate in rooms suitable for dialysis located in professional health-care facilities.

Normative and local regulations must be observed.

SIDE EFFECTS

Occasional occurrence of the following side effects is reported in current literature:

- Hypotension
- Blood loss (e.g. hypovolemia, and hypovolaemic shock, hypotension, anaemia, cardiac arrhythmia, cardiac arrest)
- Coagulation disorders (e.g. thrombocytopenia)
- Embolism (air-, thrombo-)
- Bacterial infection
- Hypothermia
- Haemolysis
- Hypersensitivity or hypersensitivity-like reactions (e.g. dyspnoea, cardiopulmonary arrest, hypotension, urticaria, flushing, itching, abdominal pain, nausea, convulsions and unconsciousness).

Additional side effects might be specific to other devices and drugs used in the therapy.

Refer to the Instructions for Use of the multiFiltratePRO device for more information on treatment related side effects.

Reporting of serious incidents

If any serious incident occurs in relation to the device, including those not listed in this leaflet, the treating physician shall be informed immediately. Within the EU the user must report any serious incident that has occurred in relation to the device to the manufacturer according to labelling (📄) and the competent authority of the EU Member State in which the treatment is performed.

A serious incident can be any incident that directly or indirectly leads to the death of a patient, user or other person; to the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health threat.

CONTRAINDICATIONS

Product related contraindications

Hypersensitivity to any of the material of the tubings as listed under composition or to ethylene oxide residues.

In severe cases the treatment must be discontinued, the blood should not be returned to the patient, and appropriate emergency medical treatment should be initiated.

Therapy related contraindication

Refer to the Instructions for Use of the multiFiltratePRO device for more information of general contraindications for extracorporeal blood purification.

PERFORMANCE CHARACTERISTICS

Refer to the Instructions for Use of the multiFiltratePRO device for essential performance parameters.

METHOD OF ADMINISTRATION

Handling Instructions

Refer to the Instructions for Use of the multiFiltratePRO device of the regarding how to handle the disposable during set-up (e.g. selection of treatment type), priming, treatment and reinfusion. Also, the Instructions for Use of other disposables such as the haemofilter must be taken into account.

Preparation

• Refer to the Instructions for Use of the multiFiltratePRO device for details on any preparatory handling of the disposable before it is ready for use or during its use and information to verify if the disposable is properly installed.

• The multiFiltratePRO tubings are intended to be used and compatible with multiFiltratePRO device (refer to label) and shall be used only after appropriate instruction or training.

• The infusion administration lines should be clamped, except when needed.

• Unpack and connect the multiFiltratePRO tubings aseptically without touching open connectors.

• Colour codes should be followed and used in line with the corresponding markings on the device.

• Handle pumps and pump clips in accordance to the Instructions for Use of the device.

• Tighten all closure caps and ensure that all connectors and caps are secured.

• Fill and rinse the disposable in accordance with the Instructions for Use of the device or supplemental Instructions for Use for additional options, and training, if applicable.

• Fresenius Medical Care disposable is designed to withstand the maximum and minimum manufacturer's recommended pressures and flow rates generated in use with the respective multiFiltratePRO device. All the other relevant specifications for the device apply.

Treatment

• In cases where the extracorporeal volume of the circuit consists of a significant amount of the patient's circulating blood volume, as well as in anaemic and/or haemodynamically unstable patients, priming the circuit with blood or blood substitution fluid may be considered.

• To ensure correct functionality of the bubble catchers, fill to about 1 cm below the top. This avoids fluid contacting the hydrophobic filter as well as air entering the extracorporeal circuit. In the event that fluid reaches and comes into contact with the hydrophobic filter within the Transducer Protector (TP) the disposable must be replaced.

• The disposable has an additional connection intended to be connected to the pressure measurement system. This can be used to connect a new pressure measurement line (available as an accessory in the Fresenius Medical Care product range).

• Do not push back the fluid in the pressure measurement line by using a syringe. This might damage the membrane in the TP and thus lead to contamination. In the event that fluid might have passed through the TP, check the device for contamination after treatment is completed. If contamination has occurred, the device must be taken out of service and disinfected according to the manufacturer's recommendations before further use.

• After the manipulation of lines or the usage of components during the treatment check and, if necessary, restore the correct position of the lines and components.

• Disinfect the corresponding access sites without protection cap prior to connection with other products.

Reinfusion

• Refer to the Instructions for Use of the multiFiltratePRO device for the termination of the treatment and also to the "Disposal" Section of the present Instructions for Use.

• Close all clamps on the multiFiltratePRO tubings before removing the disposable to reduce the risk of fluid leakage.

WARNINGS AND PRECAUTIONS

Warnings

• The disposable must only be used by individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown.

• The multiFiltratePRO tubings are intended for single use only. Re-use may be hazardous to both patient and operator (e.g. impaired performances, contamination).

• Do not use after use-by date (refer to label).

• If the carton is damaged, check the contained products carefully. Do not use if the sterile package is damaged, if the protective or closure caps are not in place or if there is any visible damage on the finished products (e.g. kinked tubes).

• Make sure that all tubes and chambers are correctly inserted into the respective holders. Avoid kinking or occluding the disposable in order to avoid mechanical and chemical damage to cellular blood constituents.

• For hygienic and functional reasons, it is recommended that tubing systems are inserted immediately before preparation only, thereby keeping the preparation and circulation times nearer the commencement of the treatment and in compliance with applicable guidelines.

• Connect the multiFiltratePRO tubings aseptically without touching open connectors.

• Disinfect with 70% alcohol access sites without protection cap prior to the connection with other products to reduce the risk of infection and let it dry before the connection.

• The compatibility of disinfectants (other than those recommended) with access sites shall be determined prior to clinical use.

• Excessive negative pressure may cause partial collapse of the pump segment resulting in an actual blood flow substantially less than indicated on the device.

• Minimum temperature of use of the disposable is 18°C (64°F).

• To ensure a secure connection between patient access and blood lines, hold and screw the coloured (blue, red) coupling nut on the blood line only. Do not apply the screwing torque to the inner part of the connector. After the connection, check that the components are firmly screwed together.

• Make sure that the components/tubes in direct interface with the device fit properly, and ensure that all caps and all screw-lock joints are secure and properly tightened (particularly those of the patient connections, the dialyzer connections, the device and the solution bags connections).

• Inspect the disposable for kinking and leaks during the priming and treatment phase and take corrective measures (e.g. tightening Luer-Lock connection) or exchanging it as necessary.

• To avoid air embolism, make sure that the return patient line is correctly inserted into the device's air bubble detector.

• Careful attention must be given to the installation and blood level of the bubble catchers. Make sure that due to pressure fluctuation neither the blood contacts the hydrophobic filter nor air enters the extracorporeal circuit.

• Change the disposable if a visible clot formation („white strip“) is seen in the blood line downstream of the return line bubble catcher, for instance at the point where the calcium infusion line merges into the blood line.

• Cleansing solutions and disinfectants may damage materials employed for the disposable. Safety and performance of use can no longer be guaranteed, and the manufacturer assumes no liability.

• The plastics used can be incompatible with drugs or disinfectants (e.g. connectors made of polycarbonate can develop cracks when in contact with aqueous solutions with pH > 10).

• If nutritional solutions are administered into the blood line, wetting of the luer-lock connection to the feeding system with lipidic fluids can weaken the properties of the plastic material used. Make sure that during the connection of the feeding system, the infusion line close to the luer-lock connection site remains completely free of nutrition solution.

Warnings: injection sites

• Do not use needles with a bevelled tip and a diameter larger than 20 gauge (outer needle diameter to be 0.9 mm or less) to puncture injection sites (if present).

Warnings: needle-free access sites

• Needleless access sites are labelled with the following:

• Use a Luer lock syringe to activate the needle-free access site. Do not use sharp needles to puncture the needle-free access site septum. Their use can lead to leaks or air entering the fluid pathway. Refer to the warning labels present on the finished product.

• Do not connect the syringe used for the needle-free access site to other parts of the finished product — if the syringe comes into contact with lubricant, it could become disconnected and lead to patient injury or death.

Operating time

• The maximum application time is 72 hours. The disposable must be replaced after the maximum usage time or respectively after a specific blood volume pumped as indicated on the primary packaging.

Particular notes on materials and substances

CMR substances and endocrine-disrupting substances

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please use this page: www.freseniusmedicalcare.com/en/svhc.

Special precautions for storage

Follow the indication of the product label. Protect from moisture, freezing and excessive heat.

Disposal

Ensure safe disposal of any unused product or waste material in accordance with local regulations. Materials that have been in contact with blood or other material of human origin may be infectious. Dispose of such materials by taking the necessary precautionary measures and in accordance with local regulations for (potentially) infectious materials.

SYMBOLS USED ON LABELS

	Medical Device
	Unique Device Identifier
	Patient information website
	Latex-free
	Sterilized using ethylene oxide
	Replace the Blood line after maximum usage time or pumped blood volume as indicated on the primary package
	Pump segment diameter / length
	Blood Priming Volume
	Units
	Temperature limit
	Do not use if package is damaged
	Sterile barrier system

DATE OF REVISION OF TEXT

07/2023

A digital version of this document is available on the following website: www.freseniusmedicalcare.com/en/product-information.

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