en

### INSTRUCTIONS FOR USE

## Tubing Systems for Continuous Blood Purification

### Models:

**multiFiltratePRO HDF**  
**multiFiltratePRO Ci-Ca HD**  
**multiFiltratePRO Ci-Ca HDF**  
**Adapter Hansen male / Luer-lock male**

Please read the following instructions carefully.

Refer to product or packaging label for:

- Type / Code of Bloodline/Tubing system
- Manufacturer and model of the machine for which the product is intended
- Applicable symbols
- Arterial / Access / Inlet Line (RED); Venous / Return / Outlet Line (BLUE); Filtrate Line (YELLOW); Dialysate Line (GREEN); Substitute Line / Plasma Line (WHITE)

#### SYMBOLS USED ON LABELS

Symbol	Meaning	Symbol	Meaning
	Medical Device		Do not use if package is damaged
	Unique Device Identifier		Consult instructions for use or consult electronic instructions for use
	Patient information website		Non-pyrogenic
	Bloodline/Tubing systems does not contain any natural Latex		Single Sterile Barrier System
	Sterile, Sterilized using ethylene oxide		Do not re-use
	Replace the Blood line after maximum usage time or pumped blood volume as indicated on the primary package		Caution (Consult the Instructions for Use for important cautionary information such as warnings and precautions)
	Pump segment diameter		Batch code
	Blood Priming Volume		Catalogue number
	Units		Date of manufacture
	Temperature limit (Indicates the temperature limits to which the medical device can be safely exposed)		Expiry date
	Manufacturer		

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Product composition

Product model	Product composition
multiFiltratePRO Ci-Ca HD	Arterial/venous tubing system, Filtrate tubing system, Dialysate tubing system, Citrate tubing (located in arterial tubing system) and calcium tubing system (located in venous tubing system)
multiFiltratePRO Ci-Ca HDF	Arterial/venous tubing system, Filtrate tubing system, Dialysate tubing system, Substitute tubing system, Citrate tubing (located in arterial tubing system) and calcium tubing system (located in venous tubing system)
multiFiltratePRO HDF*	Arterial/venous tubing system, Filtrate tubing system, Dialysate tubing system, Substitute tubing system
Adapter Hansen male/Luer-lock male	Cap, PVC tube, Connector Hansen male

\* multiFiltratePRO HDF model: dialysate tubing and substitute tubing share one feed line.

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
Structure diagram of compositions of tubing systems for continuous blood purification

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## GENERAL INFORMATION

### General description of the product

Bloodline/Tubing systems are intended to be used and compatible with an extracorporeal blood purification system multi (Version multiFiltratePRO) device.

Bloodline/Tubing systems are intended for single use.

This product is sterile and pyrogen free.

### Sterilization

This product is ethylene oxide sterilized (EO).

## COMPOSITION

Tubes: medical grade soft-PVC.

Connectors and other components: Acrylonitrile butadiene styrene (ABS), Polybutylene Terephthalate (PBT), Polycarbonate (PC), Polyamide (PA), Polyethylene (PE), Polyethylene Terephthalate Glycol-modified (PETG), Polypropylene (PP), Polytetrafluoroethylene (PTFE), Polyvinylchloride (PVC), Silicone, Polyhydroxy Ether (PHE), Thermoplastic Elastomer (TPE), Polyethylene High-Density (HDPE), Methyl Methacrylate Acrylonitrile Butadiene Styrene (MBS).

## INTENDED PURPOSE AND RELATED DEFINITIONS

### Intended purpose

Channeling of blood and fluid in an extracorporeal treatment.

### Intended Application

Bloodline/Tubing systems are intended for extracorporeal blood purification. multiFiltratePRO HDF model can be used in CVVH, CVVHD and CVVHDF treatment modes, multiFiltratePRO Ci-Ca HD model can be used in CVVHD mode, and multiFiltratePRO Ci-Ca HDF model can be used in CVVHDF treatment mode. The Adapter Hansen male/Luer-lock male model is used with the multiFiltratePRO HDF model for pre post CVVH treatment mode.

### Medical Indication

- Acute kidney injury leading e.g., to severe electrolyte disturbances and/or severe acid-base-balance disturbance requiring continuous renal replacement therapy
- Volume overload requiring continuous renal replacement therapy
- Intoxications requiring continuous renal replacement therapy

### Intended patient population

The devices have been specified by the manufacturer for the purpose of treating patients with a body weight of 40 kg and more, irrespective of their age, under consideration of the specified technical data of the device and the single-use items used (e.g., delivery rates, fill volumes).

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### Intended user group and intended environment

The disposable must only be used by medical individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown. Operate in rooms suitable for dialysis located in professional healthcare facilities.

Normative and local regulations must be observed.

### SIDE EFFECTS

Occasional occurrence of the following side effects is reported in current literature:

- Hypotension
- Blood loss (e.g. hypovolemia, and hypovolaemic shock, hypotension, anaemia)
- Coagulation disorders (e.g. thrombocytopenia)
- Embolism (air-, thrombo-)
- Bacterial infection
- Hypothermia
- Haemolysis
- Hypersensitivity or hypersensitivity-like reactions (e.g. dyspnea, cardiopulmonary arrest, hypotension, urticaria, flushing, itching, abdominal pain, nausea, convulsions and unconsciousness).

Additional side effects might be specific to other devices and drugs used in the therapy.

Refer to the Instructions for Use of the multi (Version multiFiltratePRO) device for more information on treatment related side effects.

### Reporting of serious incidents

If any serious incident occurs in relation to the device, including those not listed in this leaflet, the treating physician shall be informed immediately. Within the EU the user must report any serious incident that has occurred in relation to the device to the manufacturer according to labelling (■) and the competent authority of the EU Member State in which the treatment is performed. When the device is used outside the EU, the serious incident in relation to the device should also be reported to the local authority where the treatment is performed.

A serious incident can be any incident that directly or indirectly leads to the death of a patient, user or other person; to the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health treat.

### CONTRAINDICATIONS

#### Product related contraindications

Hypersensitivity to any of the material of the tubings as listed under composition or to ethylene oxide residues.

In severe cases the treatment must be discontinued, the blood should not be returned to the patient, and appropriate emergency medical treatment should be initiated.

#### Therapy related contraindication

Refer to the Instructions for Use of the multi (Version multiFiltratePRO) device for more information of general contraindications for extracorporeal blood purification.

### PERFORMANCE CHARACTERISTICS

Refer to the Instructions for Use of the multi (Version multiFiltratePRO) device for essential performance parameters.

### METHOD OF ADMINISTRATION

#### Handling Instructions

Refer to the Instructions for Use of the multi (Version multiFiltratePRO) device of the regarding how to handle the disposable during set-up (e.g. selection of treatment type), priming, treatment and reinfusion. Also, the Instructions for Use of other disposables such as the haemofilter, plasmafilter must be taken into account.

#### Preparation

- Refer to the Instructions for Use of the multi (Version multiFiltratePRO) device for details on any preparatory handling of the disposable before it is ready for use or during its use and information to verify if the disposable is properly installed.
- The Bloodline/Tubing systems are designed for machines for extracorporeal blood purification multiFiltrate or multi (version: multiFiltratePRO) (refer to label) and shall be used only after appropriate instruction or training. Use only Bloodline/Tubing systems approved for the selected therapy.
- The infusion administration lines should be clamped, except when needed.
- Unpack and connect the Bloodline/Tubing systems aseptically without touching open connectors.
- Colour codes should be followed and used in line with the corresponding markings on the device.
- Handle pumps and pump clips in accordance to the Instructions for Use of the device.
- Tighten all closure caps and ensure that all connectors and caps are secured.
- Fill and rinse the disposable in accordance with the Instructions for Use of the device or supplemental Instructions for Use for additional options, and training, if applicable.
- Fresenius Medical Care disposable is designed to withstand the maximum and minimum manufacturer's recommended pressures and flow rates generated in use with the respective multi (Version multiFiltratePRO) device. All the other relevant specifications of the device apply.

#### Treatment

- In cases where the extracorporeal volume of the circuit consists of a significant amount of the patient's circulating blood volume, as well as in anaemic and/or haemodynamically unstable patients, priming the circuit with blood or blood substitution fluid may be considered.
- To ensure correct functionality of the bubble catchers, fill to about

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1 cm below the top. This avoids fluid contacting the hydrophobic filter as well as air entering the extracorporeal circuit. In the event that fluid reaches and comes into contact with the hydrophobic filter within the Transducer Protector (TP) the disposable must be replaced.

- Bloodline/Tubing systems that have the possibility of an additional connection intended to be connected to the pressure measurement system. This can be used to connect a new pressure measurement line (available as an accessory in the Fresenius Medical Care product range).
- Do not push back the fluid in the pressure measurement line by using a syringe. This might damage the membrane in the TP and thus lead to contamination. In the event that fluid might have passed through the TP, check the device for contamination after treatment is completed. If contamination has occurred, the device must be taken out of service and disinfected according to the manufacturer's recommendations before further use.
- After the manipulation of lines or the usage of components during the treatment check and, if necessary, restore the correct position of the lines and components.
- Disinfect the corresponding access sites without protection cap prior the connection with other products.

#### Reinfection

- Refer to the Instructions for Use of the multi (Version multiFiltratePRO) device for the termination of the treatment and also to the "Disposal" Section of the present Instructions for Use.
- Close all clamps on the Bloodline/Tubing systems before removing the disposable to reduce the risk of fluid leakage.

#### WARNINGS AND PRECAUTIONS

##### Warnings

- The disposable must only be used by individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown.
- The Bloodline/Tubing systems are intended for single use only. The correct function of all interfaces is ensured only for single use. Re-use may be hazardous to both patient and operator (e.g. impaired performances, contamination).
- Do not use on and after expiry date (refer to label).
- If the carton is damaged, check the contained products carefully. Do not use if the sterile package is damaged, if the protective or closure caps are not in place or if there is any visible damage on the finished products (e.g. kinked tubes).
- Make sure that all tubes and chambers are correctly inserted into the respective holders. Avoid kinking or occluding the disposable in order to avoid mechanical and chemical damage to cellular blood constituents.
- For hygienic and functional reasons, it is recommended that tubing systems are inserted immediately before preparation only, thereby keeping the preparation and circulation times nearer the commencement of the treatment and in compliance with applicable guidelines.
- Connect the Bloodline/Tubing systems aseptically without touching open connectors.
- Disinfect with 70% alcohol access sites without protection cap prior to the connection with other products to reduce the risk of infection and let it dry before the connection.
- The compatibility of disinfectants (other than those recommended) with access sites shall be determined prior to clinical use.
- Excessive negative pressure may cause partial collapse of the pump segment resulting in an actual blood flow substantially less than indicated on the device.
- Minimum temperature of use of the disposable is 18°C (64°F).
- To ensure a secure connection between patient access and blood line, hold and screw the colored (blue, red) coupling nut on the blood line only. Do not apply the screwing torque to the inner part of the connector. After the connection, check that the components are firmly screwed together.
- Make sure that the components/tubes in direct interface with the device fit properly, and ensure that all caps and all screw-lock joints are secure and properly tightened (particularly those of the patient connections, the dialyzer connections, the device and the solution bags connections).
- Inspect the disposable for kinking and leaks during the priming and treatment phase and taking corrective measures (e.g. tightening Luer-Lock connection) or exchanging it as necessary.
- To avoid air embolism, make sure that the return patient line is correctly inserted into the device's air bubble detector.
- Careful attention must be given to the installation and blood level of the bubble catchers (see above). Make sure that due to pressure fluctuation neither the blood contacts the hydrophobic filter nor air enters the extracorporeal circuit.
- Change the disposable if a visible clot formation ("white strip") is seen in the blood line downstream of the return line bubble catcher, for instance at the point where the calcium infusion line merges into the blood line.
- Cleansing solutions and disinfectants may damage materials employed for the disposable. Safety and performance of use can no longer be guaranteed, and the manufacturer assumes no liability.
- The plastics used can be incompatible with drugs or disinfectants (e.g. connectors made of polycarbonate can develop cracks when in contact with aqueous solutions with pH > 10).
- If nutritional solutions are administered into the blood line, wetting of the luer-lock connection to the feeding system with lipidic fluids can weaken the properties of the plastic material used. Make sure that during the connection of the feeding system, the infusion line

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close to the luer-lock connection site remains completely free of nutrition solution.

#### Warnings: injection sites

- Do not use needles with a beveled tip and a diameter larger than 20 gauge (outer needle diameter to be 0.9 mm or less) to puncture injection sites (if present).

#### Warnings: needle-free access sites

- Needleless access sites are labelled with the following:

- Use a Luer lock syringe to activate the needle-free access site. Do not use sharp needles to puncture the needle-free access site septum. Their use can lead to leaks or air entering the fluid pathway. Refer to the warning labels present on the finished product.
- Do not connect the syringe used for the needle-free access site to other parts of the finished product — if the syringe comes into contact with lubricant, it could become disconnected and lead to patient injury or death.

#### Operating time

- The maximum application time is 72 hours. The disposable must be replaced after the maximum usage time or respectively after a specific blood volume pumped as indicated on the primary packaging.

#### Particular notes on materials and substances

##### CMR substances and endocrine-disrupting substances

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please use this page: [www.freseniusmedicalcare.com/en/svhc](http://www.freseniusmedicalcare.com/en/svhc).

#### Special precautions for storage

Follow the indication of the product label. Protect from moisture, freezing and excessive heat.

#### Disposal

Ensure safe disposal of any unused product or waste material in accordance with local regulations.

Materials that have been in contact with blood or other material of human origin may be infectious. Dispose of such materials by taking the necessary precautionary measures and in accordance with local regulations for (potentially) infectious materials.

#### WARRANTY

The manufacturer shall not be liable for any misuse, improper handling, non compliance with instructions for use and cautionary notes, and for any damage incurred subsequent to the manufacturer's delivery of the Bloodline/Tubing systems.

#### DATE OF REVISION OF TEXT

2023-07 Fifth Modified

A digital version of this document is available on the following website:

[www.freseniusmedicalcare.com/en/product-information](http://www.freseniusmedicalcare.com/en/product-information).

#### INFORMATION ON THE MANUFACTURER

Legal manufacturer

Registrant/Manufacturer: Fresenius Medical Care AG

Address of Registrant/Manufacturer: Else-Kröner-Str.1, 61352 Bad Homburg, GERMANY

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注册人/生产企业住所: Else-Kröner-Str.1, 61352 Bad Homburg, GERMANY

**Telephone of Registrant/Manufacturer:** +49 6172 609-0

**Manufacturing Address(es):** Via Crema 8 26020 Palazzo Pignano (CR), ITALY; No: 16, Liman Serbest Bölgesi Mahallesi 07070 Antalya, TURKEY; Stegne 11 Ljubljana 1000, SLOVENIA

**Agent / After-Sale Service Provider in China:** Fresenius Medical Care (Shanghai) Co., Ltd.

**Address of Agent/After-Sale Service Provider:** Section A, Building 01, No. 439 Futexi Road 1st, China (Shanghai) Pilot Free Trade Zone

**Telephone Number of Agent/After-Sale Service Provider:** 021-61152800

**Medical Device Registration Certificate No./Product Technical Requirements No.:** GuoXieZhuJin20203100035

**Date of Manufacture:** See outer package

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