


Dialysis Water Distribution Loop Dialysis Concentrate Distribution Loop



Instructions for Use

Edition: 02A-2023
Date of issue: 2023-08
Part no.: F50015018



**FRESENIUS
MEDICAL CARE**

Table of contents

1	Index	7
2	Important information	9
2.1	How to use the Instructions for Use	9
2.2	Warnings, significance	10
2.3	Notes, significance	10
2.4	Tips, significance	10
2.5	Brief description	11
2.6	Intended purpose and related definitions	12
2.6.1	Intended purpose.....	12
2.6.2	Medical indication.....	12
2.6.3	Intended patient population.....	12
2.6.4	Intended user group and intended environment.....	12
2.7	Side effects	13
2.8	Contraindications	13
2.9	Residual risks	14
2.10	Interaction with other systems	15
2.10.1	Intended combined use.....	15
2.11	Therapy restrictions	15
2.12	Considerations for working on the DWDL and DCDL	15
2.13	Expected service life	16
2.14	Duties of the responsible organization	16
2.14.1	Further aspects of the responsible organization.....	16
2.15	Operator responsibility	17
2.15.1	Reporting incidents.....	17
2.16	Disclaimer of liability	18
2.17	Technical documentation	18
2.18	Warnings	19
2.18.1	Basic warnings.....	19
2.18.2	Warnings related to hygiene and biology.....	21
2.19	SVHC (REACH)	22
2.20	Addresses	23
3	Structure and views	25
3.1	Overall view of the distribution systems	25

3.2	Overall view of DWDL.....	26
3.3	Overall view of DCDL	26
4	Operation	27
4.1	Overview of coupling systems on the media supply systems.....	27
4.1.1	MediaC	27
4.1.2	MediaP	27
4.1.3	MediaR.....	27
4.2	Observe when operating.....	28
4.3	Operating the coupling systems	28
4.3.1	To be observed before connecting the hemodialysis device	28
4.3.2	Connecting and disconnecting.....	29
4.4	Microbiological analysis at the sampling port.....	31
4.4.1	Preparation	31
4.4.2	Accessories, equipment.....	31
4.4.3	Procedure for collecting a sample at the sampling port.....	32
4.5	Microbiological analysis at the dialysis water connection	34
4.5.1	Preparation	34
4.5.2	Accessories, equipment.....	34
4.5.3	Procedure for collecting a sample at the dialysis water connection	35
4.6	Collecting a sample for chemical analysis.....	36
4.6.1	Preparation	36
4.6.2	Accessories, tools.....	36
4.6.3	Procedure for collecting a sample for chemical analysis	37
5	Alarms	39
6	Cleaning, disinfection	41
6.1	Generally applicable regulations for cleaning and disinfection	41
6.2	Precautions	41
6.2.1	Operator safety	41
6.3	Surface cleaning, surface disinfection	43
6.3.1	Surface cleaning	43
6.3.1.1	Cleaning agents for surfaces.....	43
6.3.2	Surface disinfection	43
6.3.2.1	Surface disinfectant.....	44
6.4	Disinfecting the DWDL	45
6.4.1	Reasons for disinfecting the DWDL.....	45
6.4.2	Disinfecting the system.....	45
6.5	Cleaning connectors and connection ports.....	47
6.5.1	Cleaning the connectors	47
6.5.2	Cleaning the connection ports	47

6.6	Disinfecting connectors and connection ports	48
6.6.1	Disinfecting the connectors	48
6.6.2	Disinfecting the connection ports.....	49
7	Functional description	51
7.1	Description of procedures	51
7.1.1	Functions	51
8	Consumables, accessories, and additional equipment	53
8.1	Consumables	54
8.2	Accessories	55
8.2.1	Accessories for the DCDL distribution system	55
8.2.2	Accessories for the DWDL distribution system	56
8.3	Additional equipment	58
8.3.1	Additional equipment for the DWDL distribution system	58
8.3.2	Additional equipment for the DCDL distribution system	58
9	Installation	59
9.1	Installation requirements	59
9.1.1	General information	59
9.1.2	Environmental conditions	59
9.1.3	Structural installation requirements	59
9.2	Operational qualification	60
9.2.1	To be observed before the operational qualification.....	60
9.2.2	Operational Qualification procedure	60
9.3	Removal from service, decommissioning	61
9.3.1	Removal from service.....	61
9.3.2	Decommissioning	61
10	Transport and storage	63
10.1	Transport and storage conditions	63
10.2	Environmental compatibility/disposal	64
10.2.1	Handling of disinfectants	64
11	Technical Safety Checks and maintenance	65
11.1	Important information for performing the TSC/MA	65
11.2	Maintenance procedures	66
11.2.1	Quality assurance and care measures	66

12	Specifications	67
12.1	DWDL and DCDL product data	67
12.2	Identification label (DWDL and DCDL identification)	68
12.3	Transport/storage	69
12.4	Materials used	69
13	Definitions	71
13.1	Definitions and terms	71
13.2	Abbreviations	72
13.3	Symbol	73
13.4	Certificates	74
14	Options	75
15	Appendix	77
15.1	Medical Device Register for the DWDL and DCDL	77
15.1.1	Responsible organization and identification	77
15.1.2	Contents of the Medical Device Register for the DWDL and DCDL	79
15.2	Training Record	81
15.3	Quality of dialysis water	86

1 Index

A

Abbreviations 72
Accessories 55
Addresses 23
Alarms 39
Appendix 77

B

Brief description 11

C

Certificates 74
Chemical quality of dialysis water 87
Consumables 54
Contraindications 13

D

Decommissioning 61
Definitions 71
Disclaimer of liability 18
Duties of the responsible organization 16

E

Environmental compatibility/disposal 64
Expected service life 16

F

Functional description 51

I

Important information 9
Installation 59
Intended combined use 15
Intended patient population 12
Intended purpose 12
Intended purpose and related definitions 12

Intended user group and intended environment 12
Interaction with other systems 15
International service support 23

L

Local service support 23

M

Maintenance procedures 66
Medical Device Register 79
Medical indication 12
Microbiological quality of fluids for hemodialysis 86

N

Notes, significance 10

O

Operation 27
Operational qualification 60, 71
Operator responsibility 17
Operator safety 41
Options 75

Q

Quality assurance and care measures 66
Quality of dialysis water 86

R

Removal from service 61
Residual risks 14
Restrictions 15

S

Service life 16
Side effects 13
Specifications 67
Surface cleaning 43

Surface disinfection 43
SVHC (REACH) 22
Symbol 73

T

Technical documentation 18
Technical Safety Checks and maintenance 65
Tips, significance 10
Training Record 81
Transport and storage 63

W

Warnings 19
Warnings related to hygiene and biology 21
Warnings, significance 10

2 Important information

2.1 How to use the Instructions for Use

Description	In this document, the dialysis water distribution system is referred to as Dialysis Water Distribution Loop (DWDL) and the dialysis concentrate distribution system is referred to as Dialysis Concentrate Distribution Loop (DCDL) .				
Identification	The document can be identified by the following information on the title page and on the labels, if any: <ul style="list-style-type: none"> – Document edition – Date of issue of the document – Document part number 				
Footer	The footer contains the following information: <ul style="list-style-type: none"> – Company name – Product name – The abbreviation for the document type and the international abbreviation for the document language, e.g., IFU-EN means Instructions for Use in English. – Edition information, e.g., 04A-2021 means edition 04A from the year 2021 – Page identification 				
Organization of the chapters	To facilitate the use of documents from Fresenius Medical Care, the organization of the chapters has been standardized in all manuals. There may therefore be chapters within this document without any content. Chapters without content are identified.				
Styles used in the document	The following styles may be used in the document: <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 10px;"> <thead> <tr style="background-color: #ADD8E6;"> <th style="text-align: left; padding: 5px;">Style</th> <th style="text-align: left; padding: 5px;">Description</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">➤ Instruction</td> <td style="padding: 5px;">Instructions are indicated by an arrow ➤. Instructions must be followed. Example: ➤ Carry out instruction.</td> </tr> </tbody> </table>	Style	Description	➤ Instruction	Instructions are indicated by an arrow ➤. Instructions must be followed. Example: ➤ Carry out instruction.
Style	Description				
➤ Instruction	Instructions are indicated by an arrow ➤. Instructions must be followed. Example: ➤ Carry out instruction.				
Illustrations	The illustrations used in the documents may differ from the original if this does not have any influence on the function.				
Importance of the instructions	The Instructions for Use are part of the accompanying documents and are an essential part of the DWDL and the DCDL . They include all information necessary for the use of the distribution systems. The Instructions for Use must be carefully studied before operational qualification/start-up of the DWDL and DCDL distribution systems.				
Changes	Changes to the technical documents will be released as new editions or supplements. In general, this document is subject to change without notice.				

Reproduction

Reproduction, even in part, is only permitted with written approval.

2.2 Warnings, significance

Advises the operator that failure to comply with the measures for avoiding the hazard may result in serious or fatal personal injuries.

Warning

Type and cause of the hazard

Possible consequences if the hazard arises.

➤ Measures for avoiding the hazard.

Warnings can deviate from the example above in the following cases:

- If a warning refers to several hazards.
- If a warning cannot be assigned to one particular hazard.

2.3 Notes, significance



Note

Advises the operator that failure to observe this information can result in the following:

- Damage to the **DWDL** and **DCDL** distribution systems.
 - Specific functions not being executed at all or not being executed correctly
-


2.4 Tips, significance



Tip

Information providing useful tips for easy handling.

2.5 Brief description



The distribution systems are classified as class IIa (MDR) equipment.

The distribution systems are connection units which can be extended by the responsible organization with additional components to form a complete supply system.

They transport dialysis water via the **Dialysis Water Distribution Loop** distribution system and dialysis concentrate via the **Dialysis Concentrate Distribution Loop** distribution system.

The transported dialysis water can be used for dialysis treatments or for the production of dialysis concentrates. The transported dialysis concentrate can be used for mixing common dialysates.

2.6 Intended purpose and related definitions



Note

If the distribution system is used in a manner which is not specified by the manufacturer, the quality and characteristics of the dialysis water and dialysis concentrate transported by the distribution system may be impaired.

2.6.1 Intended purpose

Dialysis Water Distribution Loop (DWDL): Central supply of dialysis water

Dialysis Concentrate Distribution Loop (DCDL): Central supply of dialysis concentrate

2.6.2 Medical indication

DWDL: Renal insufficiency requiring renal replacement therapy, supported by a central supply of dialysis water.

DCDL: Renal insufficiency requiring renal replacement therapy, supported by a central supply of dialysis concentrate.

2.6.3 Intended patient population

The distribution systems **DWDL** and **DCDL** do not have a clinical effect of their own. The devices solely transport the dialysis water (**DWDL**), respectively dialysis concentrate (**DCDL**) required for the preparation of standard dialysates. Thus, there are no limitations for an intended patient population. The intended patient population should be defined by the compatible hemodialysis device.

2.6.4 Intended user group and intended environment

The distribution systems **DWDL** and **DCDL** must only be installed, operated and used by individuals with the appropriate training, knowledge and experience, and who are certified to have been trained. The distribution systems **DWDL** and **DCDL** must be operated in rooms suitable for the operation of hemodialysis devices located in professional health-care facilities.

2.7 Side effects

There are no side effects that can solely be traced back to the use of dialysis water and dialysis concentrate as they do not have a direct clinical effect on its own. Dialysis water and dialysis concentrate are always used in combination with a hemodialysis treatment. An elevated level of calcium, magnesium and iron in the dialysis water may lead to the hard water syndrome, resulting in nausea, vomiting, asthenia, and/or hypertension.

For reference the following treatment-related side effects known for hemodialysis, as reported in current literature, are listed:

- Acute urticaria
- Anxiety
- Impaired quality of life
- Clotting
- Blood loss
- Depressive symptoms
- Dialysis disequilibrium syndrome
- Thirst
- Vomiting
- Fever
- Hemolysis
- Hypotension
- Itching
- Cardiac arrhythmia
- Headache
- Seizures
- Cramps
- Micro air embolisms
- Cardiac tamponade
- Dialyzer reactions
- Sleep disturbance
- Pain (chest and back)
- Shivering
- Falls
- Nausea
- Restlessness

2.8 Contraindications

There are no known contraindications. Dialysis water and dialysis concentrate are never used directly on the patient. Contraindications may exist in the therapy-related context of hemodialysis:

- Hyperkalemia (only with potassium-containing hemodialysis concentrates)
- Hypokalemia (only with potassium-free hemodialysis concentrates)
- Uncontrollable blood-clotting disorders

Relative contraindications (predictors for poor treatment outcome/treatment decision on an individual basis):

- Hypotensive heart failure
- Malignant disease with poor prognosis
- Severe peripheral arterial disease (no access possible)
- Severe mental illness to the extent that patient is not aware of treatment and cannot comply.

A different method of extracorporeal treatment may be indicated for hemodynamically unstable patients.

2.9 Residual risks

The central supply of dialysis water and dialysis concentrate via the distribution systems carries a residual risk of microbial contamination, germ growth and the spread of pathogens. Patients may suffer from infection or sepsis.

The use of chemical agents for cleaning and disinfecting the distribution systems entails the residual risk of chemical residues in the dialysis water and dialysis concentrate. Patients may be exposed to toxic or pyrogenic substances.

2.10 Interaction with other systems

2.10.1 Intended combined use

The medical devices **DWDL** and **DCDL** are intended for use and combination with the following devices:

- **DWDL**
 - AquaA
 - AquaBplus, AquaB LITE
 - Granumix 107S/507S
 - Granumix plus
 - MediaP
 - MediaC
 - MediaR

- **DCDL**
 - CDS3
 - Granumix plus
 - MediaP
 - MediaC
 - MediaR

2.11 Therapy restrictions

none

2.12 Considerations for working on the DWDL and DCDL

Warning

Risk of injury to patients and operators from improper service work on the distribution system

Improper servicing may impair the safe functioning of the distribution system.

- Operational qualification, extensions, maintenance procedures, modifications, or repairs may only be carried out by the manufacturer or by persons authorized by the manufacturer.

To perform the Technical Safety Checks and the maintenance procedures, contact the local service support.

Use only original spare parts. To identify and order spare parts, test equipment, and tools, always use the electronic Spare Parts Catalog.

- Further information on installation (see chapter 9 on page 59).
- Further information about maintenance procedures (see chapter 11.2 on page 66).
- Further information about transport and storage (see chapter 10 on page 63).

2.13 Expected service life

The expected service life of the distribution systems is 10 years.

If the maintenance procedures and checks are performed as prescribed and at the defined intervals, the distribution system will run safely between checks.

2.14 Duties of the responsible organization

The responsible organization must ensure that the following requirements are met:

- Compliance with the national or local regulations concerning the installation, operation, use, and maintenance of the device.
- Compliance with the accident prevention regulations.
- Ensuring that the distribution system is in proper and safe working order.
- Ensuring that the Instructions for Use are available at all times.
- The distribution system may only be operated under the operating conditions specified by the manufacturer.

2.14.1 Further aspects of the responsible organization

The responsible organization must ensure that the technical design of the distribution systems matches the requirements of the other components used to make up the complete system.

The responsible organization must prepare a plan for emergency operation to supply dialysis devices with dialysis water and dialysis concentrate based on the available system components and must make this plan available to the operators of the system.

The responsible organization should inform the local water supplier of the dialysis operation and insist on prior discussion with regards to water composition, availability, etc. This measure does not relieve the responsible organization of its obligation to regularly check the inlet water composition.

Before the responsible organization may start operating the distribution systems, the individual responsible for operation must have been verifiably instructed by the manufacturer in how to use the systems and must be thoroughly familiar with the contents of the Instructions for Use. The manufacturer provides training for distribution systems.

The local service support organization is available to answer any further questions (see chapter 2.20 on page 23).

Disinfection DCDL

The **DCDL** is not designed for disinfection.

Disinfection DWDL

Bacterial growth in the distribution system depends on the individual components, and the type and time of use. Bacterial growth in the distribution systems must be prevented by continuous operation of the system with a minimum of idle times and by preventive measures such as chemical disinfection or heat disinfection.

Samples for microbial testing must therefore be collected from the distribution system and from the individual parts of the system in accordance with the applicable regulations. As the complete system consists of a number of smaller systems, the responsible organization is responsible for the complete system.

For further information on disinfecting the **DWDL** (see chapter 6.4.1 on page 45).

2.15 Operator responsibility

Warning

Risk of injury from defects on the distribution systems

If the distribution systems have the following defects, the indicated measures must be taken.

Defects on the distribution systems:

- Mechanical damage
- Deterioration of performance characteristics: pressure drop, leakages
- Other defects

Measures:

- The distribution system must be taken out of service.
 - The responsible organization or the local service support must be notified.
-

2.15.1 Reporting incidents

Within EU Member States, the user must report any severe incidents associated with the product to the manufacturer in accordance with the identification as well as the responsible authority of the Member State in which the user is located.

2.16 Disclaimer of liability



Warning

Risks affecting the proper functioning of the device

The distribution system has been approved for use with certain consumables and accessories. Should the responsible organization wish to use other consumables and accessories than those listed in this chapter, the suitability must be checked beforehand by gathering the appropriate manufacturer information. The applicable legal regulations must be complied with.

The manufacturer does not assume any responsibility or liability for personal injury or other damage, and the use of non-approved or unsuitable consumables or accessories resulting in damage to the distribution system will void the warranty.



Tip

For further information on the topic of consumables, accessories, additional equipment (see chapter 8 on page 53).

2.17 Technical documentation

Descriptions and other technical documents can be provided by the manufacturer on request. These are intended to support appropriately trained personnel of the responsible organization in maintaining and repairing the system.

2.18 Warnings

The list of warnings and notes that follows is only an excerpt. Knowledge of all warnings mentioned in this document is required for the safe operation of the distribution systems.

2.18.1 Basic warnings



Warning

Undetected fluid leakage outside of dialysis hours

Leakage can cause damage to buildings.

- To prevent damage to buildings outside dialysis hours (unattended times without staff) caused by water leakage, a leakage monitoring system with a shut-off function, such as the **AquaDETECTOR** with leakage sensors, should be installed in every room with supply ports.
 - If no leakage monitoring system is installed, it is recommended that all supply tubes be disconnected from the distribution system outside dialysis hours (unattended times without staff).
-



Note

The responsible organization must ensure that the Technical Safety Checks (TSC) are performed.



Warning

Risk of injury to patients and operators if the TSC intervals are not complied with

Non-compliance with the TSC intervals may impair the safe functioning of the distribution system.

- The Technical Safety Checks/maintenance procedures (local service) must be performed on the distribution system at least once every 24 months.
 - The Technical Safety Checks and maintenance procedures may only be performed by certified service technicians with electrical, system-related, and medical-technical knowledge.
-



Note

Adherence to applicable laws and regulations

- Observe all applicable local laws and regulations concerning the handling of laboratory equipment and reagents.
-



Warning

Risk of burning and scalding from hot surfaces or hot dialysis water during heat disinfection

Contact with hot surfaces or hot dialysis water can cause burning or scalding.

- Do not touch the accessible components of the dialysis water distribution system while heat disinfection is in progress.
- Do not attempt to remove the dialysis water manually while heat disinfection is in progress.



Warning

Undetected fluid leakage due to insufficient checks

Leakage can cause damage to buildings.

- Regular visual inspection and leakage checks of all tubing, connectors, and piping containing fluid are required.
- Tubing must be protected against possible mechanical damage.



Note

The distribution system is not designed for supporting additional loads.

2.18.2 Warnings related to hygiene and biology



Warning

Risk of contamination from unsuitable dialysis water

There is a risk of spreading germs.

- Check the dialysis water quality at regular intervals, and, if necessary, perform the disinfection/cleaning cycles of the dialysis water supply system.



Warning

Risk of poisoning – Not drinking water

The dialysis water and dialysis concentrate transported via the distribution system do not meet the requirements for drinking water.



Warning

Risk of contamination from insufficient cleaning/disinfection

There is a risk of spreading germs.

- The distribution system may only be cleaned and disinfected by persons who have been instructed on the proper handling of the system during these procedures.
- The operator must observe and follow the general safety precautions.
- The distribution system may only be disinfected after consultation with the manufacturer of the system or persons authorized by the manufacturer.



Warning

Risk of chemical burns when working with acidic substances (disinfectant / cleaning agent)

Contact with chemicals can cause chemical burns.

- Always handle acidic fluids with care and do not spill any disinfectant concentrate.
- Wear appropriate personal protective equipment (gloves, goggles, etc.) in line with the safety precautions for the disinfectant / cleaning agent used.
- Observe the safety precautions for the disinfectant / cleaning agent used, including the relevant first aid measures.



Note

Risk of infection

Observe the applicable local laws and regulations concerning the handling of potentially infectious material.

2.19 SVHC (REACH)

Information on the subject of SVHC in accordance with Article 33 of Regulation (EC) 1907/2006 ("REACH") is available from the following website:

www.freseniusmedicalcare.com/en/svhc



2.20 Addresses

Manufacturer

Fresenius Medical Care AG & Co. KGaA
Else-Kröner-Str. 1
61352 Bad Homburg
GERMANY
Phone: + 49 6172 609-0
www.freseniusmedicalcare.com

International service support


Fresenius Medical Care
Deutschland GmbH
Technical Operations
Technical Coordination Office (TCO)
Hafenstrasse 9
97424 Schweinfurt
GERMANY

Local service support



3 Structure and views


3.1 Overall view of the distribution systems




Legend

- Green** **Dialysis Concentrate Distribution Loop (DCDL)**
- Blue** **Dialysis Water Distribution Loop (DWDL)**

3.2 Overall view of DWDL



3.3 Overall view of DCDL



4 Operation

4.1 Overview of coupling systems on the media supply systems


4.1.1 MediaC



Blue: dialysis water connection

Gray: drain connections

4.1.2 MediaP



Yellow: gas connections

Blue: dialysis water connections

Gray: drain connections

Green: dialysis concentrate connections

4.1.3 MediaR



Green: dialysis concentrate connections

Blue: dialysis water connection

Gray: drain connection

4.2 Observe when operating



Tip

- For basic information regarding **Definitions**, refer to chapter 13 (see chapter 13.1 on page 71).
-



Tip

- For basic information regarding **Product data**, refer to chapter 12 (see chapter 12.1 on page 67).
-



Tip

- For basic information regarding **Cleaning and disinfection**, refer to chapter 6 (see chapter 6 on page 41).
-

4.3 Operating the coupling systems

Warning**Patient hazard as a result of a disorder of the electrolyte balance because of a wrong dialysate composition**

Mixing up dialysis concentrates can result in a dialysis fluid that is unsuitable for the patient.

- When connecting the hemodialysis device to the media supply system, ensure that the concentrate connections are allocated correctly.
-

4.3.1 To be observed before connecting the hemodialysis device

Warning**Risk of contamination as a result of improper handling of connection sites**

There is a risk of spreading germs.

- Before connecting the hemodialysis device to the media supply system, disinfect the couplings and counter-couplings. Observe the Instructions for Use for the hemodialysis device.
 - It is absolutely essential to avoid contaminating the connections through contact with skin or other non-sterile objects.
-

Warning**Leakage due to damage on seals**

Leakage can cause damage to buildings.

- Before connecting to the corresponding counter-couplings, carefully check all dialysis concentrate couplings for salt crystals and, if necessary, clean and disinfect them. Observe the Instructions for Use for the hemodialysis device.
-

**Note****To be observed before connecting**

- The connectors must be disinfected before connecting to prevent possible contamination (see chapter 6.6 on page 48).
-

**Note****Cleaning and disinfecting the connectors**

- For information about cleaning and disinfecting the connectors and using the recommended disinfectants, (see chapter 6.5 on page 47) and (see chapter 6.6 on page 48).
-

4.3.2 Connecting and disconnecting

Warning**Risk of leakage as a result of improper handling of the coupling system**

If the lock is not released completely, this can cause increased wear and, as a consequence, fluid leakage.

- To connect and disconnect the coupling system, push the sleeve of the lock as far as it will go.
-



Recommended equipment

- Rubber gloves

**Step-by step
connecting/disconnecting**

- Grasp the sleeve of the coupling.
- Push the sleeve to the end position and hold it in this position.
- Push the coupling as far as it will go onto the nipple and release the sleeve.
- Disconnect the coupling by following the connection steps in reverse order.




Note

To be observed during disconnecting:

- Do not place disconnected tubes onto the floor and do not connect with other tubes.
 - Hang the tubes on a suitable support with the opening pointing down.
-

4.4 Microbiological analysis at the sampling port

The sampling port of the **DWDL** is the sampling valve, which can be opened by turning the valve.



4.4.1 Preparation

- Have a cooled shipping box available.
- The reverse osmosis system must be operated in the **RINSE** or **SUPPLY** mode for at least 20 minutes before the sample is collected.
- While collecting the sample, the reverse osmosis system must be in the **RINSE** or **SUPPLY** program.
- Collect the sample according to the procedure described for collecting a sample at the sampling port (see chapter 4.4.3 on page 32).

4.4.2 Accessories, equipment

The following equipment is recommended by the manufacturer:

- Rubber gloves
- Alcohol-based disinfectant (approx. 70–80 % alcohol, without lipid replenisher)

The **sampling set for the Fresenius sampling valve** (F00010382) can be used for sampling the dialysis water.

4.4.3 Procedure for collecting a sample at the sampling port








Illustration	Description
 <p data-bbox="180 741 252 772">Fig. 1</p>	<p data-bbox="778 389 1294 421">Fig. 1 – Disinfecting the sampling valve:</p> <ul style="list-style-type: none"> <li data-bbox="778 436 1426 499">➤ Disinfect the sampling valve using an alcohol-based disinfectant (without lipid replenisher). <li data-bbox="778 515 1347 546">➤ Wipe away any contamination using a swab. <li data-bbox="778 562 1378 593">➤ Then repeat the disinfection procedure (Fig. 1). <p data-bbox="778 609 1324 672">Caution: Observe the dwell time of the disinfectant!</p>
 <p data-bbox="132 1014 204 1046">Fig. 2</p> <p data-bbox="132 1270 204 1301">Fig. 3</p>	<p data-bbox="778 804 1433 835">Fig. 2, Fig. 3 – Connecting and locking the adapter:</p> <ul style="list-style-type: none"> <li data-bbox="778 851 1433 913">➤ Place the adapter of the sampling bag onto the sampling valve (Fig. 2). <li data-bbox="778 929 1426 1023">➤ Then lock the adapter (Fig. 3). The multiway valve on the sampling set must be set so as to ensure that no fluid can flow.
 <p data-bbox="247 1738 319 1769">Fig. 4</p>	<p data-bbox="778 1397 1238 1429">Fig. 4 – Opening the sampling valve</p> <ul style="list-style-type: none"> <li data-bbox="778 1444 1426 1507">➤ Turn the sampling valve counterclockwise to open it (Fig. 4).

Illustration	Description
 <p data-bbox="411 618 485 651">Fig. 5</p>	<p data-bbox="834 293 1283 327">Fig. 5 – Rinsing the sampling valve</p> <ul style="list-style-type: none"> <li data-bbox="834 338 1477 405">➤ Turn the multiway valve clockwise by 90° to ensure a flow from the adapter to the rinse tubing. <li data-bbox="834 416 1477 483">➤ Rinse the sampling valve for approx. 60 seconds via the rinse tubing.
 <p data-bbox="411 1028 485 1061">Fig. 6</p>	<p data-bbox="834 680 1110 714">Fig. 6 - Filling the bag</p> <ul style="list-style-type: none"> <li data-bbox="834 725 1477 792">➤ Then turn the multiway valve again 90° clockwise to fill the bag. <li data-bbox="834 804 1477 904">➤ Caution: Make sure to return the multiway valve to its locked position in time to prevent the bag from bursting! <li data-bbox="834 916 1477 983">➤ Close the clamp immediately, release the lock, and remove the bag.
 <p data-bbox="288 1534 362 1568">Fig. 7</p>	<p data-bbox="834 1088 1394 1122">Fig. 7 – Completing the sampling procedure</p> <ul style="list-style-type: none"> <li data-bbox="834 1133 1477 1200">➤ The sampling valve is now closed again by turning it clockwise. <li data-bbox="834 1211 1477 1312">➤ Disconnect the disposable parts after the multiway valve and close the bag immediately with the enclosed plug.
 <p data-bbox="411 1839 485 1872">Fig. 8</p>	<p data-bbox="834 1597 1430 1630">Fig. 8 – Preparing the bag for the shipping box</p> <ul style="list-style-type: none"> <li data-bbox="834 1641 1430 1675">➤ Gently squeeze the bag to check for any leaks. <li data-bbox="834 1686 1477 1753">➤ Affix the completed label on the bag and place it into the prepared shipping box. <li data-bbox="834 1765 1445 1832">➤ The bag must be delivered to the test laboratory within 24 hours.

4.5 Microbiological analysis at the dialysis water connection

The dialysis water connection on the media supply system serves as the sampling port.



4.5.1 Preparation







- Have a cooled shipping box available.
- The reverse osmosis system must be operated in the **RINSE** or **SUPPLY** mode for at least 20 minutes before the sample is collected.
- While collecting the sample, the reverse osmosis system must be in the **RINSE** or **SUPPLY** program.
- Disconnect the connection tube of the hemodialysis device from the dialysis water connection of the media supply.
- Collect the sample according to the procedure described for collecting a sample for microbiological analysis at the dialysis water connection (see chapter 4.5.3 on page 35).

4.5.2 Accessories, equipment

The following equipment is recommended by the manufacturer:

- Rubber gloves
- Alcohol-based disinfectant (approx. 70–80 % alcohol, without lipid replenisher)
- The **bag with adapter** (6030671) can be used for sampling the dialysis water.

4.5.3 Procedure for collecting a sample at the dialysis water connection

Illustration	Description
 <p>Fig. 1</p>  <p>Fig. 2</p>	<ul style="list-style-type: none"> ➤ Disinfect the dialysis water connection with an alcohol-based disinfectant (Fig. 1) and use a swab to wipe away any contamination (Fig. 2). ➤ Then repeat the disinfection procedure (Fig. 1 and 2). <p>Caution: Observe the acting time of the disinfectant!</p>
 <p>Fig. 3</p>  <p>Fig. 4</p>  <p>Fig. 5</p>  <p>Fig. 6</p>	<ul style="list-style-type: none"> ➤ The multiway valve on the sampling set (A) must be set so as to ensure that no fluid can flow (Fig. 3). ➤ The adapter of the sampling bag is placed on the dialysis water connection and locked (B) (Fig. 3). ➤ Then turn the multiway valve 90° clockwise (C) and "rinse" the dialysis water connection for approx. 60 seconds via the rinse tubing (Fig. 4). ➤ Now turn the multiway valve again 90° clockwise to fill the bag (Fig. 5). ➤ After approx. 250 ml (approx. half filled), promptly return the multiway valve to its original position (A) (Fig. 3) to prevent the bag from bursting. ➤ Close the clamp immediately, release the lock, and remove the bag. ➤ Disconnect the disposable parts after the multiway valve and immediately close the bag with the enclosed plug (Fig. 6). Gently squeeze the bag to check for any leaks. Affix the completed label on the bag and place it immediately into the prepared shipping box. The bag must be delivered to the test laboratory within 24 hours.

4.6 Collecting a sample for chemical analysis

4.6.1 Preparation

Dialysis water consumption is only possible if the reverse osmosis system is in **SUPPLY** mode, or if it is producing dialysis water during a manual rinse program in **RINSE** mode.

Before collecting the sample, the reverse osmosis system must have been in operation for at least 20 minutes. If the device is not in the **SUPPLY** mode, the manual rinse program must be started.

The sample is collected in the **SUPPLY** or **RINSE** mode.


4.6.2 Accessories, tools

The following accessories and tools are recommended by the manufacturer:



- Rubber gloves

For the chemical sampling, use the sample containers provided by the laboratory.

- For the sampling of the dialysis water, the **bag with adapter** (part no.: 6030671) can be used.



4.6.3 Procedure for collecting a sample for chemical analysis

Illustration	Description
 <p>Fig. 1 A</p>	<ul style="list-style-type: none"> ➤ The dialysis water must be collected while the system is in operation (during SUPPLY) or after extensive rinsing of the system (see above). Before collecting the sample, rinse the sampling valve (A or B) (approx. 2–10 L). <p>Caution: To avoid contamination of the sample by dirty sample containers, use only the containers supplied by the laboratory!</p> <ul style="list-style-type: none"> ➤ When collecting the sample, do not open the sample containers until directly before collecting the sample and make sure to close the containers immediately after taking the sample to prevent contamination. ➤ The sample containers must be filled up to the top. ➤ When collecting a sample the valve must first be sufficiently rinsed (approx. 5 L). Then the sample containers must be filled from the free flowing jet of water.
 <p>Fig. 2</p>	<ul style="list-style-type: none"> ➤ When collecting the sample on the media supply column, attach the bag with adapter to the coupling and fix it in place with the latch (Fig. 2 – B). ➤ Then thoroughly rinse the coupling (approx. 2 L) via the rinse tubing. ➤ Fill the sample container via the rinse tubing. <p>Caution: When collecting a sample at the coupling on the media supply column, do not use the bag with adapter as a sample container. Use the bottles supplied by the laboratory as sample containers!</p>

5 Alarms

This chapter does not contain any content applicable to this product.

6 Cleaning, disinfection

6.1 Generally applicable regulations for cleaning and disinfection



Warning**Risk of contamination from insufficient cleaning/disinfection**

There is a risk of spreading germs.

- The distribution system may only be cleaned and disinfected by persons who have been instructed on the proper handling of the system during these procedures.
 - The operator must observe and follow the general safety precautions.
 - The distribution system may only be disinfected after consultation with the manufacturer of the system or persons authorized by the manufacturer.
-

6.2 Precautions

6.2.1 Operator safety



Warning**Risk of chemical burns when working with acidic substances (disinfectant / cleaning agent)**

Contact with chemicals can cause chemical burns.

- Always handle acidic fluids with care and do not spill any disinfectant concentrate.
 - Wear appropriate personal protective equipment (gloves, goggles, etc.) in line with the safety precautions for the disinfectant / cleaning agent used.
 - Observe the safety precautions for the disinfectant / cleaning agent used, including the relevant first aid measures.
-



Warning

Risk of burning/scalding from hot surfaces or hot dialysis water during heat disinfection

Contact with hot surfaces or hot dialysis water can cause burning or scalding.

- Do not touch the accessible components of the dialysis water distribution system while heat disinfection is in progress.
 - Do not attempt to remove the dialysis water manually while heat disinfection is in progress.
-

6.3 Surface cleaning, surface disinfection

6.3.1 Surface cleaning

If the surface is contaminated by dust and dirt, clean the accessible surfaces of the **DWDL** and **DCDL**.



Note

Surface cleaning agents

The following procedures must be followed when using cleaning agents:

- If the **DWDL** and **DCDL** distribution systems are extremely dirty, wipe the affected sections with a damp cloth.
- Do not use any cleaning agents which contain acetone.
- Do not use solvents, diluting agents, or chemical cleaning sprays.
- Do not use any abrasives or aggressive cleaning agents and solvents.
- Do not use rough cleaning tools (e.g., scouring pad or similar) to clean the **DWDL** and **DCDL** distribution systems.

6.3.1.1 Cleaning agents for surfaces



Tip

The use of dialysis water is recommended for cleaning surfaces.

6.3.2 Surface disinfection



Note

The manufacturer recommends using **ClearSurf** for disinfecting the surface of the **DWDL** and **DCDL**.

- Proceed in accordance with the instructions of the manufacturer of the disinfectant to disinfect the surface.
- The manufacturer does not assume any liability for potential damage to the surfaces of the **DWDL** and **DCDL** if a different disinfectant than the one recommended is used for disinfection.

6.3.2.1 Surface disinfectant



Tip

The use of **ClearSurf** (1 % dilution) or **ClearSurf Wipes** (ready-to-use wipes) is recommended for disinfecting the surfaces.

For further information on disinfectants, refer to the chapter on consumables (see chapter 8.1 on page 54).

6.4 Disinfecting the DWDL

6.4.1 Reasons for disinfecting the DWDL

If a water supply as specified by the applicable regulations can no longer be ensured:

- After repairs to the dialysis water circuit.
- If the system has been idle for more than 72 consecutive hours.
- ISO 23500-1 "Guidance for the preparation and quality management of fluids for haemodialysis and related therapies" recommends a regular (e.g., monthly) preventive disinfection to avoid significant biofilm formation (biofouling).
- Regular disinfection according to the requirements of the responsible organization
- If a microbiological test reveals an elevated microbial count.

Recommended disinfectant

- **Puristeril 340**
- or:
- **Puristeril plus, Minncare®**

6.4.2 Disinfecting the system



Warning

Risk of contamination from insufficient cleaning/disinfection

There is a risk of spreading germs.

- The distribution system may only be cleaned and disinfected by persons who have been instructed on the proper handling of the system during these procedures.
- The operator must observe and follow the general safety precautions.
- The distribution system may only be disinfected after consultation with the manufacturer of the system or persons authorized by the manufacturer.



Note

The system may only be disinfected by trained Clinic Technicians or trained System Technicians who are trained and certified in the appropriate procedures.

- **Performing chemical disinfection**

The safety precautions, consumables, and operating steps – including performing the residual test – for chemical disinfection are described in the Service Manual of the reverse osmosis systems.

- **Performing heat disinfection**

The safety precautions and operating steps for heat disinfection are described in the Service Manual of the reverse osmosis systems.

6.5 Cleaning connectors and connection ports

The connectors and connection ports are cleaned to remove possible dialysis concentrate residues or salt crystals.

- Recommended equipment**
- Rubber gloves
 - Dialysis water
 - Small container
 - Soft, lint-free cloth or swab

6.5.1 Cleaning the connectors



- After disconnecting, immerse the connectors in a clean container filled with dialysis water and shake them gently.
- Then use a cloth or swab to pat dry the outer surfaces of the connectors.

6.5.2 Cleaning the connection ports

- Recommended equipment**
- Rubber gloves
 - Laboratory spray bottle with dialysis water
 - Soft, lint-free cloth or swab




- After disconnecting the couplings, use a laboratory spray bottle filled with dialysis water to rinse the connection ports installed on the media supply system.
- Then use a cloth or swab to carefully pat dry the connection ports.
Or: Use a damp cloth to carefully clean the connection ports and then pat them dry with a dry cloth or swab.
- After the connectors and connection ports have been cleaned, they must be disinfected (see chapter 6.6 on page 48).

6.6 Disinfecting connectors and connection ports

Recommended equipment

- Rubber gloves
- Alcohol-based disinfectant (approx. 70–80 % alcohol, without lipid replenisher)
- Soft, lint-free cloth or swab

6.6.1 Disinfecting the connectors



- After disconnecting, clean the connectors if necessary (see chapter 6.5 on page 47).
- Then wet the connectors with an alcohol-based disinfectant. After the dwell time, dry the outer surfaces of the connectors with a cloth or swab or allow the disinfectant to evaporate completely.



Note

To be observed during disconnecting:

- Do not place disconnected tubes onto the floor and do not connect with other tubes.
 - Hang the tubes on a suitable support with the opening pointing down.
-

6.6.2 Disinfecting the connection ports



- After disconnecting the connectors, use a laboratory spray bottle to rinse the connection ports installed on the media supply system and then carefully pat them dry with a cloth or swab.
- Then wet them with an alcohol-based disinfectant. After the dwell time, dry the outer surfaces of the connection ports with a cloth or swab or allow the disinfectant to evaporate completely.

7 Functional description

This chapter provides a brief functional description of the **DWDL** and **DCDL** distribution systems.

7.1 Description of procedures

7.1.1 Functions

The distribution systems **DWDL** and **DCDL** are connection units, including the tapping points for the dialysis water (**DWDL**) or dialysis concentrate (**DWDL**), between the reverse osmosis system (**DWDL**) or the concentrate supply system (**DCDL**) and the connected hemodialysis machines. The distribution systems **DWDL** and **DCDL** are used to transport dialysis water (**DWDL**) or dialysis concentrate (**DCDL**).

The trained personnel of the dialysis ward are considered users of the **DWDL** or **DCDL** distribution systems.

Installation and start-up are performed by trained technicians authorized by the manufacturer.

General functional description of the device

The **DWDL** is used to connect a reverse osmosis system with either a concentrate supply system or directly with a hemodialysis device in order to transport dialysis water. There is no direct contact with patients.

The **DCDL** distribution system is used to connect a concentrate supply system with a hemodialysis machine in order to transport dialysis concentrate. There is no direct contact with patients.

The **DWDL** or **DCDL** distribution systems are planned and installed according to the local conditions, regarding the length as well as the design and number of installation parts. The **DWDL** and **DCDL** distribution systems are to be regarded as permanently installed units.

8 Consumables, accessories, and additional equipment



Warning



Risks affecting the proper functioning of the device

The distribution system has been approved for use with certain consumables and accessories. Should the responsible organization wish to use other consumables and accessories than those listed in this chapter, the suitability must be checked beforehand by gathering the appropriate manufacturer information. The applicable legal regulations must be complied with.

The manufacturer does not assume any responsibility or liability for personal injury or other damage, and the use of non-approved or unsuitable consumables or accessories resulting in damage to the distribution system will void the warranty.




Upon request the local service support will provide information about further accessories, consumables, and other additional equipment.

8.1 Consumables

Part number	Consumables	Figure
5085851	Puristeril plus Active substance: Peracetic acid; D, GB, DK, E, FIN, I, NL, S	n/a
5085671	Puristeril 340 Active substance: Peracetic acid; D, GB, DK, E, FIN, I, NL, S	n/a
n/a	Minncare®	n/a
6030711	ClearSurf concentrate (or: ClearSurf wipes) Surface disinfectant	n/a
6299161	Peracetic acid test 5–50 mg/L	n/a
6030671	Bag with adapter Sampling set for standard configuration	
F00010382	Sampling set for the Fresenius sampling valve Sampling set for dialysis water distribution systems	



8.2 Accessories

8.2.1 Accessories for the DCDL distribution system

Part number	Description	Information
6309241	Tube PE white 8 x 12 mm white; length 100 m; Intended for concentrate type K1	Contact with dialysis concentrate 
M512671	Tube PE red 8 x 12 mm; length 100 m; Intended for concentrate type K2	Contact with dialysis concentrate 
M512681	Tube PE blue 8 x 12 mm; length 100 m; Intended for concentrate type K3	Contact with dialysis concentrate 
5453721	Reinforced tube, white D6 mm	Contact with dialysis concentrate
F40005702	Double nozzle D8 mm; PPSU	Contact with dialysis concentrate
n/a	O-ring EPDM 4 x 4 mm	Part of the MediaR, MediaC, MediaP
n/a	Nipple, closure conc. DN4, G1/4", PPSU	Part of the MediaR, MediaC, MediaP
n/a	Nipple, closure conc. DN4, G1/4", PVDF	Part of the MediaR, MediaC, MediaP
F00010505	Concentrate nipple PPSU; including seal	Contact with dialysis concentrate
F00010506	Concentrate nipple PVDF; including seal	Contact with dialysis concentrate
F40005755	Concentrate block K0 MediaR	Contact with dialysis concentrate
F00010855	Concentrate module MediaR ; including nozzles	Contact with dialysis concentrate

Part number	Description	Information
n/a	Concentrate block 0–0°	Part of the MediaC, MediaP
F00010501	Concentrate block 2 x straight MediaC, MediaP	Contact with dialysis concentrate
F00010502	Concentrate unit K1 double MediaC, MediaP ; modular system	Contact with dialysis concentrate
F00010503	Concentrate unit K2 double MediaC, MediaP ; modular system	Contact with dialysis concentrate
F00010504	Concentrate unit K3 double MediaC, MediaP ; modular system	Contact with dialysis concentrate
F40005704	Concentrate block 45–0°	Contact with dialysis concentrate
F40005705	Concentrate block 0–45°	Contact with dialysis concentrate
M026391	Tube clamp Single; 14.5 mm	No contact with dialysis concentrate

8.2.2 Accessories for the DWDL distribution system

Part number	Description	Information
6345031	XLPE tube, natural 25 x 3.5 mm (100 m distribution)	Contact with dialysis water 
6309351	XLPE tube, natural 25 x 3.5 mm (50 m distribution)	Contact with dialysis water 
6316031	Connector 90° Stainless steel for PEX tube 25 x 3.5 mm	Contact with dialysis water
6316041	Connector, straight Stainless steel for PEX tube 25 x 3.5 mm	Contact with dialysis water
6325801	PEX connector 180° U-connector for stainless steel cable duct	Contact with dialysis water
n/a	Dialysis water connection, straight MSM PEX double with connecting element	Part of the MediaC, MediaP

Part number	Description	Information
n/a	Dialysis water connection, double, straight MSM PEX double with connecting element	Part of the MediaC, MediaP
F00007306	Dialysis water supply center Double arch MediaC , single	Contact with dialysis water
F00007307	Dialysis water supply center Double arch MediaC , double	Contact with dialysis water
F00006458	Fluid Fly Loop Secondary ring main 2 x 2 m, strain relief, ring nut, device connection, nipple with sample collection	Contact with dialysis water
F00006459	Fluid Fly Loop Secondary ring main 2 x 2 m, strain relief, ring nut, device connection, nipple without sample collection	Contact with dialysis water
n/a	Flat gasket EPDM; 11.5 x 5 x 3.6	Part of the MediaR, MediaC, MediaP
n/a	Nipple closure NW 6-G 1/4" 11 mm dialysis water	Part of the MediaC, MediaP
n/a	Nipple closure NW 6-R 1/4" 11 mm dialysis water	Part of the MediaR
F00010499	Dialysis water nipple Walther incl. seal	Contact with dialysis water
F00010500	Dialysis water nipple FIDICA incl. seal	Contact with dialysis water
F00010492	Dialysis water block MediaC, MediaP single, straight	Contact with dialysis water
F00010493	Dialysis water block MediaC, MediaP double, straight	Contact with dialysis water
F00010494	Dialysis water block stainless steel MediaC, MediaP single – 2 x 90° angled – modular system	Contact with dialysis water
F00010495	Dialysis water block stainless steel MediaC, MediaP double – 2 x 90° angled – modular system	Contact with dialysis water
F00010826	Dialysis water module MediaR	Contact with dialysis water
F00010873	PEX connection, long (set), Dialysis water block MediaR – modular system	Contact with dialysis water
6309401	Sleeve 25 x 3.5 mm	No contact with dialysis water

8.3 Additional equipment

8.3.1 Additional equipment for the DWDL distribution system

Part number	Description	Information
F40001179	Sampling set PE-Xa	Contact with dialysis water
F00010866	Lock washer 19 Dialysis water MediaC , MediaP , with fixing	No contact with dialysis water
F00010507	Lock washer 19 Dialysis water/waste water MediaR , with fixing	No contact with dialysis water

8.3.2 Additional equipment for the DCDL distribution system

Part number	Description	Information
F00010869	Lock washer 17 Dialysis concentrate MediaC , MediaP , with fixing	No contact with dialysis concentrate
F00010508	Lock washer 17 Dialysis concentrate MediaR , with fixing	No contact with dialysis concentrate

9 Installation

9.1 Installation requirements

9.1.1 General information

Follow the applicable installation guidelines

For new installations, the applicable installation guidelines must be followed.

To be observed before the operational qualification

The operational qualification of the distribution systems is performed together with the operational qualification of the supplying devices.

Observe national and local regulations

The national or local installation, operation, use, and maintenance regulations must be complied with.

State of components

Prior to installation, check the components of the distribution systems for damage that may have occurred during transport.



Note

Protecting components containing fluid

- Prior to operational qualification, check the components of the distribution system for damage that may have occurred during transport. If there are signs of damage to the hydraulic components, do not use the distribution system.
-

Access to the connection points

The connection points of the distribution systems must be accessible for checks.

Performance parameters of connected devices

The performance parameters of the connected device must be taken into account when designing the distribution system. For further information, see the Service Manuals of the supplying devices.

9.1.2 Environmental conditions

Observe the local conditions

- The installation site must be free from frost and dust.
- The components must not be exposed to continuous, direct sunlight.

9.1.3 Structural installation requirements

Leakage sensor

The use of a leakage sensor is recommended.

Pilot holes

The hole pattern with hole sizes and distances must be applied. For further information, the Installation Guideline (IGL) must be followed.

9.2 Operational qualification

9.2.1 To be observed before the operational qualification

Tester's qualification	Operational Qualification must be performed by the technical service department of Fresenius Medical Care or a person authorized by them. The Operational Qualification may only be performed by persons qualified to properly perform the specified checks based on their educational background, training, knowledge, and experience. Furthermore, the persons performing the checks must not be bound by any directives when performing this activity.
For operational qualification only	The following information is only intended for the operational qualification. It is not applicable to the operational requalification of DWDL and DCDL distribution systems that have been decommissioned or temporarily decommissioned.
Specifications	<ul style="list-style-type: none"> – Observe the information on the specifications. – Specific connection and performance data must be taken from the Specifications chapter.
Use of spare parts	Any installation, modification, or repair work requiring the DWDL and DCDL to be opened may only be performed by manufacturer-authorized persons and is permitted only when using original spare parts.
Identification of supply ports	Identification of the distribution system supply ports: <ul style="list-style-type: none"> – Dialysis water – Dialysis concentrate (C1, C2, C3)
Test equipment and tools	The activities described in this document require the availability of the necessary technical test equipment and tools.
Maintenance procedures (MA)	For further information (see chapter 11.2 on page 66).



Note

All local regulations regarding technical safety must be observed.

Precautions Repair any visible damage.

9.2.2 Operational Qualification procedure



Note

➤ When performing the operational qualification on the distribution systems, the descriptions in the Service Manual must be followed.

● After Operational Qualification

Warning

Risk of contamination from insufficient cleaning/disinfection

Before start-up, a chemical disinfection must be performed on the **DWDL**. The successful disinfection must be verified by means of a microbiological analysis.



Note

- The senior physician must be informed about the results of the microbiological analysis. The Technical Safety Checks must be performed and reported.
-

9.3 Removal from service, decommissioning



Note

- For information regarding the removal from service or decommissioning of the **DWDL** and **DCDL** distribution system, contact the local service support.
-

9.3.1 Removal from service



Note

If the **DWDL** and **DCDL** distribution systems are removed from service after the operational qualification, the following has to be observed:

- On operational requalification, the water supply pressure must be checked against the prescribed minimum pressure.
-

9.3.2 Decommissioning



Note

- For information regarding decommissioning of the **DWDL** and **DCDL** distribution systems, contact the local service support.
-

10 Transport and storage



Note

The following transport and storage conditions and further information regarding transport and storage affect the **DWDL** and **DCDL** distribution systems.

10.1 Transport and storage conditions

Storage temperature range 5 to 45 °C

**Note**

Protect the **DWDL** and **DCDL** distribution systems from frost.

Relative humidity 20 to 80 % at 20 °C (non-condensing)

Atmospheric pressure 700 hPa to 1150 hPa

**Note****Protection from exposure to UV light**

Do not expose the components of the distribution system to direct sunlight (UV rays may cause premature aging of the materials).

➤ Do not store outdoors!

10.2 Environmental compatibility/disposal

In the EU member states, the **DWDL** and **DCDL** can be returned to the manufacturer. Please also observe the applicable local legal regulations.

Before the device is sent off for disposal, the responsible organization must ensure that all consumables attached to the **DWDL** and **DCDL** are removed and the **DWDL** and **DCDL** are disinfected as specified by the manufacturer (see chapter 6 on page 41).

The responsible organization must also inform the disposal company responsible for dismantling and disposing of the **DWDL** and **DCDL** of the following before the start of the disposal measures:

- It is possible that the **DWDL** and **DCDL** distribution systems may be contaminated when returned. Therefore, it is vital to take suitable precautions when dismantling it, such as wearing personal protective equipment.
- Further information will be provided if requested by disposal companies.

10.2.1 Handling of disinfectants

It is absolutely essential to observe the manufacturer's specifications of the disinfectants used (relating to protective clothing, storage, dosing, expiration date).

The local regulations with regard to the disposal of waste water, if defined, must be clarified and observed before using disinfectants.

11 Technical Safety Checks and maintenance

11.1 Important information for performing the TSC/MA

Checks	The Technical Safety Checks (TSC) must be performed every 24 months.
Tester's qualification	<p>The checks must be performed by the manufacturer's service support or a manufacturer-authorized person.</p> <p>The checks may only be performed by persons qualified to properly perform the specified checks owing to their educational background, training, knowledge, and experience. Furthermore, the persons performing the checks must not be bound by any directives when performing this activity.</p>
Specifications	Observe the information on the specifications.
Documentation	<ul style="list-style-type: none">➤ To perform the Technical Safety Checks (TSC) and the maintenance procedures (MA), contact the local service support.➤ Performance of the Technical Safety Checks must be entered in the Medical Device Register.

11.2 Maintenance procedures

The following procedures must be performed by the operator in accordance with the specifications regarding intervals.

11.2.1 Quality assurance and care measures

Action	Accessories/target state	Interval	Comment
Cleaning and disinfection measures	<ul style="list-style-type: none"> ➤ Surface cleaning or surface disinfection ➤ Cleaning or disinfecting the connectors and connection ports 	Recommendation: When required	(see chapter 6.3 on page 43), (see chapter 6.5 on page 47), (see chapter 6.6 on page 48)
Visual inspection and leakage check	<ul style="list-style-type: none"> ➤ Perform a visual inspection of all connectors and tubing containing fluid. 	Daily	(see chapter 4.3.1 on page 28)


12 Specifications

12.1 DWDL and DCDL product data

Specifications	DWDL	DCDL
Dimensions Internal diameter (mm) Wall thickness (mm) External diameter (mm)	25 x 3.5 mm 18.0 3.5 25.0	12 x 2 mm 8.0 2.0 12.0
Weight (kg)	depending on the project planning	
Material	Polyethylene, cross-linked	Low-density polyethylene (LDPE)
Medium	Dialysis water according to ISO 23500-3	Dialysis concentrate according to ISO 23500-4: – Acetate-based, acidic dialysis concentrates – Citrate-based, acidic dialysis concentrates
Operating conditions		
Operating pressure	0–6 bar	
Operating temperature range	5–35 °C	5–30 °C
Atmospheric pressure	700 hPa to 1150 hPa	
Relative humidity	20 to 80 % at 20 °C (non-condensing)	
Medium operating temperature	5–35 °C	5–30 °C
Maximum temperature during heat disinfection	95 °C	--
Transport and storage conditions		
Storage temperature range	5 to 45 °C	
Atmospheric pressure	700 hPa to 1150 hPa	
Relative humidity	20 to 80 % at 20 °C (non-condensing)	
ISO standards	ISO 23500-1 ISO 23500-2 ISO 23500-3	ISO 23500-1 ISO 23500-4
Materials used	According to ISO 10993-1	
The DWDL and DCDL distribution systems are planned and installed according to the local conditions regarding the length as well as the design and number of installation parts. The DWDL and DCDL distribution systems are to be regarded as permanently installed units.		

12.2 Identification label (DWDL and DCDL identification)

The identification label shown is only an example. The actual data is the data specified on the identification label of the **DWDL** and the **DCDL**.



- 1 Type identification
- 2 Serial number
- 3 Barcode, Code 39
- 4 Max. inlet pressure
- 5 Manufacturer: date of manufacture and manufacturer's address
- 6 (11) Date of manufacture YYMMDD, 6 digits
- 7 (21) Serial number, 8 digits
- 8 (01) *GTIN (**SAP: EAN/UPC code), 13 digits plus digit 0
- 9 ***UDI scan code
- 10 UDI identification
- 11 Medical device identification
- 12 CE mark
- 13 Part number and identification label edition
- 14 REF = SAP material number

*GTIN = **G**lobal **T**rade **I**tem **N**umber

SAP: EAN/UPC code = SAP product code: **European **A**rticle **N**umber/**U**niversal **P**roduct **C**ode

***UDI = **U**nique **D**evice **I**dentification

12.3 Transport/storage

For further information, (see chapter 10 on page 63).

12.4 Materials used

For further information, (see chapter 12.1 on page 67).

13 Definitions



















13.1 Definitions and terms

Coupling system	Hydraulic, mechanical connections establish the connection between the hemodialysis system and the dialysis media supply system.
Dialysate, dialysis fluid	The exchange fluid used in dialysis.
Dialysis concentrate	Dialysis concentrate is a highly concentrated solution consisting of solid and/or liquid components and dialysis water. Dialysis concentrate is used in the hemodialysis device with other components for the production of dialysis fluid.
Dialysis concentrate distribution system	Transport line which provides dialysis concentrate for use in dialysis units.
Dialysis water	Water suitable for dialysis treatments (water processed by a reverse osmosis system which meets the requirements of ISO 23500-3). A high-pressure pump, membrane module, and appropriate monitoring equipment are used to produce dialysis water from drinking water.
Dialysis water distribution system	Transport line which provides dialysis water for use in dialysis units.
Initial operational qualification	Initial start-up
Media supply system	Main interface and connection unit between the hydraulic lines such as dialysis water supply line, dialysis concentrate supply line, drainage line, and the hemodialysis system.
Operational qualification	Start-up
Operational requalification	Recommissioning
Permeate	This term is used as a synonym for dialysis water. This term must only be used in a technical context.

13.2 Abbreviations

C1	Dialysis concentrate 1
C2	Dialysis concentrate 2
C3	Dialysis concentrate 3
DCDL	Dialysis Concentrate Distribution Loop
DWDL	Dialysis Water Distribution Loop
Fig.	Figure (diagram)
MA	Maintenance
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RO	Reverse osmosis system
SVHC	Substance of Very High Concern
TSC	Technical Safety Checks

13.3 Symbol

Symbols	Description
	The CE mark documents compliance with the MDR; notified body: TÜV PRODUCT SERVICE 0123
	Year/month/day of manufacture
	Serial number
	Medical device
	Material number
	Unique Device Identification
	Follow the Instructions for Use!
	eIFU, www.freseniusmedicalcare.com/en/product-information-fme
	Caution! Do not stack!
	Do not use hooks!
	Store upright! Do not tilt!
	Caution: Fragile
	Keep dry!
	Warning: Hot surface
	Permissible temperature range
	Atmospheric pressure operating conditions range
	Relative air humidity operating conditions range
	Protect against sunlight (UV light)!

13.4 Certificates

Upon request the local service support will provide the currently valid versions of the certificates.

14 Options

This chapter does not contain any content applicable to this product.

15 Appendix

15.1 Medical Device Register for the DWDL and DCDL

15.1.1 Responsible organization and identification

The following page shows the master copy of the address of the responsible organization and the product identification.

Dialysis Water Distribution Loop, Dialysis Concentrate Distribution Loop	Address of the responsible organization & identification	 FRESENIUS MEDICAL CARE
--	---	---

Address of the responsible organization
Name: _____
Address: _____
City: _____
Phone: _____
Site of installation _____

Internal medical consultant
Name, phone: _____
Name, phone: _____
Name, phone: _____
Name, phone: _____
Name, phone: _____

Identification
DWDL, DCDL
Type: Dialysis water and dialysis concentrate
Classification: IIa
Registration number: _____
Identification number of the certified body: 0123
Serial number: _____
Flexible insulation available; yes <input type="checkbox"/> , no <input type="checkbox"/>
Manufacturer: Fresenius Medical Care & Co. KGaA, 61352 Bad Homburg

Tests and checks	
Type	Intervals
Technical Safety Check (TSC)	Every 24 months
_____	Every _____ months
_____	Every _____ months

Contracts regarding tests and checks:
Technical Safety Checks:
Company name: _____
Address: _____
Phone: _____

15.1.2 Contents of the Medical Device Register for the DWDL and DCDL

The following page shows the contents of the Medical Device Register for the **DWDL** and **DCDL** distribution systems.

Dialysis Water Distribution
Loop, Dialysis Concentrate
Distribution Loop

Contents of the Medical Device Register




1	Instructions for Use
Monitoring	
2	Microbiological and chemical monitoring for DWDL – Results of the microbiological analysis – Results of the chemical analysis – Sampling plans
3	Settings reports
4	Service reports, device training, malfunctions – Device Training Record – Service reports and documentation of modifications to the device equipment – Reporting of incidences – Documentation of malfunctions and repeated, identical operating errors
5	Technical Safety Checks (TSC) and revalidation
Validation phase	
6	Installation qualification (IQ) – Installation report – Validation plan
7	Operational qualification (OQ) – Training record OQ – Sampling plan OQ – Disinfection plan OQ – Start-up report OQ
8	Performance qualification (PQ) – Results of the microbiological analysis PQ – Results of the chemical analysis PQ

15.2 Training Record

Applicability, purpose	The responsible organization must ensure that operators receive proper instruction. Operator instruction is based on the Instructions for Use and, if present, the enclosed Additional Sheets. The manufacturer recommends using this Training Record to document the operator instruction performed.
Significance of the warnings	Observing all warnings in the Instructions for Use is essential for using the device safely. Instruction must be given on all warnings in the Instructions for Use.

● Explanations on the Training Record report

General information	<ul style="list-style-type: none"> – The report heading records the circumstances of operator instruction. – The report footer records the trainer and participants. – The chapters of the Instructions for Use are listed on separate lines up to the second level.
Y/N/NA	<ul style="list-style-type: none"> – <input type="checkbox"/>/–/– Chapter required for proper operator instruction. – <input type="checkbox"/>/□/– Chapter recommended for proper operator instruction. – <input type="checkbox"/>/–/□ If option is available: Chapter required for proper operator instruction. – <input type="checkbox"/>/□/□ If option is available: Chapter recommended for proper operator instruction. <p>➤ Record the instruction of the relevant content and warnings completed by marking ✓ in field Y.</p> <p>➤ Record chapters or options that have not been instructed with ✓ in field N.</p> <p>➤ Record unavailable options with ✓ in field NA.</p>

 FRESENIUS MEDICAL CARE		Training Record	DWDL and DCDL
Customer name:		Start date:	
Address:		End date:	
Customer name:		Software version: n/a	
Serial number DWDL 1:			
Serial number DWDL 2:			
Serial number DWDL 3:			
Serial number DCDL 1:			
Serial number DCDL 2:			
Serial number DCDL 3:			
Description			Y/N/NA
1	Index		
2	Important information		
2.1	How to use the Instructions for Use	<input type="checkbox"/> /□/–	
2.2	Warnings, significance	<input type="checkbox"/> /–/–	
2.3	Notes, significance	<input type="checkbox"/> /□/–	
2.4	Tips, significance	<input type="checkbox"/> /□/–	

Description	Y/N/NA
2.5 Brief description	<input type="checkbox"/> /0/-
2.6 Intended purpose and related definitions	<input type="checkbox"/> /0/-
2.7 Side effects	<input type="checkbox"/> /0/-
2.8 Contraindications	<input type="checkbox"/> /0/-
2.9 Residual risks	<input type="checkbox"/> /0/-
2.10 Interaction with other systems	<input type="checkbox"/> /0/-
2.11 Therapy restrictions	<input type="checkbox"/> /0/-
2.12 Considerations for working on the DWDL and DCDL	<input type="checkbox"/> /-/-
2.13 Expected service life	<input type="checkbox"/> /0/-
2.14 Duties of the responsible organization	<input type="checkbox"/> /0/-
2.15 Operator responsibility	<input type="checkbox"/> /-/-
2.16 Disclaimer of liability	<input type="checkbox"/> /-/-
2.17 Technical documentation	<input type="checkbox"/> /0/-
2.18 Warnings	<input type="checkbox"/> /-/-
2.19 SVHC (REACH)	<input type="checkbox"/> /0/-
2.20 Addresses	<input type="checkbox"/> /0/-
3 Structure and views	
3.1 Overall view of the distribution systems	<input type="checkbox"/> /0/-
3.2 Overall view of DWDL	<input type="checkbox"/> /0/-
3.3 Overall view of DCDL	<input type="checkbox"/> /0/-
4 Operation	
4.1 Overview of coupling systems on the media supply systems	<input type="checkbox"/> /0/-
4.2 Observe when operating	<input type="checkbox"/> /0/-
4.3 Operating the coupling systems	<input type="checkbox"/> /-/-
4.4 Microbiological analysis at the sampling port	<input type="checkbox"/> /0/-
4.5 Microbiological analysis at the dialysis water connection	<input type="checkbox"/> /0/-
4.6 Collecting a sample for chemical analysis	<input type="checkbox"/> /0/-
5 Alarms	
6 Cleaning, disinfection	
6.1 Generally applicable regulations for cleaning and disinfection	<input type="checkbox"/> /-/-
6.2 Precautions	<input type="checkbox"/> /-/-
6.3 Surface cleaning, surface disinfection	<input type="checkbox"/> /0/-
6.4 Disinfecting the DWDL	<input type="checkbox"/> /-/-
6.5 Cleaning connectors and connection ports	<input type="checkbox"/> /0/-
6.6 Disinfecting connectors and connection ports	<input type="checkbox"/> /0/-
7 Functional description	
7.1 Description of procedures	<input type="checkbox"/> /0/-
8 Consumables, accessories, and additional equipment	
8.1 Consumables	<input type="checkbox"/> /0/-
8.2 Accessories	<input type="checkbox"/> /0/-

Description	Y/N/NA
8.3 Additional equipment	<input type="checkbox"/> / <input type="checkbox"/> /-
9 Installation	
9.1 Installation requirements	<input type="checkbox"/> / <input type="checkbox"/> /-
9.2 Operational qualification	<input type="checkbox"/> /-/-
9.3 Removal from service, decommissioning	<input type="checkbox"/> / <input type="checkbox"/> /-
10 Transport and storage	
10.1 Transport and storage conditions	<input type="checkbox"/> / <input type="checkbox"/> /-
10.2 Environmental compatibility/disposal	<input type="checkbox"/> / <input type="checkbox"/> /-
11 Technical Safety Checks and maintenance	
11.1 Important information for performing the TSC/MA	<input type="checkbox"/> / <input type="checkbox"/> /-
11.2 Maintenance procedures	<input type="checkbox"/> / <input type="checkbox"/> /-
12 Specifications	
12.1 DWDL and DCDL product data	<input type="checkbox"/> / <input type="checkbox"/> /-
12.2 Identification label (DWDL and DCDL identification)	<input type="checkbox"/> / <input type="checkbox"/> /-
12.3 Transport/storage	<input type="checkbox"/> / <input type="checkbox"/> /-
12.4 Materials used	<input type="checkbox"/> / <input type="checkbox"/> /-
13 Definitions	
13.1 Definitions and terms	<input type="checkbox"/> / <input type="checkbox"/> /-
13.2 Abbreviations	<input type="checkbox"/> / <input type="checkbox"/> /-
13.3 Symbol	<input type="checkbox"/> / <input type="checkbox"/> /-
13.4 Certificates	<input type="checkbox"/> / <input type="checkbox"/> /-
14 Options	
15 Appendix	
15.1 Medical Device Register for the DWDL and DCDL	<input type="checkbox"/> / <input type="checkbox"/> /-
15.2 Training Record	<input type="checkbox"/> / <input type="checkbox"/> /-
15.3 Quality of dialysis water	<input type="checkbox"/> / <input type="checkbox"/> /-
Comments:	

**Note**

- Observe index, important information, and all warnings in the Instructions for Use!

Trainer			
Date	Name		Signature
Participant			
Date	Function	Name	Signature

Trainer			
Date		Name	Signature

15.3 Quality of dialysis water

The microbiological and chemical purity of the dialysis fluid prepared in the dialysis clinic is of critical importance for the quality of the patient's treatment. The quality of the dialysis water should comply with local regulations. If no local regulations apply, compliance with the applicable requirements of ISO 23500-3 "Water for haemodialysis and related therapies" is necessary.

The quality of the dialysis water should be monitored regularly for listed chemical and microbiological contaminants. The monitoring schedule should be based on the results of the system validation. In an existing water treatment system operated under stable conditions, the chemical contaminants in the dialysis water should be monitored at least once every year. This excludes total chlorine which, if present in the feed water, should be monitored at the beginning of each treatment day.

Compliance with requirements for the chemical parameters according to ISO 23500-3 may necessitate additional water pretreatment stages or a change in the yield on the device. The composition of the dialysis water must be checked as part of the performance qualification (PQ), and the water pretreatment must be adjusted as needed.

● Microbiological quality of fluids for hemodialysis

Reference	Medium	Permissible maximum values	
		Total viable microbial count [CFU/ml]	Endotoxin concentration [EU/ml]
ISO 23500-3 Water for haemodialysis and related therapies	Dialysis water	< 100 (AL* 50)	< 0.25 (AL* 0,125)
ISO 23500-5 Quality of dialysis fluid for haemodialysis and related therapies	(Standard) dialysis fluid **	< 100 (AL* 50)	< 0.5 (AL* 0.25) (Ph.Eur: < 0.25)

*AL = Action Level. Starting at this concentration, steps must be taken to stop the trend from reaching higher, unacceptable values. This value is typically 50 % of the maximum allowable level.

**Tests for bacterial growth and endotoxins are not required if the dialysis device's fluid pathway is fitted with a bacteria-retentive and endotoxin-retentive filter that has an appropriate capacity, has been validated by the manufacturer, and is operated and monitored according to the manufacturer's instructions (e.g., DIASAFE plus).

● **Chemical quality of dialysis water**

ISO 23500-3					
Contaminants with proven toxicity in dialysis	Maximum allowable level [mg/L]	Electrolytes	Maximum allowable level [mg/L]	Trace elements	Maximum allowable level [mg/L]
Aluminum	0.01	Calcium	2	Antimony	0.006
Lead	0.005	Potassium	8 (*2)	Arsenic	0.005
Fluoride	0.2	Magnesium	4 (*2)	Barium	0.1
Total chlorine	0.1	Sodium	70 (*50)	Beryllium	0.0004
Copper	0.1			Cadmium	0.001
Nitrate as (N)*	2			Chrome	0.014
Sulfate	100 (*50)			Mercury	0.0002 (*0.001)
Zinc	0.1			Selenium	0.09
				Silver	0.005
				Thallium	0.002

* Values according to the European Pharmacopoeia (Ph. Eur.); applicable regulations must be observed. Other deviations in the Ph.Eur. are: Nitrate: Alarm limit = 2 mg/L nitrate in proportion to the total nitrate molecule NO₃. Other harmful substances that are only listed in the Ph.Eur.: Ammonium (NH₄): 0.2 mg/L; heavy metals (e.g., Pb): 0.1 mg/L; chloride: 50 mg/L.

For continued compliance with quality standards, checks and disinfections of the dialysis water system must be performed regularly.

Recommended chemical surveillance

Annual inspection	The dialysis water should be checked for chemical contamination at least once a year.
Offline tests	If the feed water or pretreated water is chlorinated and offline tests are used, the total chlorine test should be performed downstream of the activated carbon filter at the beginning of each treatment day, before the first patient treatment. If chloramine is used at a concentration of 1 mg/L or more to disinfect the drinking water supply, the test should be repeated before the start of each patient session. If no patient sessions are scheduled, the test should be carried out approximately every 4 hours during operation.
Online tests	For online tests in the water pretreatment system, the chlorine and total hardness parameters, for example, can be monitored using AquaSENS .

