DALI 500 / 750 Adsorber

de	Gebrauchsanweisung
en	Instructions for Use
bg	Инструкции за работа
bs	Upute za uporabu
cs	Návod k použití
da	Brugervejledning
el	Οδηγίες χρήσεως
es	Instrucciones de uso
et	Kasutusjuhend
fi	Käyttöohje
fr	Mode d'emploi
hr	Upute za uporabu
hu	Használati utasítás
it	Istruzioni per uso
lt	Naudojimo instrukcijos
lv	Lietošanas instrukcija
nl	Gebruiksaanwijzing
no	Bruksanvisning
pl	Instrukcja użycia
pt	Instruções de uso
ro	Instrucțiuni de utilizare
ru	Упаковочный вкладыш
sk	Návod na použitie
sl	Navodila za uporabo
sr	Uputstvo za upotrebu
sv	Bruksanvisning
tr	Kullanma talimatları

de DALI 500 / 750 ADSORBER

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INSTRUCTIONS FOR USE

Please read the following instructions carefully.

⚠️ Caution

Before using the DALI 500 / 750 Adsorber, read also the Instructions for Use (IFU) of the other medical devices dedicated for DALI treatment (DALI Kit components, DALI Priming Solution and suitable active device (hemoadsorption system)). Before treating the patient, the physician must make an individual benefit-risk analysis.

GENERAL INFORMATION

General description of the product

The DALI Adsorbers are sterile and intended for single-use only for lipoprotein removal during whole blood apheresis. The DALI Adsorbers must be used together with the other components from the DALI Kit, the DALI Priming Solution and with an active device (hemoadsorption system). The DALI Adsorbers consist of a polycarbonate housing filled with adsorbent. The adsorbent is made of porous polyacrylamide beads specified with diameters of 150 – 230 µm coated with negatively charged polyacrylic acid. This serves as a ligand that binds positively charged apoprotein B of Low-Density Lipoprotein (LDL), Lipoprotein(a) (Lp(a)) and Very Low-Density Lipoprotein (VLDL) based on charge interaction. The extracorporeal volume of the adsorbers is approximately 285 mL for DALI 500 and 381 mL for DALI 750.

Sterilization method: INLINE steam

Be aware, since steam is used in the sterilization process, water drops may be formed in the adsorbers.

COMPOSITION

Model	DALI 500	DALI 750
Catalogue number	9795301	9795311
ADSORBENT		
Matrix	Methacrylate	
Ligand	Polyacrylic acid	
HOUSING		
Corpus, funnel and header	Polycarbonate	
Sieves	Polyester	
Membrane	Acrylic Copolymer	

INTENDED PURPOSE AND RELATED DEFINITIONS

Intended purpose

Removal of specific substances from blood in an extracorporeal treatment.

Medical indication

Hypercholesterolemia and/or isolated increased Lp(a) leading to cardiovascular and/or peripheral vascular diseases.

Clinical benefits

Lipoprotein apheresis leads immediately to a high reduction of LDL-c and Lp(a) with additional benefits to rheology by pleiotropic effects through reduction of further atherogenic factors, e.g. oxidized low-density lipoprotein (oxLDL). The effective reduction of lipoproteins with lipoprotein apheresis leads to a reduction of the incidence rate of cardiovascular events. Especially in patients with homozygous familiar hypercholesterolemia, a reduction of their xanthelasmas and xanthomas can be observed over time. In conclusion, the clinical benefit has been evaluated to outweigh the risks of use.

Intended patient population

The DALI 500 / 750 Adsorber is intended to be used in adult patients with a minimum weight of 40 kg.

There is no data for the use of DALI in children (not even > 40 kg), adult patients < 40 kg, pregnant patients and nursing mothers. Therefore, we cannot recommend treating such patients.

The choice of the selected system configuration and the blood volume to be treated depend on the patient's clinical situation and on the prescription of the treating physician. For the selection of the system configuration see also the table regarding recommended treatment parameters.

Intended user group and intended use environment

The DALI Adsorbers must only be used by health care professionals who are familiar with aseptic techniques and who have been trained in the proper operation and handling of the medical device. The DALI system including the DALI Adsorbers shall only be used in professional healthcare institutions under medical supervision. A monitoring of the patient's vital parameters during treatment by the attending staff has to occur. Equipment and drugs for emergency treatment of hypotension/shock or anaphylactic/anaphylactoid reaction must be available.

SIDE EFFECTS

Extracorporeal treatment modality itself places a burden on the patient and can result in side effects which may limit the performance of the treatment and may continue for some time after the treatment. This stress varies from patient to patient and may result from factors specific to the patient's condition, operating parameters, equipment, anticoagulation, medication, etc. During apheresis treatment, the blood and plasma comes in contact with the extracorporeal "foreign" surfaces which may result in effects on the blood components, such as activation of the coagulation and complement system, and could potentially trigger an anaphylactic reaction.

Physicians have a duty to inform the patient about the possibility of side effects and that the patient's capacity may be limited for some time after the treatment.

Organ class (IMDRF****)	Side effects
Nervous system	Dizziness, speech disorders, loss of consciousness, movement disorders** (shivering, paralytic symptoms), numbness**
Blood and lymphatic system	Anemia, leucopenia, pancytopenia, thrombocytopenia, coagulation disorders (binding of coagulation factors from the intrinsic pathway), hypovolemia (through blood loss), hypervolemia and oedema, hemolysis (symptoms see hypersensitivity reactions), electrolyte imbalance (especially hypocalcemia**)
Immune system	Hypersensitivity/allergic reactions (leading to, e.g. flush, hypotension, nausea, headache, back and abdominal pain, local reaction (pruritus, erythema, eczemas), urticaria, dyspnea, angioedema (mouth area, larynx, pharynx), loss of consciousness, shock), DALI-induced thrombocytopenia, heparin-induced thrombocytopenia (HIT)*
Vascular system	Bleeding/blood loss (also: access-related), embolism (air-, foreign body-, thrombo-embolism), phlebitis*, cardiovascular insufficiency, hematoma*, thrombosis*
Heart	Arrhythmia (bradycardia, tachycardia), angina pectoris
Gastrointestinal system	Nausea, vomiting, abdominal pain, diarrhea
Metabolism	Protein loss, iron loss, hypoglycemia/hyperglycemia*
Musculoskeletal system	Tingling**, cramps/muscle spasms**
Other/generalized disorders	Fatigue, discomfort, hypothermia, hyperthermia (heat or cold sensation), pain (abdominal, chest, eye, teeth, extremities, headache), hypotension, hypertension, blurred vision, pallor, chills, sweating

* Treatment related ** Anticoagulation (Citrate) related *** International Medical Device Regulators Forum Additional side effects might be specific for other devices or drugs used in the therapy, e.g. catheters.

Reporting of serious incidents

A serious incident means any incident that directly or indirectly leads to any of the following: the death of a patient, user or other person; the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; a serious public health threat. If any serious incident occurs in relation to the device, inform the treating physician immediately.

Report this serious incident to the manufacturer according to labelling (☞) and if required to the corresponding local authority. Within the EU the user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

Medical warnings/precautions (to prevent side effects)

- To avoid hypersensitivity or hypersensitivity-like reactions, a slow start into the DALI therapy in the first treatments is necessary (please refer to chapter "Method of administration/Recommended treatment parameters") to adapt the patient on the extracorporeal treatment. These reactions usually disappear after the adaptation phase. Appropriate precautionary measures (emergency treatment equipment) must be available. If necessary, the treatment must be discontinued without returning the blood of the extracorporeal circuit to the patient and medical treatment steps initiated immediately.
- To ensure a safe and effective treatment, comply with the variables treated blood volume, blood flow and Acid-Citrate-Dextrose A (ACD-A) flow as summarized in the recommendation in chapter "Method of Administration/Recommended treatment parameters".
- To prevent any severe change in blood cells, take complete blood counts regularly. Acute or chronic thrombocytopenia during DALI treatment have been reported and therefore close controls of the blood count (especially from thrombocytes) are necessary, especially when the therapy is intensified, e.g. when the therapy frequency or the adsorber configuration is increased. The DALI apheresis procedure may result in temporary mild leucocyte and erythrocyte loss. Chronic long-term therapy may result in a decrease in hemoglobin or hematocrit values.
- Use of insufficient anticoagulation can lead to thrombocyte activation and in the consequence to loss of thrombocytes. Nevertheless, despite sufficient anticoagulation, if too many thrombocytes are lost between two treatments, the treatment with the DALI system must be ceased.
- Temporary protein loss of up to 16 % per treatment may occur. This is caused by an acute reduction in fibrinogen, C-Reactive Protein (CRP) and immunoglobulins (which can lead to temporary increased susceptibility to infection) and by a dilution effect from the blood reinfusion process. Anticoagulation or reduced coagulation factors may cause a change in prothrombin time (PT) or partial thromboplastin time (PTT).
- To prevent mild anemia induced by iron/ferritin loss, regular monitoring of iron and/or other parameters of the iron metabolism (e.g. transferrin saturation, ferritin) is recommended. If necessary, iron can be supplemented after the treatment.
- The treatment frequency may influence the achievement of the treatment goals. When increasing the treatment frequency, a potentially increased blood loss and more frequent contact of the blood with foreign material must be considered.
- Ensure that patients who underwent coronary intervention are allowed to recover before the next apheresis.

Specific for citrate-containing anticoagulation

A pre-existing ionized hypocalcemia may have to be treated before initiating the DALI procedure to reduce the risk of any clinically relevant hypocalcemia during the first hour of treatment.

Fast infusion of citrate-containing anticoagulants may result in symptoms such as "tingling", usually in the patient's mouth area. In addition, muscular tension, shivering, paralytic symptoms, tightness in the chest or cold sensation may occur. This type of reaction can be caused by mild hypocalcemia because free ionic calcium has been complexed by non-metabolized citrate. Depending on the degree of the clinical signs of hypocalcemia, measures must be taken to counteract further decompensation of the calcium balance. Measures are adaptation of parameters of the active device according to the recommendations given in the table in chapter "Method of administration/Anticoagulation", and/or administration of calcium (orally or in the form of calcium gluconate (10% i.v.)) under the supervision of the physician.

CONTRAINDICATIONS

Product-specific contraindications

- Patients who currently take Angiotensin Converting Enzyme (ACE) inhibitors (e.g. enalapril, ramipril, lisinopril, etc.), as well as combination drugs which contain ACE inhibitors, or other drugs which influence the regulation of bradykinin (e.g. Angiotensin-Receptor-Nephrilysin-Inhibitors).
- Patients who suffer from thrombocytopenia ($<40 \times 10^9/\mu\text{l}$).
- Patients with known hypersensitivity to the materials used in the DALI 500 / 750 Adsorbers.
- Patients who suffer from: acute systemic infection, acute cold or influenza infections and/or severe acute cardiovascular diseases which forbid any extracorporeal therapy like DALI.

Treatment-related contraindications

- Patients who are not able to tolerate the stress of an extracorporeal treatment due to their physical status or their clinical condition.
- Patients with a history of severe coagulopathy issues such as thrombosis or hemorrhagic diathesis, are predisposed to an increased risk of bleeding potentially caused by the extracorporeal apheresis procedure and the accompanying anticoagulation.
- Patients in whom adequate anticoagulation cannot be achieved (e.g. due reduced citrate metabolism, HIT-II).

Relative treatment-related contraindications

(require a particularly careful benefit-risk evaluation before the application of apheresis):

- Patients who suffer from clinically relevant hypotension where any further drop in blood pressure may be harmful.
- Patients with impaired kidney function and reduced cardiac capacity, who may react sensitively to acute fluid and electrolyte shifts.
- Patients with confirmed or assumed intracranial disease in whom a fluid imbalance or changes in blood pressure should be avoided.
- Patients with bleeding tendency that has not been clarified.
- Patients tending to hypocalcemia due to a concomitant disease.

Note

Additional contraindications might apply with other devices used in the therapy. Please observe the Instructions for Use enclosed with the active device, tubings and solutions, etc.

PERFORMANCE CHARACTERISTICS

Performance data depend on the system configuration used. Several studies evaluating the performance of the DALI Adsorbers have demonstrated mean relative reduction values, e.g. 70 % for LDL-c, 65 % for Lp(a) and 15 % for HDL-c. These levels were obtained with the treatment system following the physician's prescription.

Please note that, besides the chosen treatment configuration, the actual reduction values also depend on the patient's blood composition, clinical settings (refer also to paragraph "Limitations of the medical device") as well as on the measurement method. For the treatment, consider that a rebound effect will occur after treatment due to the underlying hyperlipidemia.

Interactions

The DALI apheresis may alter the blood concentration of drugs by adsorption. If medically justifiable, the application time of concomitant drug therapy should therefore be adjusted to the DALI therapy (e.g. administration of drugs after apheresis treatment). The type of concomitant therapy as well as the dosage and time of drug administration must be determined by the attending physician.

All medical devices with blood-contacting, negatively charged surfaces like DALI 500 / 750 Adsorbers activate the bradykinin synthesis. Normally bradykinin is quickly degraded by ACE, but if ACE inhibitors are used in a therapy, this degradation is delayed. This may lead to increased bradykinin levels and anaphylactic reactions and even shock, which may threaten the patient's life (see Contraindications). This also applies for other drugs, which influence bradykinin regulation (synthesis or inhibition of metabolism). Patients who were previously treated with an ACE inhibitor or other drugs, which influence bradykinin regulation, must stop taking the drug early enough. In this case, the terminal half-life period of the active pharmaceutical ingredient must be taken into account (e.g. after 5-fold half-life period, a reduction to below 5 % residual concentration can be expected).

METHOD OF ADMINISTRATION

Before use

The adsorber must only be taken out of its individual packaging immediately prior to use. Properly dispose of products that have not been used.

Use the adsorbers only if:

- they are undamaged (i.e. no damage to packaging, housing incl. protective caps),
- the label is readable and expiry date not exceeded,
- the protective caps on the connectors have not become loose or dislodged,
- they have reached room temperature (20 – 25°C) to avoid unnecessary cooling down of the patient.

Anticoagulation

For the DALI treatment, citrate anticoagulation is mandatory. The ACD-A solution contains citrate, which prevents coagulation. In addition, it prevents thrombocytes and the complement system from being activated in the extracorporeal circuit. Any insufficient anticoagulation may result in an increased residual blood volume in the adsorber (after reinfusion) and thrombocyte activation. This may lead to an increased thrombocytes loss in the course of time. In case of thrombocytes loss, consider increasing the ACD-A flow to avoid thrombocytopenia. The anticoagulant requirements may vary with the patient's condition and treatment modality. It has to be prescribed by the treating physician.

The following anticoagulation parameters are recommended:

Preparation (priming/rinsing)	DALI Standard ACD-A ratio 1:20	or	DALI Low Citrate* ACD-A ratio 1:40
Before starting the treatment	Initial i.v. (unfractionated) heparin bolus: 20 IU/kg BW		
During treatment ACD-A ratio acc. to blood flow	1:(blood flow/2.7 mL/min)** Note: When changing the blood flow, adjust the ACD-A ratio according to the formula above. If less ACD-A (ACD-A flow < 2.7 mL/min) is chosen, the initial heparin bolus must be increased.		

* DALI Low Citrate can be chosen during preparation to mitigate hypersensitivity, especially during first uses of DALI Adsorbers

** Set manually by user or (if available) set automatically by ACD-Automatic mode

In case of signs of hypocalcaemia symptoms:

To take countermeasures, the ACD-Automatic mode, if selected, must be switched off. Initially, the ACD-A flow should be reduced by adapting the ACD-A flow/blood flow ratio.

⚠ Caution

The attending staff must check anticoagulation and blood flow rates at regular intervals.

The following treatment parameters are recommended:

	First treatments	Routine treatments
Blood flow	1 st to 3 rd treatment: start with low rates: 20 – 30 mL/min Depending on the patient's response, the blood flow can be slowly increased by 10 mL/min every two minutes to a max. 50 mL/min.	80 mL/min
Treated blood volume	1 st treatment: 1/3 of target blood volume 2 nd treatment: 2/3 of target blood volume 3 rd treatment: total target blood volume	Target volume: Patient's 1.6-fold blood volume
Configuration	1 st to 3 rd treatment: a single adsorber (DALI 500 or DALI 750) must be used.	Depends on the patient's clinical situation and on the prescription of the treating physician. It is recommended to use the smallest DALI Kit configuration possible required to achieve the treatment goals. In consideration of the extracorporeal volume, we do not recommend the DALI 1250 configuration for patients weighing less than 50 kg.

Note

Patients may experience hypersensitivity reactions during treatment with the DALI system. The first three treatments must be particularly monitored.

In some cases, it might be necessary to shortly interrupt the treatment if the patient shows reactions similar to hypersensitivity when the first treated blood coming from the adsorber returns to the patient. In this case, restart using low rates (20 – 30 mL/min) as described above. In some patients, hypersensitivity reactions are not mitigated after restarting. In such cases, the treatment must be terminated without returning the blood.

Patients who are excessively stressed by the blood-letting process into the extracorporeal circuit at the start of the treatment, should be supplied with a solution for fluid volume replacement through a venous indwelling cannula on their contralateral arm.

Blood flows below the recommended blood flow range can be set. The coagulation must be monitored, and the decrease in clearance values might be steeper than expected.

Intended application time

In general, the duration of the treatment depends on the blood volume of the patient and the prescription. Usually one treatment per week is carried out, resulting in approximately 50 treatments per year.

Handling

Priming and patient connection

Pre-treat the adsorber with the Priming/ACD-A solution to saturate the ligand with electrolytes and to set the pH. Add a bolus of 25,000 IU unfractionated heparin to the first primer bag. Excess heparin is rinsed out when the 2nd primer bag is connected. To connect the adsorber to the tubing lines, open the caps of the connectors and connect them to the corresponding DALI tubing lines. Refer to the IFU of the active device for the detailed instruction of the set-up. Check whether all connections have been properly established to prevent accidental separation and leaks. Fill the adsorber with DALI Priming Solution and then close the adsorber pressure line clamp directly at the pressure connector of the adsorber. Shake the adsorber to achieve a homogeneous gel bed and an adequate adsorption performance. Set the level in the drip chamber to relieve pressure from the system. Open the pressure line clamp again and rinse the adsorber with DALI Priming Solution according to the program of the active device, including set pump stops (flow rate: 400 mL/min; rinsing volume: 4800 mL).

⚠ Caution

Do not shake the adsorber after priming to avoid the gel bed from being disturbed and to prevent particles from being released. Since the system is unstable both chemically and physically, treatment must be started within 2 hours after priming. Any valid hygiene directives of the country and/or hospital in question that specify a shorter maximum durability below 2 hours after setup must be followed. Do not rinse with physiological saline solution after priming because, otherwise, the electrolytes may be desorbed. This would result in hypocalcemia in the patient. A sieve is installed in the return line downstream the DALI Adsorber(s), which serves as a second safety barrier for particles in case the housing sieve is ruptured. The same mesh size is used as in the DALI housing, therefore particles above 40 µm are retained.

Termination

Refer to the IFU of the active device for the termination of the treatment and also to the Disposal section of the present IFU.

WARNINGS AND PRECAUTIONS

Technical warnings and precautions

- Use aseptic handling methods. Improper handling of the connection sites can lead to contamination, and pathogens can enter the extracorporeal blood circuit.
- Blood and plasma may be contaminated with pathogenic agents of transmittable diseases. After having contact with the patient's blood or plasma, the adsorbers are considered to be potentially infectious. For proper disposal, please refer to the Disposal section of the present IFU.
- DALI 500 / 750 Adsorbers are intended for single use only. Re-use may be hazardous to both the patient and operator. Cleansing solutions and disinfectants should not come into contact with the product as they may damage materials employed for the housing and thus impair the product's safety and performance. DALI 500 / 750 Adsorbers are not intended for re-sterilization. If the products are re-sterilized, biological contamination may not be prevented completely. In addition, re-sterilisation may cause malfunctions.

Limitations of the medical device

The adsorbers used in the DALI system possess a limited adsorption capacity. In individual cases, efficiency may be limited in patients with an extremely high level of LDL cholesterol (e.g. patients suffering from a homozygous familial hypercholesterolemia) or Lp(a) before apheresis. In these cases, it is possible to use a larger configuration, having a higher capacity, e.g. DALI 1000 or DALI 1250.

Particular notes on materials and substances

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please use this page:
www.freseniusmedicalcare.com/en/svhc

SPECIAL PRECAUTIONS FOR STORAGE

Do not store below +5°C or above +25°C.
 Keep away from sunlight.
 Fragile, handle with care.
 The shelf life is three years (see also use-by date on the product label).

DISPOSAL

Ensure safe disposal of any unused product or waste material in accordance with local regulations. Materials that have been in contact with blood or other material of human origin may be infectious. Dispose of such materials by taking the necessary precautionary measures and in accordance with local regulations for (potentially) infectious materials.
 Contains synthetic polymer microparticles. Avoid release of the adsorber content into the environment; do not open or destroy the housing even after use.

SYMBOLS USED ON LABELS

	Sterile fluid path Sterilized using steam		Temperature limit
	Sterile barrier system with protective packaging outside		Catalog number
	Batch code		Do not reuse
	Medical Device		Unique Device Identifier
	Use-by date		Manufacturer
	Date of manufacture		Consult Instructions for Use
	Keep away from sunlight		Patient information website
	Caution		Fragile, handle with care
	Do not use if the product sterile barrier system or its packaging is compromised		

A digital version of this document is available on the following website:
www.freseniusmedicalcare.com/en/product-information

INFORMATION ON THE MANUFACTURER

Legal manufacturer

**FRESENIUS
MEDICAL CARE**

Fresenius Medical Care AG
 Else-Kröner-Str. 1, 61352 Bad Homburg, GERMANY
 +49 6172 609-0

CE 0123

The CE marking documents compliance with the applicable requirements set out by the European medical device regulations and other applicable EU legislation providing for its affixing.

Notified body with regards to the European medical device regulation:
 TÜV SÜD Product Service GmbH (0123)

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