DAI Priming Solution 3000 ml

- | | | | |
|-----------|---------------------------------|-----------|----------------------------------|
| de | Gebrauchsanweisung | tr | Kullanım Kılavuzu |
| en | Instructions for Use | lt | Naudojimo instrukcijos |
| fr | Notice d'utilisation | et | Kasutusjuhend |
| it | Istruzioni per l'uso | hu | Alkalmazási előírás |
| es | Instrucciones de uso | hr | Upute za uporabu |
| pt | Instruções de Utilização | ro | Instrucțiuni de utilizare |
| da | Brugervejledning | sk | Návod na používanie |
| sv | Bruksanvisning | sr | Uputstvo za upotrebu |
| nl | Gebruiksaanwijzing | pl | Instrukcja używania |
| cs | Návod k použití | no | Bruksanvisning |
| el | Οδηγίες Χρήσης | ru | Инструкция по применению |

Fresenius Medical Care AG

Eise-Kröner-Str. 1, 61352 Bad Homburg, GERMANY

+49 6172 609-0

www.freseniusmedicalcare.com/en/product-information

DALI Priming Solution 3000 ml

de

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en

INSTRUCTIONS FOR USE

Please read the following instructions carefully. Before use, read also through the operating instructions for the haemoadsorption device used and for the other components of the DALI system as outlined in the DALI Kit Instructions for Use (IFU), especially from the DALI Adsorber.

GENERAL INFORMATION

General description of the product

DALI Priming Solution (also called priming solution in the following) is a sterile, bicarbonate buffered solution for the initial priming of Fresenius Medical Care whole-blood adsorber system DALI.

The product is a double chamber bag of 3000 ml volume, with an injection port, a plug and a Luer-Lock connector with closure cap. The Luer-Lock connector is supposed to be connected to a tube set of a haemoadsorption device. DALI Priming Solution has been sterilized by steam.

COMPOSITION

Ion concentration in each liter of the ready-to use solution:

Na ⁺	134	mmol/l
K ⁺	4	mmol/l
Ca ²⁺	1.75	mmol/l
Mg ²⁺	0.5	mmol/l
Cl ⁻	106.5	mmol/l
HCO ₃ ⁻	36	mmol/l
pH	8.0 – 8.7	
Theoretical osmolarity	283	mOsm/l

INTENDED PURPOSE AND RELATED DEFINITIONS

Intended purpose

Rinsing of extracorporeal blood circuit.

Medical indication

Hypercholesterolemia and/or isolated increased Lp(a) leading to cardiovascular and/or peripheral vascular diseases.

Intended patient population

Adult patients with a minimum body weight of 40 kg.

Pregnancy and lactation:

There are limited clinical data available for pregnancy and no clinical data for lactation with use of this device in these populations. The device should only be used during pregnancy and lactation when the benefit to the mother clearly outweighs the potential risks to the fetus and to the infant, respectively.

Intended user group and intended use environment

Only to be used by health care professionals who are familiar with aseptic techniques and who have been trained in the proper operation and handling of the medical device. The product shall only be used in professional healthcare institutions under medical supervision.

SIDE EFFECTS

Product-specific side effects:

Possibility of electrolyte disorder in case of failure in adjustment of pH and/or electrolyte saturation of adsorber.

Treatment-related side effects:

Refer to the Instructions for Use of the DALI Adsorber or the haemoadsorption device for more information on treatment-related side effects.

Reporting of serious incidents

If any serious incident occurs in relation to the device, including those not listed in this leaflet, the treating physician shall be informed immediately.

Within the EU, the user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer according to labeling (■) and the competent authority of the EU Member State in which the user and/or patient is established.

A serious incident can be any incident that directly or indirectly leads to the death of a patient, user or other person; to the temporary or permanent serious deterioration in the health of a patient, user or other person, or a serious public health threat.

MEDICAL WARNINGS AND PRECAUTIONS (to prevent side effects)

- Mix the two compartments of each bag prior to use to prevent a wrong pH within the DALI system and thus the possibility of haemolysis.
- Do not connect the 2 DALI Priming solutions via a Y-connector to prevent infusion of unbound heparin to the patient.
- Use the primed DALI system within 2 h to avoid any shift in electrolytes or pH.
- Never rinse the primed DALI system with saline prior to use in order to keep the electrolyte saturation of the adsorber.

CONTRAINDICATIONS

Product-specific contraindications:

- Do not use for infusion / volume substitution.

Treatment-related contraindications:

Refer to the Instructions for Use of the DALI Adsorber or the haemoadsorption device for more information on treatment-related contraindications.

PERFORMANCE CHARACTERISTICS AND CLINICAL BENEFITS

Performance characteristics

DALI Priming Solution is used in combination with ACD-A solution for priming of DALI Adsorber, monitored by the haemoadsorption device for lipoprotein whole blood adsorption. The priming solution supports the pH adjustment and electrolyte saturation of DALI Adsorber before the start of therapy. For the essential performance parameters of the DALI system, please refer to the Instructions for Use of the DALI Adsorber.

Clinical benefits

DALI Priming Solution has no clinical benefit on its own. The clinical benefit of the whole blood adsorption system DALI for patients with hypercholesterolemia/isolated Lp(a) increase includes the reduction of the incidence rate of cardiovascular events due to reduction of low-density lipoprotein-cholesterol (LDL-c) and Lp(a). Please refer to the Instructions for Use of the DALI Adsorber for further information.

INTERACTIONS

- No drug other than heparin shall be injected into the priming solution.
- For interactions of the DALI system, please refer to the IFU of the DALI Adsorber.

METHOD OF ADMINISTRATION

Dosing

For preparation of DALI Adsorber, usually the mixed solution from two bags of DALI Priming Solution is applied.

For information about flow speed and mixing ratio with ACD-A, please refer to the instructions for use of the adsorber kit or the haemoadsorption device.

Handling Instructions

Before use

- Inspect the bag
- Only use a bag if the solution is clear, the peel seam is intact and the packaging is undamaged.

- Labeling must be clear and the expiry date must not be exceeded.
- Recommended bag temperature for application is +20 °C to +25 °C maximum (to avoid hypothermia).

Handling


Opening of the overwrapping

- The double-chamber bag containing the solution must only be taken out of its overwrapping immediately prior to its use. Before opening, the overwrapping must be checked for damage. Bags with damaged overwrapping must be discarded.
- Remove the outer packaging from the priming bag and put the priming bag on a firm surface.
- The peel seam that separates the two compartments must be intact after unpacking, otherwise compliance of the solution specification cannot be guaranteed and the bag must be discarded.


Mixing of both compartments

Both compartments must be mixed immediately before use. Do not use unmixed. Mixing of the chambers before use is mandatory, as application of solution from the left chamber only can lead to haemolysis.

Preparation of mixed solution: Please follow steps 1– 4.

- 1 Unfold middle peel seam.
- 2 Roll the bag diagonally towards the bag connector. Middle peel seam will open.

- 3 Mix the solution by pushing bag with both hands twice.
- 4 Make sure that the peel seam is completely open.

After mixing both compartments, a check must be carried out to verify that the peel seam is completely open, that the solution is clear and colorless and that the bag is not leaking. Use only if the solution is clear and the container is not damaged. A leaking bag must be discarded immediately.

Addition of heparin:

Add a bolus of 25,000 I.U unfractionated heparin to the first priming bag after mixing of both compartments using the injection port.

Injection port

Use of the injection port is the responsibility of the treating physician and it can be used with established cannula not exceeding 23 Gauge.

Disinfect the injection port with 70% ethanol and let it dry before application of cannula.

Connection of the bag to extracorporeal circuit

The bag with the mixed solution is connected to the tube set of the haemoadsorption device according to the instructions for the devices. The bag is sterile before use and the Luer-Lock connector with cap must not be disinfected before use. Contamination of the solution or any parts in contact with the solution must be avoided. Any bag of DALI Priming Solution that has been disconnected from the tube set must be

DALI Priming Solution 3000 ml

en

discarded. Depending on the haemoadsorption device in use, the change of the first bag to the second one must either be performed manually by the user or it is initiated by the device.

Termination

After priming, disconnect the tubings from the bag, then the DALI system can be connected to the patient according to the Instructions for Use of the haemoadsorption device.

WARNINGS AND PRECAUTIONS

Technical Warnings and Precautions

- Use aseptic techniques during preparation and in connection to extracorporeal blood circuit.
- Use the mixed priming solution within 30 minutes after mixing of compartments, otherwise compliance of the solution specification cannot be guaranteed and the bag must be discarded.
- Application of disinfectants can lead to damages to the Luer-Lock connector with cap.
- Do not re-use.
- Any unused portion of solution is to be discarded.

Particular notes on materials and substances

There is no hypersensitivity to the bag system or any of the materials known.

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH"), please use this page: www.freseniusmedicalcare.com/en/svhc

SPECIAL PRECAUTIONS FOR STORAGE

Keep away from sharp objects such as knives and scissors. Do not store below +5 °C or above +25 °C. Keep dry.

DISPOSAL

Ensure safe disposal of any unused product or waste material in accordance with local regulations.

Materials that have been in contact with blood or other material of human origin may be infectious. Dispose of such materials by taking the necessary precautionary measures and in accordance with local regulations for (potentially) infectious materials.

SYMBOLS USED ON LABELS

	Medical Device
	Manufacturer
	Catalog number
	Batch code
	Date of manufacture
	Use-by date
	Unique Device Identifier

	Units
	Caution
	Consult instructions for use
	Patient information website
	Do not use if package is damaged and consult instructions for use
	Sterilized using steam
	Single sterile barrier system with protective packaging outside
	Do not re-use
	Non-pyrogenic
	No presence of natural rubber latex
	Temperature limit

The CE marking documents compliance with the applicable requirements set out by the European medical device regulations and other applicable EU legislation providing for its affixing. Notified body with regards to the European medical device regulation: TÜV SÜD Product Service GmbH (0123)

CE 0123

DATE OF REVISION OF TEXT

2023-09

A digital version of this document is available on the following website:

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