# Filtrate bag 10 L (CN)

## INSTRUCTIONS FOR USE

Please read the following instructions carefully.

## GENERAL INFORMATION

### General description of the product

The filtrate bag is an unsterile product; it consists of tubes, connectors, cap, clamp and fluid bag.

The product is intended to be used in combination with an extracorporeal acute blood purification system multiFiltrate and multiFiltratePRO.

## COMPOSITION

Materials: Tube: medical grade soft-PVC; Connectors and other components: PP, PVC, ABS, PE.

## INTENDED PURPOSE

### Medical purpose

This product is used as a container to collect the clinical waste such as dialysate/body fluids at hospitals/clinics in combination with Fresenius Medical Care CRRT (Continuous Renal Replacement Therapy) devices multiFiltrate and multiFiltratePRO during CRRT treatment, haemoperfusion procedures or therapeutic plasma exchange.

This product is used as ancillary device for collecting waste fluid during HD/F treatment.

### Intended patient population

The bags have been specified by the manufacturer for the purpose of treating patients irrespective of their age and body weight, under consideration of the specified technical data of the medical device and the single-use items used (e.g., delivery rates, fill volumes).

### Intended user/environment

The disposable may only be used by individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown.

Operation in suitable rooms in professional health care facilities. Normative and local regulations must be observed.

## SIDE EFFECTS

No specific side effect due to the use of Filtrate bag is known.

Side effects related to the multiFiltrate and multiFiltratePRO must be taken into account. Please refer to the Instructions for Use (IFU) of multiFiltrate and multiFiltratePRO for more information on the side effects.

### Reporting of serious incidents

If any serious incident occurs in relation to the device, including those not listed in this leaflet, the treating physician shall be informed immediately. Within the EU the user must report any serious incident that has occurred in relation to the device to the manufacturer according to labelling (📄) and the competent authority of the EU Member State in which the user is established.

A serious incident can be any incident that directly or indirectly leads to the death of a patient, user or other person; to the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health treat.

## CONTRAINDICATIONS

There are no contraindications known for Filtrate bag 10L.

For general contraindications on extracorporeal blood purification therapies, please refer to the Instructions for Use of multiFiltrate and multiFiltratePRO.

## PERFORMANCE CHARACTERISTICS AND CLINICAL BENEFITS

Filtrate Bag 10L has to be used in combination with multiFiltrate and multiFiltratePRO during CRRT treatments haemoperfusion procedures and therapeutic plasma exchange.

## METHOD OF ADMINISTRATION

### Handling Instructions

#### Before use

Refer to the Instructions for Use of the multiFiltratePRO and multiFiltrate device.

## Handling

The product is intended to be used in combination with an extracorporeal acute blood purification system multiFiltrate and multiFiltratePRO.

Open the package to get the product, hang filtrate bag onto the hooks of the lower scales, remove the luer cap, connect the male Luer-Lock with the Filtrate Line (multiFiltrate filtrate line "F") according to Instructions for Use for multiFiltrate and multiFiltratePRO device.

Upon completion of solution collection, enable the function "bag change" and close the clamp, according to Instructions for Use for multiFiltrate and multiFiltratePRO System.

Discard the used product and material waste as described in the "Disposal" Section of this IFU.

Where present, the available stopcock can be used to empty the bag after disconnection from the multiFiltrate and multiFiltratePRO device.

### Intended application time

The maximum application time is 72 hours (equal to the maximum application time of multiFiltrate and multiFiltratePRO tubing systems for extracorporeal acute blood purification).

### Termination

Refer to the Instructions for Use of the multiFiltratePRO and multiFiltrate device for the termination of the treatment and also to the "Disposal" Section of the present Instructions for Use.

## WARNINGS AND PRECAUTIONS

### Warnings

Do not use the product if there is damage to the bag or components.

If the carton is damaged, check the products contained carefully. Do not use if the package is damaged, if the protective or closure caps are not in place, or if there is any visible damage on the finished products (e.g. kinked tubes).

Do not use one filtrate bag for different patients.

Do not use after use-by date (refer to label).

Inspect the extracorporeal circuit for kinking and leaks during the priming phase and treatment, taking corrective measures or exchanging the disposable as necessary.

The system contains large-size packaging and foils that can cause choking, and which should be kept away from children.

### Particular notes on materials and substances

There is no hypersensitivity to the bag system or any of the materials known.

One or more components contain in a concentration > 0.1 mass% according to Article 33 and 59 (1, 10) of Regulation (EC) No. 1907/2006 ("REACH"): – 2-ethylhexyl 10-ethyl-4,4-dioctyl-7-oxo-8-oxa-3,5-dithia-4-stannatetradecanoate (DOTE)

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please use this page: [www.freseniusmedicalcare.com/en/svch](http://www.freseniusmedicalcare.com/en/svch).

## SPECIAL PRECAUTIONS FOR STORAGE

Follow the indication of the product label. Protect from moisture, freezing and excessive heat.

## DISPOSAL

The dialysate waste may be contaminated with pathogenic agents of transmittable diseases; thus, the bag must be considered to be potentially infectious.

The accredited organization is liable for the safe disposal of the product. Any unused product or waste material shall be disposed of in accordance with local requirements and in accordance with local requirements for potentially contaminated materials.

## SYMBOLS USED ON LABELS

	Medical Device
	Unique Device Identifier
	Patient information website
	Latex-free
	Exchange the Bloodline/Tubing System after maximum usage time or pumped blood volume as indicated on the primary package
	Units

## WARRANTY

The manufacturer shall not be liable for any misuse, improper handling, non-compliance with instructions for use and cautionary notes and for any damage incurred subsequent to the manufacturer's delivery of the bloodlines.

## DATE OF REVISION OF TEXT

01/2021

For the electronic version of instruction for use (e-IFU), please use this page:

[www.freseniusmedicalcare.com/en/product-information](http://www.freseniusmedicalcare.com/en/product-information).

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Fresenius Medical Care AG & Co. KGaA