# Therapeutic Plasma Exchange Tubing

## MODELS:

multiFiltratePRO TPE

## INSTRUCTIONS FOR USE

Please read the following instructions carefully.

## GENERAL INFORMATION

### General description of the product

The therapeutic plasma exchange tubing (multiFiltrate PRO TPE) is intended to be used and compatible with an extracorporeal blood purification system multiFiltratePRO device.

multiFiltratePRO TPE is intended for single use.

multiFiltratePRO TPE is designed to provide extracorporeal blood treatment of patients who need therapeutic plasma exchange due to various diseases.

multiFiltratePRO TPE is composed by the arterial line / access line (RED line), venous line / return line (BLUE line), filtrate line (YELLOW line) and plasma line (WHITE line).

The product is pyrogen free.

### Sterilization

multiFiltratePRO TPE is ETO sterilized.

## COMPOSITION

Tubes: medical grade soft-PVC.

Connectors and other components: Acrylonitrile butadiene styrene (ABS), Polybutylene Terephthalate (PBT), Polycarbonate (PC), Hydrophobic Acrylic Copolymer membrane, Polypropylene (PP), Polytetrafluoroethylene (PTFE), Polyvinylchloride (PVC), Silicone, Polyisoprene (PHE), Thermoplastic elastomer (TPE).

## INTENDED PURPOSE AND RELATED DEFINITIONS

### Intended purpose

Channeling of blood and fluid in an extracorporeal treatment.

### Medical indication

Conditions requiring the exchange of blood plasma via TPE.

### Intended patient population

The devices have been specified by the manufacturer for the purpose of treating patients with a body weight of 40 kg and more, irrespective of their age, under consideration of the specified technical data of the device and the single-use items used (e.g., delivery rates, fill volumes).

### Intended user group and intended environment

The disposable must only be used by individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown.

Operate in rooms suitable for dialysis located in professional healthcare facilities.

Normative and local regulations must be observed.

## SIDE EFFECTS

Occasional occurrence of the following side effects is reported in current literature:

- Hypotension
- Blood loss (e.g. hypovolemia, and hypovolaemic shock, hypotension, anaemia, cardiac arrhythmia, cardiac arrest)
- Coagulation disorders (e.g. thrombocytopenia)
- Embolism (air-, thrombo-)
- Bacterial infection
- Hypothermia
- Haemolysis
- Hypersensitivity or hypersensitivity-like reactions (e.g. dyspnea, cardiopulmonary arrest, hypotension, urticaria, flushing, itching, abdominal pain, nausea, convulsions and unconsciousness).

Additional side effects might be specific for other devices and drugs used in the therapy.

Refer to the multiFiltratePRO device Instructions for Use for more information on treatment related side effects.

### Reporting of serious incidents

If any serious incident occurs in relation to the device, including those not listed in this leaflet, the treating physician shall be informed immediately. Within the EU the user must report any serious incident that has occurred in relation to the device to the manufacturer according to labelling (M) and the competent authority of the EU Member State in which the treatment is performed.

A serious incident can be any incident that directly or indirectly leads to the death of a patient, user or other person; to the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health threat.

## CONTRAINDICATIONS

### Product related contraindications

Hypersensitivity to any of the material of the tubings as listed under composition or to ETO residues.

In severe cases the treatment must be discontinued, the blood should not be returned to the patient, and appropriate emergency medical treatment should be initiated.

### Therapy related contraindication

Refer to the multiFiltratePRO device Instructions for Use for more information of general contraindications for extracorporeal blood purification.

## PERFORMANCE CHARACTERISTICS

Refer to the multiFiltratePRO device Instructions for Use for essential performance parameters.

## METHOD OF ADMINISTRATION

### Handling Instructions

Refer to the multiFiltratePRO device Instructions for Use of the regarding how to handle the disposable during set-up, priming, treatment and reinfusion. Also, the Instructions for Use of other disposables such as the plasma filter must be taken into account.

### Preparation

Refer to the multiFiltratePRO device Instructions for Use for details on any preparatory handling of the disposable before it is ready for use or during its use and information to verify if the disposable is properly installed.

multiFiltratePRO TPE is intended to be used and compatible with multiFiltratePRO device (refer to label) and shall be used only after appropriate instruction or training.

The infusion administration lines should be clamped, except when needed.

- Unpack and connect the multiFiltratePRO TPE aseptically without touching open connectors.
- Colour codes should be followed and used in line with the corresponding markings on the device.
- Handle pumps and pump clips in accordance to the Instructions for Use of the device.
- Tighten all closure caps and ensure that all connectors and caps are secured.
- Fill and rinse the disposable in accordance with the Instructions for Use of the device or supplemental Instructions for Use for additional options, and training, if applicable.

Fresenius Medical Care disposable is designed to withstand the maximum and minimum manufacturer's recommended pressures and flow rates generated in use with the respective multiFiltratePRO device. All the other relevant specifications of the device apply.

### Treatment

To ensure correct functionality of the bubble catchers, fill to about 1 cm below the top. This avoids fluid contacting the hydrophobic filter as well

as air entering the extracorporeal circuit. In the event that fluid reaches and comes into contact with the hydrophobic filter within the Transducer Protector (TP) the disposable must be replaced.

- The disposable has an additional connection intended to be connected to the pressure measurement system. This can be used to connect a new pressure measurement line (available as an accessory in the Fresenius Medical Care product range).
- Do not push back the fluid in the pressure measurement line by using a syringe. This might damage the membrane in the TP and thus lead to contamination. In the event that fluid might have passed through the TP, check the device for contamination after treatment is completed. If contamination has occurred, the device must be taken out of service and disinfected according to the manufacturer's recommendations before further use.
- After the manipulation of lines or the usage of components during the treatment check and, if necessary, restore the correct position of the lines and components.
- Disinfect the corresponding access sites without protection cap prior the connection with other products.

### Reinfusion

- Refer to the multiFiltratePRO device Instructions for Use for the termination of the treatment and also to the "Disposal" Section of the present Instructions for Use.
- Close all clamps on the multiFiltratePRO TPE before removing the disposable to reduce the risk of fluid leakage.

## WARNINGS AND PRECAUTIONS

### Warnings

- The disposable must only be used by individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown.
- The multiFiltratePRO TPE are intended for single use only. The correct function of all interfaces is ensured only for single use. Re-use may be hazardous to both patient and operator (e.g. impaired performances, contamination).
- Do not use on and after expiry date (refer to label).
- If the carton is damaged, check the products contained carefully. Do not use if the sterile package is damaged, if the protective or closure caps are not in place, or if there is any visible damage on the finished products (e.g. kinked tubes).
- Make sure all tubes and chambers are correctly inserted into the respective holders. Avoid kinking or occluding the disposable in order to avoid mechanical and chemical damage to cellular blood constituents.
- For hygienic and functional reasons, it is recommended that tubing systems are inserted immediately before preparation only, thereby keeping the preparation and circulation times nearer the commencement of the treatment and in compliance with applicable guidelines.
- Connect the multiFiltratePRO TPE aseptically without touching open connectors.
- Disinfect with 70% alcohol access sites without protection cap prior to the connection with other products to reduce the risk of infection and let it dry before the connection.
- The compatibility of disinfectants (other than those recommended) with access sites shall be determined prior to clinical use.
- Excessive negative pressure may cause partial collapse of the pump segment resulting in an actual blood flow substantially less than indicated on the device.
- Minimum temperature of use of the disposable is 18°C (64°F).
- Do not use needles with a beveled tip and a diameter larger than 20 gauge (outer needle diameter to be 0.9 mm or less) to puncture injection sites (if present).
- To ensure a secure connection between patient access and blood line, hold and screw the colored (blue, red) coupling nut on the blood line only. Do not apply the screwing torque to the inner part of the connector. After the connection, check that the components are firmly screwed together.
- Make sure that the components/tubes in direct interface with the device are fitting properly and ensure that all caps and all screw-lock joints are secure and properly tightened (particularly those of the patient connections, the dialyzer connections, the device and the solution bags connections).
- Inspect the disposable for kinking and leaks during the priming and treatment phase, taking corrective measures (e.g. tightening Luer-Lock connection) or exchanging it as necessary.
- To avoid air embolism, make sure that the return patient line is correctly inserted into the device's air bubble detector.
- Careful attention must be given to the installation and blood level of the bubble catchers. Make sure that due to pressure fluctuation neither the blood contacts the hydrophobic filter nor air enters the extracorporeal circuit.
- Change the disposable if a visible clot formation („white strip“) is seen in the blood line downstream of the return line bubble catcher, for instance at the point where the calcium infusion line merges into the blood line.
- Cleansing solutions and disinfectants may damage materials employed for the disposable. Safety and performance of use can no longer be guaranteed, and the manufacturer assumes no liability.
- The plastics used can be incompatible with drugs or disinfectants (e.g. connectors made of polycarbonate can develop cracks when in contact with aqueous solutions with pH > 10).
- If nutritional solutions are administered into the blood line, wetting of the luer-lock connection to the feeding system with lipidic fluids can weaken the properties of the plastic material used. Make sure that during the connection of the feeding system, the infusion line close to the luer-lock connection site remains completely free of nutrition solution.

### Operating time

- The maximum application time is 72 hours. The disposable must be replaced after the maximum usage time or respectively after a specific blood volume pumped as indicated on the primary packaging.

### Particular notes on materials and substances

#### CMR substances and endocrine-disrupting substances

For SVHC information according to Article 33 of Regulation (EC) No. 1907/2006 ("REACH") please use this page:  
www.freseniusmedicalcare.com/en/svhc

### Special precautions for storage

Follow the indication of the product label. Protect from moisture, freezing and excessive heat.

### Disposal

Ensure safe disposal of any unused product or waste material in accordance with local regulations. Materials that have been in contact with blood or other material of human origin may be infectious.

Dispose of such materials by taking the necessary precautionary measures and in accordance with local regulations for (potentially) infectious materials.

## SYMBOLS USED ON LABELS

	Medical Device
	Unique Device Identifier
	Patient information website
	Latex-free
	Sterile, Sterilized using ethylene oxide
	Replace the Bloodline after maximum usage time or pumped blood volume as indicated on the primary package
	Pump segment diameter / length
	Blood Priming Volume
	Units
	Temperature limit
	Do not use if package is damaged
	Refer to the instructions
	Non-pyrogenic
	Sterile barrier system
	Do not re-use
	Caution (Consult the Instructions for Use for important cautionary information such as warnings and precautions)
	Batch code
	Catalogue number
	Date of manufacture
	Expiry date

## WARRANTY

The manufacturer shall not be liable for any misuse, improper handling, non compliance with instructions for use and cautionary notes, and for any damage incurred subsequent to the manufacturer's delivery of the multi-FiltratePRO TPE.

## DATE OF REVISION OF TEXT

2023/07

For the electronic version of instruction for use (e-IFU), please use this page:  
www.freseniusmedicalcare.com/en/product-information.

## INFORMATION ON THE MANUFACTURER

### Legal manufacturer

**FRESENIUS  
MEDICAL CARE** CE 0123

**Fresenius Medical Care AG**  
Else-Kröner-Str. 1, 61352 Bad Homburg

GERMANY ☎ +49 6172 609-0

www.freseniusmedicalcare.com/en/product-information

**Manufacturer/Registrant:** Fresenius Medical Care AG

**Address of Manufacturer/Registrant:** Else-Kröner-Str. 1, 61352 Bad Homburg, GERMANY

**Manufacturing Address(es):** No: 16, Liman Serbest Bölgesi Mahallesi, 07070 Antalya, TURKEY

**Agent / After-Sale Service Provider in China:** Fresenius Medical Care (Shanghai) Co., Lt

**Address of Agent/After-Sale Service Provider:** Section A, Building 01, No. 439 Futexi Road 1st, China (Shanghai) Pilot Free Trade Zone

**Telephone Number of Agent/After-Sale Service Provider:**

021-61152800

**Registration Certificate No./Product Technical Requirements**

No.: GuoXieZhuJin20223100504

**Date of Manufacture:** See outer package

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