

# **Eco Chair Serie II – a good way to take care**

Actualway Product Portfolio



**FRESENIUS  
MEDICAL CARE**

# Actualway Eco Serie II Therapy Chair

## The manufacturer

Actualway is a Portuguese manufacturer of medical equipment for care & treatment such as Therapy, Sampling/Blood Donation, Gynaecology, Multi-functional Chairs, Couches and many accessories.

Fresenius Medical Care cooperates since 2018 for offering functional, comfortable and economic dialysis treatment chairs with high stakeholder satisfaction in Europe and abroad.

Actualway has been increasing its distribution all over the world thanks to their promotion efforts. Actualway's sales strategy is based on building a long term relationship with all their partners, providing quality and service.

## Actualway certification

Actualway products are manufactured in Portugal, tested according to applicable norms.

Actualway is certified

- QM-certification according to ISO 9001
- QM- certification according to EN ISO 13485

**The Actualway vision:** The best way to take care is to get our partner's and patient's satisfaction meeting their needs.

**The Actualway promise:** We produce reliable, strong and comfortable products, because we do care.

We provide quality, after-sales service, reliable products adapted to the needs of patients and users, and an excellent delivery time. Come and join us, so that we can expand business together!



## Features and advantages

- 3 motors
- Trendelenburg position ( $-12^\circ$ ) in 12 sec
- Electric height adjustment from 53 to 78 cm
- Back rest and leg rest adjustment
- Cervical pillow
- Lockable handset
- Manual foot rest adjustment
- Vertical and lateral arm rest adjustment
- 4 castors with 10 cm diameter and individual brake
- Push handle
- Equipotential bonding cable



Lockable handset



Vertical and lateral arm adjustment



Manual foot rest adjustment




Vertical and lateral arm adjustment

## Options

- Memory foam upholstery
- Paper roll holder
- IV serum pole (Stainless/Epoxy)
- ABS table
- Foot rest motor

| Technical specifications:          |                               |
|------------------------------------|-------------------------------|
| Power supply                       | 100 V–240 V; 50/60 Hz         |
| Protection                         | IPX4                          |
| Handset and control box protection | IPX6                          |
| Maximum load capacity              | 250 kg + optional accessories |
| Logistic data:                     |                               |
| Net weight                         | 80 kg                         |
| Gross weight                       | 110 kg                        |
| Pallet                             | 80 x 120 x 160 cm             |
| Dimensions:                        |                               |
| Width                              | 55 cm                         |
| Total width                        | 88 cm                         |
| Height                             | 53–78 cm                      |
| Length                             | 198–215 cm                    |
| Frame                              | 80 x 19 x 117 cm              |

# Colour options and Article numbers




Snow




Olive




Canard




Myosotis



Sapphire




Sand




Pepper



Brown




Macadam




Mist




Light Grey




Pottery



Cherry



Garnet



Coal

## Upholstery colour

### Actualway Dialysis Chair Eco

F00008506

#### Actualway Therapy Chair (Serie II) Eco

(3 independent motors: height, back rest and leg rest adjustment, foot rest manual adjustment, lockable handset, Trendelenburg position (handset operated), vertical and horizontal armrest adjustment, 4 castors with individual locking brake (10 cm diameter), weight capacity: 250 kg + extra accessories, equipotential bonding cable, pillow and foot rest in upholstery colour)

Myosotis

### Actualway Therapy Chair Eco Options (new orders only)

F00008520

#### Actualway Eco IV serum pole

F00008521

#### Actualway Eco ABS foldable table

F00008522

#### Actualway Eco additional 4th motor for foot rest adjustment

F00008523

#### Actualway Eco paper roll holder

### Actualway Therapy Chair Eco Instruction for Use (IFU)

Please note: If not otherwise indicated dialysis chairs will be delivered with IFU in English language

F00008570

#### Actualway Therapy Chair Eco IFU English

F00008571

#### Actualway Therapy Chair Eco IFU Portuguese

F00008572

#### Actualway Therapy Chair Eco IFU French

F00008573

#### Actualway Therapy Chair Eco IFU Spanish

#### Actualway Therapy Chair Eco IFU in other language (Please indicate language)

# Article numbers

## Actualway Therapy Chair Colour Options (new orders only)

|           |   |
|-----------|---|
| F00008527 | <b>Actualway upholstery colour Snow</b>       |
| F00008528 | <b>Actualway upholstery colour Sand</b>       |
| F00008529 | <b>Actualway upholstery colour Brown</b>      |
| F00008530 | <b>Actualway upholstery colour Olive</b>      |
| F00008531 | <b>Actualway upholstery colour Macadam</b>    |
| F00008532 | <b>Actualway upholstery colour Pottery</b>    |
| F00008534 | <b>Actualway upholstery colour Cherry</b>     |
| F00008538 | <b>Actualway upholstery colour Garnet</b>     |
| F00008540 | <b>Actualway upholstery colour Myosotis</b>   |
| F00008541 | <b>Actualway upholstery colour Canard</b>     |
| F00008544 | <b>Actualway upholstery colour Sapphire</b>   |
| F00008546 | <b>Actualway upholstery colour Light Grey</b> |
| F00008548 | <b>Actualway upholstery colour Pepper</b>     |
| F00008550 | <b>Actualway upholstery colour Coal</b>       |
| F00008625 | <b>Actualway upholstery colour Mist</b>       |
|           |   |
| F00008556 | <b>Actualway pillow colour Snow</b>           |
| F00008557 | <b>Actualway pillow colour Sand</b>           |
| F00008558 | <b>Actualway pillow colour Brown</b>          |
| F00008559 | <b>Actualway pillow colour Olive</b>          |
| F00008560 | <b>Actualway pillow colour Macadam</b>        |
| F00008561 | <b>Actualway pillow colour Pottery</b>        |
| F00008562 | <b>Actualway pillow colour Cherry</b>         |
| F00008562 | <b>Actualway pillow colour Garnet</b>         |
| F00008563 | <b>Actualway pillow colour Myosotis</b>       |
| F00008565 | <b>Actualway pillow colour Canard</b>         |
| F00008566 | <b>Actualway pillow colour Sapphire</b>       |
| F00008567 | <b>Actualway pillow colour Light Grey</b>     |
| F00008568 | <b>Actualway pillow colour Pepper</b>         |
| F00008569 | <b>Actualway pillow colour Coal</b>           |
| F00008626 | <b>Actualway pillow colour Mist</b>           |
|           |   |
| F00008525 | <b>Actualway footrest = pillow colour</b>     |
| F00008526 | <b>Actualway foam armrest = pillow colour</b> |

**Manufacturer:**

Via Actual – Comércio Internacional, Lda.  
CE Castelo da Maia  
Rua Manuel Assunção Falcão, N°337, Armazém 13  
4475-041 Santa Maria de Avioso – Maia  
Portugal  
Tlf./Phone: +351 229 829 663/4  
Fax: +351 229 829 665  
E-mail: [comercial@actualway.pt](mailto:comercial@actualway.pt)  
[www.actualway.pt](http://www.actualway.pt)



**FRESENIUS  
MEDICAL CARE**

Head office: Fresenius Medical Care Deutschland GmbH · 61346 Bad Homburg v. d. H. · Germany  
Phone: +49 (0) 6172-609-0 · Fax: +49 (0) 6172-609-2191  
[www.FreseniusMedicalCare.com](http://www.FreseniusMedicalCare.com)