



# multiFiltratePRO HL7 connectivity Efficient and secure device interoperability

In the critical care setting, device connectivity to electronic medical records is designed to enhance the quality of patient care and healthcare operations.

# Access the benefits of data integration with multiFiltratePRO HL7 connectivity:

#### **Data-based clinical decision making**

Accurate and faster electronic documentation<sup>1</sup> makes information available for data-driven clinical decisions

#### **Streamline workflows**

Staff spends less time manually documenting<sup>1</sup>, which may help focus on what matters most – patient care


### **Reduce errors**

Electronic documentation reduces the risk of transcription errors<sup>1</sup>





Head office: Fresenius Medical Care Deutschland GmbH 61346 Bad Homburg v. d. H. · Germany · Phone: +49 (0) 6172-609-0 www.freseniusmedicalcare.com



#### **One connection**

Connect with leading Patient Data Management Systems (PDMS) and middleware providers using HL7 standard communication protocol



#### Secure

Transmit information through current data protection and security standards, including modern encryption protocols



#### **Support and training**

Count on our support, from planning and installation to future updates

#### **Network specifications**

Ethernet: 10/100 Mbps, RJ-45 port		
802.1X Ethernet port security	Device operating in client and host mode	TLS1.2* encryption
Static IP address or DHCP*	Supporting NTP*	HL7 2.6

\*Abbreviations: IP: Internet Protocol, DHCP: Dynamic Host Configuration Protocol, NTP: Network Time Protocol, TLS: Transport Layer Security

#### multiFiltratePRO-HL7-connector

Description	Article no.
Built-in connectivity option to support HL7	On request

## Connect with us To connect multiFiltratePRO to your PDMS, contact your local Fresenius Medical Care representative